

REF	500-3	6×3mL	EXP	2025-05-01	LOT	2FB6E4C	更新日期：	2024-05
1 / 3								

批号：2FB6E4C					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(25-OH-Vit D)25羟基维生素D					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	34.6	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	86.6	12%	<5	<100
(Cortisol)皮质醇					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	430	35.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	388	6.6%	13	293
(C-P)C-肽					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.2	6.8%	6	136
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.66	9.8%	<5	<100
(E2)雌二醇					
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	2228	15.7%	6	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1462	4.9%	6	137
Beckman Access Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	228	4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1760	4.5%	13	300
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	3238	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1286	12%	<5	<100
(FA)叶酸					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.97	12%	<5	<100
(FSH)促卵泡刺激素					
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	57.2	4.7%	6	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	62.5	5.9%	6	136
Beckman Access Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	7.49	6.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	82.7	6.3%	13	316
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	42	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	37.4	12%	<5	<100
(FT3)游离T3					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	14.6	4.3%	6	136
OTHER TELLGEN \ Chemiluminescence Immunoassay (CLIA)	pmol/L	21.9	7.9%	5	139
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	16.5	18.8%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	2.38	9.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	9.94	5.8%	13	300
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.1	12%	<5	<100
GeteinBiotech MAGICL Series Getein Biotech \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.46	1.1%	<5	<100
(FT4)游离T4					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	41.9	4.5%	6	136
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	33.3	11.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	36.5	3.7%	13	307
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	40.5	12%	<5	<100
GeteinBiotech MAGICL Series Getein Biotech \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.8	2.1%	<5	<100
(hCG)人绒毛膜促性腺激素					
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	690	4.5%	6	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	508	12%	<5	<100
(INS)胰岛素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	89.2	16.3%	<5	<100

REF	500-3	6×3mL	EXP	2025-05-01	LOT	2FB6E4C	更新日期：	2024-05
2 / 3								

		批号：2FB6E4C			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(LH)促黄体生成素					
Roche cobas e 801	IU/L	56.9	4.2%	7	136
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Ortho Vitros 3600/5600	IU/L	63.8	6.6%	6	<100
OTHER \ Chemiluminescence Immunoassay (CLIA)					
Beckman Access	IU/L	2.04	6.2%	<5	<100
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman UniCel Dxl 600/800	IU/L	48.8	3.1%	13	307
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Autolumo A2000/A2000 PLUS	IU/L	48	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	IU/L	18.6	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(PRL)催/ 泌乳素					
Ortho Vitros 3600/5600	ng/mL	35.5	6.3%	6	<100
OTHER \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	48.9	6.5%	6	136
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman Access	ng/mL	6.53	9.9%	<5	<100
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman UniCel Dxl 600/800	ng/mL	35.8	3.9%	13	308
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Autolumo A2000/A2000 PLUS	ng/mL	67.1	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	ng/mL	45.3	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(PROG)孕酮					
Ortho Vitros 3600/5600	ng/mL	16.2	5.7%	6	<100
OTHER \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	14.1	6%	6	137
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman Access	ng/mL	0.929	6.8%	<5	<100
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman UniCel Dxl 600/800	ng/mL	14.6	8.5%	13	323
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Autolumo A2000/A2000 PLUS	ng/mL	26.7	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	ng/mL	22.9	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(TES)睾酮					
Ortho Vitros 3600/5600	ng/mL	13.5	6%	6	<100
OTHER \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	8.51	7.4%	6	136
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman Access	ng/mL	0.512	8.9%	<5	<100
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman UniCel Dxl 600/800	ng/mL	7.45	5.6%	13	302
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Abbott ARCHITECT i series	ng/mL	7.79	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	12	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(TSH)促甲状腺素					
Roche cobas e 801	mIU/L	26.1	5.4%	6	131
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
OTHER	mIU/L	34.6	3.7%	5	138
TELLGEN \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	mIU/L	21.5	28.6%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	mIU/L	27.8	6.3%	13	300
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Autolumo A2000/A2000 PLUS	mIU/L	22.9	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(TT3)总三碘甲状原氨酸					
Autolumo A2000/A2000 PLUS	ng/mL	4.01	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
GeteinBiotech MAGICL Series	ng/mL	1.86	2.6%	<5	<100
Getein Biotech \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	ng/mL	2.99	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(TT4)总T4					
Abbott ARCHITECT i series	ng/mL	105	12%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	126	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
GeteinBiotech MAGICL Series	ng/mL	45.5	2.9%	<5	<100
Getein Biotech \ Chemiluminescence Immunoassay (CLIA)					

REF	500-3	6×3mL	EXP	2025-05-01	LOT	2FB6E4C	更新日期：	2024-05
3 / 3								

		批号：2FB6E4C			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(VB12)维生素B12					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	737	4.5%	7	137
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	608	5.4%	13	315
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	729	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。