

REF	530-3	6×3mL	EXP	2026-10-03	LOT	I293576	更新日期：	2024-05
1 / 1								

		批号: I293576			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CK-MB mass)肌酸激酶-MB质量					
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	36.7	10%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	57.1	4.9%	<5	<100
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	ng/mL	93	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	33.7	3.1%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	mg/L	0.973	10%	<5	<100
Roche ACCU-CHEK series Roche \ Nephelometry	mg/L	1.53	9.8%	<5	<100
(hs-CRP)超敏C反应蛋白					
Beckman AU series Beckman Coulter \ Immunturbidimetry (ITA)	mg/L	7.78	10.3%	<5	<100
(Myo / Mb)肌红蛋白					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μg/L	366	2.4%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	397	10%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	376	4.1%	<5	<100
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	μg/L	382	10%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	pg/mL	3230	10%	<5	<100
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	pmol/L	380	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1320	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	155	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	1382	2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	162	2%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	1262	5.8%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	148	5.8%	<5	<100
(Tnl)肌钙蛋白I					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	73	2.3%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.8	3.1%	<5	<100
Radiometer AQTSeries Radiometer \ Time-resolved fluorescence immunoassay (TRFIA)	ng/mL	2.1	10%	<5	<100
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	76.6	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。