

REF	530-1	6×3mL	EXP	2023-05-27	LOT	T328E61	更新日期：	2024-05
1 / 1								

批号：T328E61					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CK-MB mass)肌酸激酶-MB质量					
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.01	10.9%	<5	<100
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.1	9.4%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.61	6.2%	5	105
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.01	10%	<5	<100
(CRP)C反应蛋白					
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	mg/dL	0.176	9.6%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Roche ACCU-CHEK series Roche \ Nephelometry	mg/L	0.123	10%	<5	<100
(hs-CRP)超敏C反应蛋白					
Beckman AU series Beckman Coulter \ immunoturbidimetric	mg/L	1.42	10%	<5	<100
(Myo / Mb)肌红蛋白					
Roche cobas c 501/502/701/702 DENUO \ immunoturbidimetric	μg/L	71.3	6.9%	8	123
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	143	11.9%	8	115
Roche cobas c 501/502/701/702 Reebio \ OTHER	μg/L	66.7	2.3%	<5	<100
Roche cobas c 501/502/701/702 Reebio \ immunoturbidimetric	μg/L	58.2	7.5%	13	<100
Ortho Vitros 3600/5600 Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	125	4.7%	11	314
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	57.4	10%	<5	<100
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	μg/L	19.2	9%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	97	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	11.4	10%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.021	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。