

REF	500-3 6x3mL	EXP 2025-10-20	LOT 2FB6E52	更新日期: 2025-03 1 / 3
------------	----------------	-----------------------	--------------------	------------------------

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(TT3)总三碘甲状原氨酸	(FT3)游离T3	(TT4)总T4
(FT4)游离T4	(TSH)促甲状腺素	(Cortisol)皮质醇
(FSH)促卵泡刺激素	(LH)促黄体生成素	(PROG)孕酮
(PRL)催/泌乳素	(TES)睾酮	(E2)雌二醇
(C-P)C-肽	(FA)叶酸	(INS)胰岛素
(VB12)维生素B12	(hCG)人绒毛膜促性腺激素	(β-hCG)人绒毛膜促性腺激素β亚基
(25-OH-Vit D)25羟基维生素D		

REF	500-3	6x3mL	EXP	2025-10-20	LOT	2FB6E52	更新日期:	2025-03
								2 / 3

批号: 2FB6E52			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(25-OH-Vit D)25羟基维生素D			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	35	26.6-43.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	87.4	66.4-108.4
(Cortisol)皮质醇			
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	243	227-259
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	254	193-315
(C-P)C-肽			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.97	9.43-10.51
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.37	7.8-8.94
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.11	6.16-10.06
(E2)雌二醇			
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1948	1784-2112
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1299	987-1611
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1347	1255-1439
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	3270	2485-4055
(FA)叶酸			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3	2.28-3.72
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	2.87-4.17
(FSH)促卵泡刺激素			
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	56.8	53.6-60
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	37.8	28.7-46.9
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	42.4	32.2-52.6
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	IU/L	40.2	38.4-42
(FT3)游离T3			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.7	15.3-18.1
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.4	20.1-32.7
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23	21-25
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20	15.2-24.8
(FT4)游离T4			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	49.4	44.2-54.6
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	40.9	31.1-50.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.8	24.9-40.7
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	pmol/L	37.4	35.2-39.6
(hCG)人绒毛膜促性腺激素			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	268	137-399
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	647	585-709
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	514	391-637
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	IU/L	522	505-539
(INS)胰岛素			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	74.9	69.5-80.3
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	100	76-124
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	mIU/L	97.9	89.3-106.5
(LH)促黄体生成素			
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	63.5	57.3-69.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.8	14.3-23.3
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	48.5	36.9-60.1
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	IU/L	44.7	41.9-47.5

REF	500-3	6x3mL	EXP	2025-10-20	LOT	2FB6E52	更新日期:	2025-03
								3 / 3

(PRL)催泌乳素			
Ortho Vitros 3600/5600			
OTHER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	35.1	32.5-37.7
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	67.7	51.5-83.9
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	62	57.8-66.2
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	45.8	34.8-56.8
(PROG)孕酮			
Ortho Vitros 3600/5600			
OTHER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.8	13.2-16.4
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.2	22.3-26.1
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23.1	17.6-28.6
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.9	20.4-33.4
Beckman AU series			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	22.4	20-24.8
(TES)睾酮			
Ortho Vitros 3600/5600			
OTHER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.4	11.4-13.4
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.1	9.2-15
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.5	9.6-11.4
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.86	5.97-9.75
(TSH)促甲状腺素			
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	19.4	16.6-22.2
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	mIU/L	22.1	20.3-23.9
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	23.1	17.6-28.6
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	16.1	12.2-20
(TT3)总三碘甲状腺原氨酸			
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.75	2.29-5.21
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.02	2.3-3.74
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.06	3.09-5.03
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.33	1.73-2.93
(TT4)总T4			
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	59.7	54.8-64.6
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	127	97-157
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	106	81-131
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	128	124-132
(VB12)维生素B12			
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	736	559-913
Lifotronic ecl Series			
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	pg/mL	735	722-748
(β-hCG)人绒毛膜促性腺激素β亚基			
OTHER			
Wondfo \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.07	0.88-1.26
Wondfo shine i2900			
Wondfo \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.17	0.92-1.42

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。