

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期: 2025-03
							1 / 11

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(CK)肌酸激酶	(CK-MB)肌酸激酶-MB (活性)	(Myo / Mb)肌红蛋白
(TnI)肌钙蛋白I	(TnT)肌钙蛋白T	(NT-proBNP)N-末端脑钠肽前体
(DGXN)地高辛	(TT3)总三碘甲状腺氨酸	(FT3)游离T3
(TT4)总T4	(FT4)游离T4	(TSH)促甲状腺素
(Cortisol)皮质醇	(FSH)促卵泡刺激素	(LH)促黄体生成素
(PROG)孕酮	(PRL)催乳素	(TES)睾酮
(E2)雌二醇	(C-P)C-肽	(FA)叶酸
(INS)胰岛素	(VB12)维生素B12	(IgE)免疫球蛋白 IgE
(AFP)甲胎蛋白	(CEA)癌胚抗原	(hCG)人绒毛膜促性腺激素
(PSA)前列腺特异性抗原	(CA 125)糖类抗原 125	(CA 15-3)糖类抗原 15-3
(CA 19-9)糖类抗原 19-9	(β2M)β-2-微球蛋白	(Ferritin)铁蛋白
(β-hCG)人绒毛膜促性腺激素β亚基	(F-PSA)游离前列腺特异性抗原	(HCY)同型半胱氨酸
(17-OHP)17-羟孕酮	(SCC)鳞状细胞癌抗原	(ASD)雄烯二酮
(ALD)醛固酮	(Anti-Tg)抗甲状腺球蛋白抗体	(Anti-TPO)抗甲状腺过氧化物酶抗体
(DHEA-S)硫酸脱氢表雄酮	(EPO)促红细胞生成素	(HGH)人生长激素
(IGF1)胰岛素生长因子1	(PTH)甲状旁腺素	(SHBG)性激素结合球蛋白
(Tg)甲状腺球蛋白	(25-OH-Vit D)25羟基维生素D	(25-OH-VD3)25羟基维生素D3
(PCT)降钙素原	(ACTH)促肾上腺皮质激素	(Renin)肾素
(PINP)总I型胶原氨基端延长肽	(β-CTx)β胶原特殊序列	(TRAb)促甲状腺素受体抗体
(CT)降钙素	(IL6)白介素6	(HE4)人附睾蛋白4
(proGRP)胃泌素释放肽前体	(AMH)抗缪勒氏管激素	(BGP (OCN))骨钙素
(Digitoxin)洋地黄毒苷	(HS-CTnI) (高) 超敏肌钙蛋白I	(CK-MB mass)肌酸激酶-MB质量

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								2 / 11

批号: 3384FDF			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Roche cobas E601/602	ng/mL	1.99	1.73-2.25
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Siemens Immulite1000/2000	ng/mL	1.78	1.41-2.15
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			
(25-OH-VD3)25羟基维生素D3			
Roche cobas E601/602	ng/mL	56.7	54.4-59
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
(25-OH-Vit D)25羟基维生素D			
Roche cobas E601/602	ng/mL	31.5	28.3-34.7
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas E601/602	nmol/L	78.7	70.7-86.7
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas 4000/6000/8000	ng/mL	38.2	31.9-44.5
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas 4000/6000/8000	nmol/L	95.6	79.7-111.5
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	ng/mL	42.9	37.7-48.1
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)			
Abbott ARCHITECT i series	nmol/L	107	94-120
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)			
(ACTH)促肾上腺皮质激素			
Mindray CL-2000i/1000i/6000i series	pg/mL	1	1-1
Mindray \ Chemiluminescence Immunoassay (CLIA)			
(AFP)甲胎蛋白			
Siemens ADVIA Centaur XP/XPT	ng/mL	123	112-134
Dedicated \ Electro chemi luminescence (ECL)			
Abbott ARCHITECT ci series	ng/mL	129	122-136
Dedicated \ Chemiluminescence Immunoassay (CLIA)			
Abbott ARCHITECT i series	ng/mL	119	105-133
Abbott \ Chemiluminescence Immunoassay (CLIA)			
Roche cobas E601/602	ng/mL	56.1	0-157
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	ng/mL	122	117-127
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			
Siemens Atellica	ng/mL	122	100-144
Siemens \ Chemiluminescence Immunoassay (CLIA)			
Beckman UniCel Dxl 600/800	ng/mL	121	107-135
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)			
Abbott ARCHITECT ci series	ng/mL	121	121-121
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)			
Caris200	ng/mL	105	80-130
WanTai \ Chemiluminescence Immunoassay (CLIA)			
Tosoh AIA Series	ng/mL	117	95-139
TOSOH \ Chemiluminescence Immunoassay (CLIA)			
(ALD)醛固酮			
Mindray CL-2000i/1000i/6000i series	pg/mL	174	92-256
Mindray \ Chemiluminescence Immunoassay (CLIA)			
DiaSorin Liaison /XL	pg/mL	167	135-199
Dedicated \ Chemiluminescence Immunoassay (CLIA)			
Autolumo A2000/A2000 PLUS	pg/mL	3164	2531-3797
Autobio \ Chemiluminescence Immunoassay (CLIA)			
(AMH)抗缪勒氏管激素			
Beckman UniCel Dxl 600/800	ng/mL	2.6	2.29-2.91
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			
(Anti-Tg)抗甲状腺球蛋白抗体			
YHLO iFlash Analyzer	IU/mL	147	127-167
YHL \ Chemiluminescence Immunoassay (CLIA)			
Roche cobas E601/602	IU/mL	186	108-264
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	IU/mL	43.5	21.7-65.3
Abbott \ Chemiluminescence Immunoassay (CLIA)			
Roche cobas 4000/6000/8000	IU/mL	79.4	67.6-91.2
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	IU/mL	20.5	16-25
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)			
Siemens IM1600	IU/mL	56.4	51.4-61.4
Siemens \ Chemiluminescence Immunoassay (CLIA)			
Mindray CL-8000i	IU/mL	6.46	5.68-7.24
Mindray \ Chemiluminescence Immunoassay (CLIA)			
Siemens ADVIA Centaur XP/XPT	IU/mL	64.6	51.8-77.4
Siemens \ Chemiluminescence Immunoassay (CLIA)			
(Anti-TPO)抗甲状腺过氧化物酶抗体			
YHLO iFlash Analyzer	IU/mL	29.2	27.3-31.1
YHL \ Chemiluminescence Immunoassay (CLIA)			
Roche cobas E601/602	IU/mL	34.3	24.6-44
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								3 / 11

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	38.2	33.7-42.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	25.3	20.2-30.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	41.3	33.1-49.5
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	47.8	36.3-59.3
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.5	18.9-26.1
(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.12	2.43-3.81
(BGP (OCN))骨钙素			
Roche cobas E601/602 Roche \ Electro chemi luminescence (ECL)	ng/mL	84.7	76.9-92.5
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.6	20.3-58.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	52.6	46.5-58.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	68.3	63.2-73.4
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	57.4	50.9-63.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	54.1	47.5-60.7
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	60.9	60.2-61.6
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.9	52.3-77.5
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	85.7	69.4-102
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	70.8	53.8-87.8
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	24.1	0-62.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	47.6	41.4-53.8
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	59.8	52.9-66.7
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	50.3	43.3-57.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.6	20.5-24.7
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	52.1	51.4-52.8
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.1	45.9-70.3
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	50.6	40.4-60.8
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	85.9	0-267.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	46.8	41.7-51.9
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	309	261-357
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	131	105-157
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	85.7	74.6-96.8
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	270	269-271
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.4	31.6-47.2
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	90.5	68.8-112.2
(CEA)癌胚抗原			
Siemens ADVIA Centaur XP/XPT Dedicated \ Electro chemi luminescence (ECL)	ng/mL	8.97	8.38-9.56
Abbott ARCHITECT ci series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.2	9.3-11.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.45	7.27-9.63
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.3	9.5-15.1
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	10.8	10-11.6
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.27	4.43-14.11
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.47	7.76-9.18

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								4 / 11

Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA) Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.4	10.1-10.7
Abbott ARCHITECT i series Abbott \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	11.5	9.3-13.7
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.3	8.6-14
(CK)肌酸激酶			
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	89.4	71.3-107.5
(CK-MB mass)肌酸激酶-MB质量			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	13.4	11.4-15.4
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.3	9.8-14.8
Beckman UniCel Dx C 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.4	11.5-13.3
(CK-MB)肌酸激酶-MB (活性)			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.1	11.9-16.3
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	13.1	10.4-15.8
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	13.4-20.2
(Cortisol)皮质醇			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	160	136-184
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	245	217-273
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	205	176-234
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	174	149-199
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	225	181-269
(C-P)C-肽			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.46	1.01-1.91
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.95	1.3-2.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.4	1.96-2.84
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.13	1.96-2.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.48	1.3-1.66
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.72	1.44-2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.29	1.67-2.91
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.74	1.34-2.14
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.46	1.18-1.74
(CT)降钙素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	88.8	76-101.6
(DGXN)地高辛			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.1	1.9-2.3
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.47	1.98-2.96
(DHEA-S)硫酸脱氢表雄酮			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4939	3882-5996
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6647	5318-7976
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5711	4535-6887
(Digitoxin)洋地黄毒苷			
Roche cobas 4000/6000/8000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	2.66	2.51-2.81
(E2)雌二醇			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	965	896-1034
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	887	720-1054
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	976	791-1161
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	843	782-904

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								5 / 11

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	975	903-1047
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	946	893-999
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1042	836-1248
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	808	645-971
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1022	818-1226
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1072	922-1222
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	870	661-1079

(EPO)促红细胞生成素

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	44.6	40.6-48.6
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	39.7	31.8-47.6
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	45.9	37.3-54.5

(FA)叶酸

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.06	6.85-9.27
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.53	4.88-8.18
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.75	5.83-7.67
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.8	6.66-8.94
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.37	5.78-8.96

(Ferritin)铁蛋白

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	129	114-144
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	97.3	24.5-170.1
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	208	185-231
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	196	174-218
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	128	113-143
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	137	109-165
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	153	121-185
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	210	160-260
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	189	151-227
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	247	217-277

(F-PSA)游离前列腺特异性抗原

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.79	4.82-6.76
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.96	5.9-8.02
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.41	5.92-6.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.29	6.66-7.92
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.65	6.53-6.77
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.74	4.36-7.12

(FSH)促卵泡刺激素

Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	31.4	27.9-34.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	26.1	22.8-29.4
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.2	29.8-38.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	25	21.6-28.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	25.1	17.7-32.5
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	31.4	25.3-37.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.8	19.3-26.3

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								6 / 11

Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	25.9	24.2-27.6
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	25.6	22.3-28.9
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	38.9	31.3-46.5
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.8	21-32.6
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	23.6	17.9-29.3
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.6	17.7-23.5

(FT3)游离T3

Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.7	11-12.4
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.8	10.7-12.9
Abbott ARCHITECT ci series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.06	8.1-10.02
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.7	10-11.4
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	12.7	10.9-14.5
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.4	6.5-18.3
Mindray BS series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.1	9.9-12.3
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	12.1	10.8-13.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	9.59	8.3-10.88
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	12.2	11.3-13.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.67	7.54-9.8
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	10.5	10.2-10.8
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.36	6.35-10.37
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	9.1-13.5
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.9	10-13.8
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.4	9.4-15.4

(FT4)游离T4

Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.1	25.7-30.5
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	31.8	29.3-34.3
Abbott ARCHITECT ci series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	21.5	19.6-23.4
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17.9	16.7-19.1
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	25.5	21.1-29.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	34.3	30.4-38.2
Mindray BS series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	41.8	38.6-45
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33.5	28.4-38.6
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	24	20.9-27.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33	29.2-36.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	27.9	24.5-31.3
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	25.5	24.3-26.7
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	32.2	25.8-38.6
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	25	19-31
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	34.2	27-41.4
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	42.9	36.9-48.9
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	38.2	29-47.4

(hCG)人绒毛膜促性腺激素

Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	121	113-129
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	110	72-148
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	147	119-175

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								7 / 11

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	152	124-180
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	234	178-290
(HCY)同型半胱氨酸			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	21.8	16.7-26.9
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	21.7	17.4-26
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.1	17.3-26.9
(HE4)人附睾蛋白4			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	39.4	25.1-53.7
(HGH)人生长激素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.34	3.67-5.01
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	3.97	3.17-4.77
(HS-CTnl) (高) 超敏肌钙蛋白I			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.445	0.407-0.483
(IgE)免疫球蛋白 IgE			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.6	64.2-97
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	82.6	66.1-99.1
(IGF1)胰岛素生长因子1			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	237	190-284
(IL6)白介素6			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	245	204-286
(INS)胰岛素			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.2	25.7-30.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	24	19.5-28.5
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	33.6	31.3-35.9
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32	29.4-34.6
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	33.3	30.6-36
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29	26.6-31.4
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.2	24.8-37.6
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.1	25.3-38.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.7	26.2-39.2
(LH)促黄体生成素			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.9	16.5-21.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	15.5	14.1-16.9
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	15.3	11.1-19.5
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.1	15.4-20.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.4	18.5-24.3
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.5	16.7-20.3
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	15	13.6-16.4
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	15.8	14.3-17.3
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	20.6	18.7-22.5
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	17.9	14.1-21.7
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	28.5	22.9-34.1
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.9	18-23.8
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.1	24.4-39.8
(Myo / Mb)肌红蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	124	108-140

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								8 / 11

Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	µg/L	142	13-271
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	53.9	49.7-58.1
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	µg/L	143	115-171
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	139	111-167
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	µg/L	180	155-205
(NT-proBNP)N-末端脑钠肽前体			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1273	1110-1436
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	150	131-169
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	995	913-1077
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	117	107-127
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1184	947-1421
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	139	111-167
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pg/mL	1203	1075-1331
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	141	126-156
(PCT)降钙素原			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.01	8.36-9.66
(PINP)总I型胶原氨基端延长肽			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	153	122-184
(PRL)催/泌乳素			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.3	10.2-12.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.6	9.6-23.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.9	15.4-20.4
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.3	9.4-13.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	17.1	14.5-19.7
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	13.4	12.4-14.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.3	14.6-18
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14	13.1-14.9
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.9	13.1-18.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.2	12.2-18.2
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16	12.2-19.8
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.2	13.1-17.3
(PROG)孕酮			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.44	2.22-2.66
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.25	6.92-11.58
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.27	5.94-10.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.6	9.8-13.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.03	8.02-10.04
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.16	6.87-9.45
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.42	7.95-10.89
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.7	10.1-11.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.21	7.79-10.63
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	9-13.4
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.6	9.1-14.1
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.3	11.4-15.2
Siemens Immulite1000/2000 Siemens Dade Actin FS(APTT) \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.23	6.58-9.88

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								9 / 11

(proGRP)胃泌素释放肽前体			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/L	8758	6446-11070
(PSA)前列腺特异性抗原			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.5	12.5-16.5
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12	9.6-14.4
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	11.4	10.4-12.4
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.5	10.6-12.4
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	15.6	13.8-17.4
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	11.1	10.7-11.5
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.9	12.1-19.7
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.5	8.9-14.1
(PTH)甲状旁腺素			
Beckman UniCel DxI 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	8.39	7.63-9.15
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	8.69	5.65-11.73
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.72	5.93-7.51
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	9.16	8.17-10.15
(Renin)肾素			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	9.86	7.02-12.7
(SCC)鳞状细胞癌抗原			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.67	1.05-2.29
(SHBG)性激素结合球蛋白			
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	78.8	59.9-97.7
(TES)睾酮			
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.091	0.077-0.105
Beckman UniCel DxI 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.694	0.597-0.791
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.37	0.72-2.02
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.306	0.228-0.384
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.27	0-18.65
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.541	0.46-0.622
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.501	0.437-0.565
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.354	0-0.848
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.695	0.594-0.796
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.76	0.679-0.841
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.455	0.346-0.564
(Tg)甲状腺球蛋白			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.4	3.2-35.6
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.7	0-11
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.29	8.57-10.01
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.83	5.27-6.39
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.37	1.91-2.83
(TnI)肌钙蛋白I			
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.37	0.296-0.444
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.47	0.35-0.59
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.005	0.004-0.006

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								10 / 11

(TnT)肌钙蛋白T

Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		0.985					0.806-1.164
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		1					0.8-1.2

(TRAb)促甲状腺素受体抗体

Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L		16.2					14.5-17.9

(TSH)促甲状腺素

Siemens IM1600								
Siemens \ Chemiluminescence Immunoassay (CLIA)	miU/L		9.92					9.01-10.83
maglumi system								
Snibe \ Chemiluminescence Immunoassay (CLIA)	miU/L		10					8.8-11.2
Abbott ARCHITECT ci series								
Dedicated \ Chemiluminescence Immunoassay (CLIA)	miU/L		7.91					7.1-8.72
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA)	miU/L		12.1					11.2-13
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA)	miU/L		8.32					7.5-9.14
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	miU/L		11.2					9.7-12.7
Mindray BS series								
Mindray \ Chemiluminescence Immunoassay (CLIA)	miU/L		11.5					10.5-12.5
Roche cobas E601/602								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	miU/L		11					10.3-11.7
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	miU/L		8.47					7.47-9.47
Roche cobas 4000/6000/8000								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	miU/L		10.9					10.3-11.5
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	miU/L		9.3					8.67-9.93
Abbott ARCHITECT ci series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	miU/L		8.39					8.29-8.49
Siemens Immulite1000/2000								
Siemens \ Chemiluminescence Immunoassay (CLIA)	miU/L		10.2					7.8-12.6
Caris200								
WanTai \ Chemiluminescence Immunoassay (CLIA)	miU/L		10					7.6-12.4
Mindray CL-8000i								
Mindray \ Chemiluminescence Immunoassay (CLIA)	miU/L		8.38					7.21-9.55
Tosoh AIA Series								
TOSOH \ Chemiluminescence Immunoassay (CLIA)	miU/L		9.78					8.06-11.5

(TT3)总三碘甲状原氨酸

Siemens IM1600								
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL		2.41					2.1-2.72
Abbott ARCHITECT ci series								
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL		1.9					1.7-2.1
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL		2.51					2.38-2.64
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL		2.09					1.9-2.28
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL		2.01					1.88-2.14
Roche cobas E601/602								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		2.39					2.19-2.59
Roche cobas 4000/6000/8000								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		2.44					2.26-2.62
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL		2.28					2.01-2.55
Siemens ADVIA Centaur XP/XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL		2.82					2.26-3.38
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		2.55					2.04-3.06
Caris200								
WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL		1.64					1.25-2.03

(TT4)总T4

Siemens IM1600								
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL		146					131-161
Abbott ARCHITECT ci series								
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL		120					108-132
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL		131					124-138
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL		131					114-148
Roche cobas E601/602								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		127					110-144
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL		133					118-148
Roche cobas 4000/6000/8000								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		136					126-146
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL		135					117-153
Caris200								
WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL		103					78-128
Siemens ADVIA Centaur XP/XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL		136					110-162

REF	540-2	6x3mL	EXP	2025-11-16	LOT	3384FDF	更新日期:	2025-03
								11 / 11

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	139	111-167
(VB12)维生素B12			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	509	456-562
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	778	703-853
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	983	904-1062
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	491	424-558
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	659	523-795
(β2M)β-2-微球蛋白			
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	ng/mL	1240	1064-1416
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1440	1120-1760
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	1180	944-1416
(β-CTx)β胶原特殊序列			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	632	565-699
(β-hCG)人绒毛膜促性腺激素β亚基			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	0.126	0.116-0.136
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.116	0.099-0.133
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.125	0.113-0.137
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.136	0.109-0.163
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.24	0.206-0.274
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.166	0.126-0.206

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。