

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期: 2025-03
							1 / 11

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(CK)肌酸激酶	(CK-MB)肌酸激酶-MB (活性)	(Myo / Mb)肌红蛋白
(TnI)肌钙蛋白I	(TnT)肌钙蛋白T	(NT-proBNP)N-末端脑钠肽前体
(DGXN)地高辛	(THEO)茶碱	(TT3)总三碘甲状腺原氨酸
(FT3)游离T3	(TT4)总T4	(FT4 )游离T4
(TSH)促甲状腺素	(Cortisol)皮质醇	(FSH)促卵泡刺激素
(LH)促黄体生成素	(PROG)孕酮	(PRL)催/泌乳素
(TES)睾酮	(E2)雌二醇	(C-P)C-肽
(FA)叶酸	(INS)胰岛素	(VB12)维生素B12
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(CA 125)糖类抗原 125
(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9	(β2M)β-2-微球蛋白
(Ferritin)铁蛋白	(β-hCG)人绒毛膜促性腺激素β亚基	(F-PSA)游离前列腺特异性抗原
(HCY)同型半胱氨酸	(17-OHP)17-羟孕酮	(NSE)神经元特异性烯醇化酶
(Cyfra 21-1)细胞角蛋白19片段	(SCC)鳞状细胞癌抗原	(CA 72-4)糖类抗原 72-4
(ASD)雄烯二酮	(ALD)醛固酮	(Anti-Tg)抗甲状腺球蛋白抗体
(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮	(EPO)促红细胞生成素
(HGH)人生长激素	(IGF1)胰岛素生长因子1	(PTH)甲状旁腺素
(SHBG)性激素结合球蛋白	(Tg)甲状腺球蛋白	(25-OH-Vit D)25羟基维生素D
(25-OH-VD3)25羟基维生素D3	(PCT)降钙素原	(Renin)肾素
(CT)降钙素	(proGRP)胃泌素释放肽前体	(AMH)抗缪勒氏管激素
(HS-CTnI) (高) 超敏肌钙蛋白I	(HS-CTnT)超敏肌钙蛋白T	(CK-MB mass)肌酸激酶-MB质量

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								2 / 11

批号: 3384F73			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
<b>(17-OHP)17-羟孕酮</b>			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.78	1.38-2.18
<b>(25-OH-VD3)25羟基维生素D3</b>			
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.5	15.2-23.8
<b>(25-OH-Vit D)25羟基维生素D</b>			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	34.2	26.6-41.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	85.2	66.3-104.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	30.3	25.7-34.9
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	75.7	64.2-87.2
<b>(AFP)甲胎蛋白</b>			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	121	99-143
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	160	144-176
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	133	104-162
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	134	114-154
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	123	120-126
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	117	92-142
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	106	83-129
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	122	122-122
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	131	123-139
<b>(ALD)醛固酮</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pg/mL	211	174-248
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA )	pg/mL	172	135-209
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	pg/mL	2940	2358-3522
<b>(AMH)抗缪勒氏管激素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.73	2.38-3.08
<b>(Anti-Tg)抗甲状腺球蛋白抗体</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	64.7	53.7-75.7
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	71.9	61.4-82.4
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	23.9	23.3-24.5
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	IU/mL	6.41	5.44-7.38
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	64.6	50.4-78.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	62.4	45.7-79.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	67.2	59.1-75.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	6.94	4.09-9.79
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	58.9	47-70.8
<b>(Anti-TPO)抗甲状腺过氧化物酶抗体</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	51.2	34-68.4
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	IU/mL	24.5	20.3-28.7
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	38.1	37.4-38.8
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	38.8	31.4-46.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	24.1	19.7-28.5
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	27.9	9.5-46.3
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	35.8	35.8-35.8

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								3 / 11

(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.12	2.41-3.83
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	65.5	51.5-79.5
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	55.8	49.9-61.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	52	41.2-62.8
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	69.7	68.6-70.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	61.4	47.6-75.2
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	IU/mL	85.7	67.5-103.9
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	64.9	51.1-78.7
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	73.4	71.8-75
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	52.1	51.7-52.5
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	59.8	57.5-62.1
(CA 15-3)糖类抗原 15-3			
Siemens ADVIA 1800/2400 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	60.3	50.3-70.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	51.5	40.9-62.1
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	68	52.5-83.5
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	71.8	70.5-73.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	45.3	35.5-55.1
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	50.6	39.3-61.9
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	123	33-213
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	IU/mL	58.1	45.3-70.9
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	78.1	71.5-84.7
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	52	46.5-57.5
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	56.9	56.9-56.9
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	321	263-379
Siemens ADVIA 1800/2400 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	122	102-142
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	51.6	40.7-62.5
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	351	348-354
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	IU/mL	39.4	30.7-48.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	46.8	36.3-57.3
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	360	358-362
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	55.9	54.4-57.4
(CA 72-4)糖类抗原 72-4			
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	2.26	2.23-2.29
(CEA)癌胚抗原			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.1	9.8-12.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	12.6	12-13.2
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.8	8.3-13.3
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	9.4-11.4
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	14.3	13.8-14.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	10.5	8.1-12.9
Abbott ARCHITECT i series Abbott \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	11.5	9-14
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	13.6	12.7-14.5
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	10.8	10.1-11.5

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								4 / 11

<b>(CK)肌酸激酶</b>			
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	57.6	44.8-70.4
<b>(CK-MB mass)肌酸激酶-MB质量</b>			
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.91	6.71-9.11
<b>(CK-MB)肌酸激酶-MB (活性)</b>			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.2	7.9-12.5
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	16.8	13-20.6
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.72	7.46-11.98
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.8	9.5-12.1
<b>(Cortisol)皮质醇</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	192	134-250
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	161	126-196
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	234	178-290
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	225	178-272
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	219	195-243
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	184	166-202
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	180	163-197
<b>(C-P)C-肽</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.63	1.42-1.84
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.18	2.1-2.26
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.22	1.74-2.7
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.23	1.95-2.51
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.35	1.84-2.86
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.74	1.33-2.15
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.46	1.15-1.77
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.28	1.81-2.75
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.87	1.81-1.93
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.25	2.15-2.35
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.28	2.19-2.37
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.82	2.72-2.92
<b>(CT)降钙素</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	0.706	0.641-0.771
<b>(Cyfra 21-1)细胞角蛋白19片段</b>			
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.92	1.81-2.03
<b>(DGXN)地高辛</b>			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.48	1.92-3.04
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.97	1.5-2.44
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.01	1.6-2.42
<b>(DHEA-S)硫酸脱氢表雄酮</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5849	3556-8142
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5711	4409-7013
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6006	4685-7327
<b>(E2)雌二醇</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	933	825-1041
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	978	790-1166

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								5 / 11

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	1042	817-1267
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	pmol/L	808	629-987
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	985	952-1018
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	807	797-817
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1048	947-1149
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	698	549-847
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	959	819-1099
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	pmol/L	1010	889-1131
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	936	885-987
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	966	966-966
<b>(EPO)促红细胞生成素</b>			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA )	mIU/mL	55.5	42.8-68.2
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/mL	45.9	36.1-55.7
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	44.8	38.2-51.4
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	mIU/mL	45	40.9-49.1
<b>(FA)叶酸</b>			
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.18	7.02-9.34
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.32	7.29-9.35
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.73	6.04-9.42
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.37	5.67-9.07
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.83	4.54-7.12
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.23	7.47-8.99
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.32	4.82-7.82
<b>(Ferritin)铁蛋白</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	207	167-247
Siemens ADVIA 1800/2400 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	171	156-186
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	162	147-177
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	153	119-187
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	137	106-168
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	129	101-157
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	178	138-218
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	130	113-147
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	130	114-146
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	204	198-210
<b>(F-PSA)游离前列腺特异性抗原</b>			
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.95	6.02-7.88
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.51	4.18-6.84
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.34	5.3-5.38
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.67	6.51-6.83
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.75	3.96-7.54
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.92	5.36-8.48
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.21	4.96-5.46

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								6 / 11

**(FSH)促卵泡刺激素**

Siemens ADVIA 1800/2400	IU/L	43.6	40.8-46.4
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	IU/L	31.5	23.6-39.4
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	IU/L	34.8	27.5-42.1
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Siemens Immulite1000/2000	IU/L	38.9	30.6-47.2
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Tosoh AIA Series	IU/L	26.8	20.7-32.9
TOSOH \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	IU/L	33.8	32.5-35.1
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	IU/L	27.2	24.1-30.3
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			
Ortho Vitros 3600/5600	IU/L	20.3	17.5-23.1
OTHER \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/L	28.5	22.7-34.3
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/L	27.7	23.2-32.2
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e 801	IU/L	34.1	32.7-35.5
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	IU/L	23.5	22.6-24.4
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(FT3)游离T3**

Ortho Vitros 3600/5600	pmol/L	16.6	15.1-18.1
OTHER \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	pmol/L	9.38	8.07-10.69
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT ci series	pmol/L	11.5	11.3-11.7
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT i series	pmol/L	12.8	7.3-18.3
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Siemens ADVIA Centaur XP/XPT	pmol/L	11.3	8.9-13.7
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	pmol/L	15.1	14-16.2
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott Alinity i Series	pmol/L	11.2	10.1-12.3
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	pmol/L	14.6	13.9-15.3
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas E601/602	pmol/L	14.1	11.9-16.3
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	pmol/L	9.46	7.38-11.54
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	pmol/L	14.8	13.8-15.8
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e 801	pmol/L	14.8	14.1-15.5
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	pmol/L	12.6	12.2-13
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(FT4 )游离T4**

Siemens ADVIA 1800/2400	pmol/L	27.2	24.9-29.5
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Ortho Vitros 3600/5600	pmol/L	37	34.7-39.3
OTHER \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	pmol/L	27.4	24.1-30.7
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT ci series	pmol/L	22.8	21.9-23.7
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT i series	pmol/L	26.5	20.8-32.2
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Siemens ADVIA Centaur XP/XPT	pmol/L	32.2	25.1-39.3
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)			
Tosoh AIA Series	pmol/L	34.2	26.5-41.9
TOSOH \ Chemiluminescence Immunoassay (CLIA )			
Abbott Alinity i Series	pmol/L	24.3	24.1-24.5
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	pmol/L	34.6	32.2-37
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	pmol/L	27.7	23.4-32
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	pmol/L	31.2	23.1-39.3
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas 4000/6000/8000	pmol/L	36.2	27.5-44.9
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e 801	pmol/L	33.1	33.1-33.1
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	pmol/L	37.4	33.7-41.1
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e411	pmol/L	33	33-33
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(hCG)人绒毛膜促性腺激素**

Siemens ADVIA 1800/2400	IU/L	147	131-163
Siemens \ Chemiluminescence Immunoassay (CLIA )			

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								7 / 11

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	147	115-179
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	152	120-184
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	IU/L	212	191-233
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	172	161-183
<b>(HCY)同型半胱氨酸</b>			
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	24.8	18.8-30.8
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	μmol/L	24.3	19-29.6
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	μmol/L	22.1	17-27.2
TOSHIBA TBA Bio-assay \ Enzymatic Cycling Assays Method	μmol/L	20.6	19.9-21.3
<b>(HGH)人生长激素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	μg/L	3.53	2.78-4.28
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	μg/L	3.53	2.76-4.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.59	3.6-5.58
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	μg/L	3.97	3.1-4.84
<b>(HS-CTnI) (高) 超敏肌钙蛋白I</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.447	0.4-0.494
<b>(HS-CTnT)超敏肌钙蛋白T</b>			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.916	0.916-0.916
<b>(IgE)免疫球蛋白 IgE</b>			
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	84.9	67.1-102.7
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	IU/mL	80.6	62.7-98.5
<b>(IGF1)胰岛素生长因子1</b>			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	211	165-257
Siemens Immulite 1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA )	nmol/L	22.9	17.8-28
<b>(INS)胰岛素</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	29.5	26.1-32.9
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	35.5	29.2-41.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.3	24.5-38.1
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	mIU/L	28.7	26.1-31.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	34.7	27-42.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	25.1	19.6-30.6
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	31.2	24.5-37.9
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	mIU/L	32.1	24.9-39.3
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	22.8	21.2-24.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	34.6	32-37.2
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	32.4	31.1-33.7
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.8	31.3-34.3
<b>(LH)促黄体生成素</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	20.1	18.3-21.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.3	18.2-26.4
Siemens Immulite 1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	28.5	22.3-34.7
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	IU/L	17.9	13.8-22
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	19.3	15.1-23.5
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.2	21.8-22.6
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	16.4	15.7-17.1
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	IU/L	19	16.8-21.2

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03
								8 / 11

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	16.7	14.6-18.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	17.3	13.3-21.3
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	IU/L	16.7	15.1-18.3
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	22	21.6-22.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.2	20.4-22
<b>(Myo / Mb)肌红蛋白</b>			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	149	117-181
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	µg/L	143	113-173
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	µg/L	126	99-153
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	µg/L	127	115-139
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	131	108-154
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA )	µg/L	196	173-219
TOSHIBA TBA OTHER \ OTHER	µg/L	88.2	88.2-88.2
<b>(NSE)神经元特异性烯醇化酶</b>			
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.24	1.23-1.25
<b>(NT-proBNP)N-末端脑钠肽前体</b>			
J&J Vitros 250/350 Dedicated \ Chemiluminescence Immunoassay (CLIA )	pg/mL	2146	2077-2215
J&J Vitros 250/350 Dedicated \ Chemiluminescence Immunoassay (CLIA )	pmol/L	252	244-260
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	pg/mL	151	144-158
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	17.8	17-18.6
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1172	1132-1212
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	138	133-143
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	520	0-1632
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	61.1	0-191.7
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA )	pg/mL	2230	2034-2426
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA )	pmol/L	262	239-285
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1310	1244-1376
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	154	146-162
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1164	1117-1211
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	137	132-142
<b>(PCT)降钙素原</b>			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.28	8.05-8.51
<b>(PRL)催/泌乳素</b>			
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	12.4	11.7-13.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.4	11.9-20.9
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	15.9	12.5-19.3
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	15.2	11.8-18.6
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.9	19.1-20.7
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	14.4	12.2-16.6
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	13.2	11.4-15
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	15.1	12.2-18
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	14.2	12.8-15.6
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	18.3	14.6-22
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.8	16.4-17.2



<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03	9 / 11
------------	-------	-------	------------	------------	------------	---------	-------	---------	--------

(PROG)孕酮				
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	10.1		7.6-12.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.6		8.3-12.9
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.07		0.062-0.078
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.3		9.5-11.1
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.6		9-14.2
Siemens Immulite1000/2000 Siemens Dade Actin FS(APTT) \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.23		6.47-9.99
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.2		8.7-13.7
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.69		7.6-9.78
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.7		6.51-8.89
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.87		8.39-11.35
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	9.19		7.87-10.51
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.2		9.4-11
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	10.3		9.7-10.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.5		10.1-10.9
(proGRP)胃泌素释放肽前体				
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA )	ng/L	22000		21736-22264
(PSA)前列腺特异性抗原				
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	13.6		12.2-15
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.5		10.7-16.3
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.2		10.8-11.6
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.5		8.8-14.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	12.5		9.6-15.4
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	12.7		10-15.4
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.8		11.6-12
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.3		11.7-12.9
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	16.8		16.8-16.8
(PTH)甲状旁腺素				
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	10.8		10.2-11.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	8.46		6.6-10.32
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	8.77		6.82-10.72
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	10.6		9.3-11.9
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	pmol/L	9.49		7.21-11.77
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.33		6.04-6.62
(Renin)肾素				
DiaSorin Liaison /XL Diasorin \ Electro-Chemiluminescence Immunoassay (ECLIA)	μIU/mL	43.5		19.8-67.2
(SCC)鳞状细胞癌抗原				
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.18		1.11-1.25
(SHBG)性激素结合球蛋白				
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	80		68.6-91.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	nmol/L	78.8		60-97.6
(TES)睾酮				
Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.341		0.291-0.391
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.448		0.323-0.573
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.42		0.328-0.512

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03	10 / 11
------------	-------	-------	------------	------------	------------	---------	-------	---------	---------

Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.59	0.458-0.722
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.692	0.554-0.83
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.622	0.514-0.73
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.672	0.535-0.809
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.63	0.459-0.801
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.482	0.378-0.586
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.515	0.493-0.537
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.45	0.45-0.45

**(Tg)甲状腺球蛋白**

Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12	11-13
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.13	7.35-8.91
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.2	9.5-14.9
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.37	1.85-2.89
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.1	10.4-11.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.32	6.66-9.98

**(THEO)茶碱**

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	24.1	18.9-29.3
---	-------	------	-----------

**(Tni)肌钙蛋白I**

Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.668	0.635-0.701
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.42	0.309-0.531
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.47	0.353-0.587
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.345	0.287-0.403

**(TnT)肌钙蛋白T**

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.03	0.8-1.26
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.18	0.89-1.47
Beckman UniCel Dxl 600/800 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.912	0.821-1.003

**(TSH)促甲状腺素**

Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	9.67	8.88-10.46
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	mIU/L	8.77	6.95-10.59
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA )	mIU/L	10.2	9.4-11
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	8.77	8.61-8.93
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA )	mIU/L	9.78	7.73-11.83
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.7	8.2-13.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	8.76	6.82-10.7
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	8.26	7.58-8.94
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	10.5-11.3
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.4	10.9-11.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.7	9.8-13.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	10.1	8.9-11.3
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11	9.8-12.2
Roche cobas e 801 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	10.9	10.6-11.2

**(TT3)总三碘甲状腺原氨酸**

Siemens ADVIA 1800/2400 Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.004	0.004-0.004
Ortho Vitros 3600/5600 OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.34	2.89-3.79

<b>REF</b>	540-2	6x3mL	<b>EXP</b>	2024-04-02	<b>LOT</b>	3384F73	更新日期:	2025-03	11 / 11
------------	-------	-------	------------	------------	------------	---------	-------	---------	---------

Roche cobas 4000/6000/8000	ng/mL	2.86	2.61-3.11
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	ng/mL	2.1	1.83-2.37
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Siemens ADVIA Centaur XP/XPT	ng/mL	2.82	2.2-3.44
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT i series	ng/mL	2.15	1.68-2.62
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Abbott Alinity i Series	ng/mL	3.16	2.95-3.37
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	ng/mL	2.31	2.18-2.44
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas 4000/6000/8000	ng/mL	2.65	2.45-2.85
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	ng/mL	2.14	1.61-2.67
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	ng/mL	2.5	2.08-2.92
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e411	ng/mL	2.78	2.47-3.09
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e 801	ng/mL	2.82	2.68-2.96
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)			

**(TT4)总T4**

Siemens ADVIA 1800/2400	ng/mL	135	124-146
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Ortho Vitros 3600/5600	ng/mL	147	114-180
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas 4000/6000/8000	ng/mL	128	119-137
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	ng/mL	105	92-118
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	ng/mL	124	114-134
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Siemens ADVIA Centaur XP/XPT	ng/mL	136	106-166
Siemens \ Chemiluminescence Immunoassay (CLIA )			
Abbott ARCHITECT i series	ng/mL	128	95-161
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Abbott Alinity i Series	ng/mL	125	122-128
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	ng/mL	107	88-126
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	ng/mL	114	96-132
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Tosoh AIA Series	ng/mL	0.009	0.007-0.011
TOSOH \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	ng/mL	125	115-135
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e 801	ng/mL	110	99-121
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	ng/mL	144	140-148
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(VB12)维生素B12**

Beckman UniCel Dxl 600/800	pg/mL	511	391-631
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
YHLO iFlash Analyzer	pg/mL	486	428-544
YHL \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	pg/mL	521	451-591
Dedicated \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	pg/mL	680	517-843
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Tosoh AIA Series	pg/mL	659	517-801
TOSOH \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	pg/mL	507	440-574
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			

**(β2M)β-2-微球蛋白**

Tosoh AIA Series	ng/mL	1440	1106-1774
TOSOH \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	ng/mL	1390	1081-1699
Roche \ Immunoturbidimetry (ITA)			
Beckman AU series	ng/mL	1434	1319-1549
Beckman Coulter \ Immunoturbidimetry (ITA)			

**(β-hCG)人绒毛膜促性腺激素β亚基**

Beckman UniCel Dxl 600/800	IU/mL	0.148	0.107-0.189
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	IU/mL	0.176	0.166-0.186
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	IU/mL	0.168	0.155-0.181
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	IU/mL	0.156	0.126-0.186
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Beckman UniCel Dxl 600/800	IU/mL	0.173	0.163-0.183
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Ortho Vitros 3600/5600	IU/mL	0.215	0.182-0.248
Dedicated \ Chemiluminescence Immunoassay (CLIA )			

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。