

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
1 / 9								

		批号：3384F15			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.8	10.3%	<5	<100
(25-OH-Vit D)25羟基维生素D					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.8	12.3%	9	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	59.5	12.3%	9	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	34.1	10.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	85.1	10.5%	<5	<100
(AFP)甲胎蛋白					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	126	11.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	137	7.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	105	7.1%	31	859
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	7.1%	<5	<100
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	4%	25	700
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	4.1%	11	107
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	105	10.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	9.5%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	116	5%	<5	<100
(ALD)醛固酮					
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pg/mL	181	10.6%	10	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	179	9.3%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	2900	10%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	58.9	4.9%	9	217
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	52.5	12.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	6.31	11.4%	30	757
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.3	10.1%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	75.9	4.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	11.5	7.7%	24	650
Abbott ARCHITECT ci series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	19.2	7%	11	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	27.7	4.7%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	38.1	10.6%	5	116
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	22.6	9.7%	31	797
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	26.7	8.1%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.2	8.4%	24	645
Abbott ARCHITECT ci series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.6	2.9%	11	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	20.9	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.7	10.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	25.4	6.3%	<5	<100
(ASD)雄烯二酮					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.92	11.8%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
2 / 9								

批号：3384F15					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CA 125)糖类抗原 125					
Roche cobas E601/602	IU/mL	51.9	8.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens ADVIA Centaur XP/XPT	IU/mL	62.8	3.5%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	56.7	4.7%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	IU/mL	50.6	10.3%	20	559
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	IU/mL	68.2	6.2%	11	107
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	IU/mL	61.3	1.9%	10	204
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	86.6	9.4%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	58.1	3%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	IU/mL	62.1	10.1%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
OTHER	IU/mL	52	0%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
(CA 15-3)糖类抗原 15-3					
Abbott ARCHITECT i series	IU/mL	50.4	9.8%	8	159
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	55.6	11.4%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens ADVIA Centaur XP/XPT	IU/mL	51.4	5.9%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	26.4	6%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	IU/mL	65.1	15.2%	23	539
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	IU/mL	90.5	4.4%	13	120
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	58.7	10.4%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	45.8	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	47.7	8.4%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
OTHER	IU/mL	27.5	0%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
(CA 19-9)糖类抗原 19-9					
Roche cobas E601/602	IU/mL	55.6	9.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	IU/mL	46.5	11%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	39.8	9.8%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	46.2	4.8%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
OTHER	IU/mL	76.4	0%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
(CEA)癌胚抗原					
Roche cobas E601/602	ng/mL	10.1	12.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens ADVIA Centaur XP/XPT	ng/mL	10.4	7.2%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	10.6	6.6%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	ng/mL	12.6	4.8%	22	586
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	ng/mL	12.5	4.8%	11	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	10.7	11.2%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	9.62	3.3%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(CK)肌酸激酶					
Roche cobas 4000/6000/8000	U/L	57.1	11.9%	<5	<100
Roche \ Urease-UV-rate					
(CK-MB)肌酸激酶-MB (活性)					
Beckman UniCel Dxl 600/800	ng/mL	12.9	5.2%	7	127
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	11.9	8.5%	31	863
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman AU series	ng/mL	11.6	8.4%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	ng/mL	16.2	10.5%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
3 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.78	11.2%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	11.5%	<5	<100
(Cortisol)皮质醇					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	190	8.9%	31	869
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	191	6.4%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	199	5.5%	15	210
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	175	8.1%	12	336
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	10.3%	10	102
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	233	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	207	10.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	231	10.8%	<5	<100
(C-P)C-肽					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.98	2.6%	9	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.55	10.6%	6	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.2	3%	5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.5	7.3%	<5	<100
Abbott ARCHITECT ci series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.36	9.3%	13	105
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2	13.9%	10	208
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.23	11.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.33	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.76	11.4%	<5	<100
(DGXN)地高辛					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.42	11.6%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.09	10.4%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5785	10.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6411	10.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6080	9.6%	<5	<100
(E2)雌二醇					
Abbott ARCHITECT c series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1388	3.1%	8	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	778	3.4%	6	101
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1044	6.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1022	7.8%	32	809
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1072	11.3%	15	378
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1024	11%	11	129
Abbott ARCHITECT i series Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1062	4.8%	10	205
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	812	10.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	998	9.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	967	0.6%	<5	<100
(EPO)促红细胞生成素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	44.8	8.1%	29	750

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
4 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	46.1	5.2%	19	406
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	56.1	10.3%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	45.2	2.4%	<5	<100
(FA)叶酸					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.45	8%	6	103
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.22	6.4%	32	733
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.98	5.9%	31	837
Beckman AU series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.69	4.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.46	10.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.58	0.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.92	11.9%	<5	<100
(Ferritin)铁蛋白					
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	197	4.7%	9	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	194	7.1%	7	125
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	128	6.8%	31	833
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	135	5.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	176	7.1%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	155	7.1%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	129	4.7%	20	470
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	180	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	139	10.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	187	0.2%	<5	<100
OTHER Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	141	0%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.04	8.6%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.23	7.3%	31	863
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.47	5.4%	22	565
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.86	5.6%	10	204
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	10.4%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.84	3.3%	<5	<100
(FSH)促卵泡刺激素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	29	3.4%	6	101
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	27.8	8.2%	31	778
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	36.9	5.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.9	7.9%	15	380
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	30.4	9.9%	11	126
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.9	4.7%	10	201
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	31.8	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	27.1	10.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.9	1.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	40.1	1.2%	<5	<100

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
5 / 9								

		批号：3384F15			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(FT3)游离T3					
Roche cobas E601/602	pmol/L	14.1	8.6%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	pmol/L	8.96	7.2%	33	1577
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Siemens ADVIA Centaur XP/XPT	pmol/L	11.8	4.4%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	pmol/L	9.15	10.9%	27	698
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	pmol/L	11	6%	21	545
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	10.6	9.7%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	pmol/L	13.4	4.7%	15	199
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	pmol/L	14.5	4%	13	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	pmol/L	13.2	1.9%	10	232
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	pmol/L	11.7	4.1%	10	202
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	pmol/L	11.5	10.4%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	9.83	2.7%	<5	<100
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)					
(FT4)游离T4					
Roche cobas E601/602	pmol/L	36.1	9.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens ADVIA Centaur XP/XPT	pmol/L	31.6	5.6%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	pmol/L	26.9	6%	31	1639
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Abbott ARCHITECT i series	pmol/L	23.6	2.8%	<5	<100
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman UniCel Dxl 600/800	pmol/L	29.4	6.9%	27	699
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	26.8	9.5%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT ci series	pmol/L	21.4	4.2%	19	549
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	pmol/L	33.7	7.2%	15	198
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	pmol/L	34.6	6.6%	11	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	pmol/L	32.3	3.6%	10	233
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	pmol/L	22.9	4.3%	10	205
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	pmol/L	29.8	10.1%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	pmol/L	34.6	10.4%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(HCY)同型半胱氨酸					
Abbott ARCHITECT i series	μmol/L	22.9	10%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	μmol/L	22.4	10.7%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	μmol/L	25.1	12%	<5	<100
Roche \ Enzymatic Cycling Assays Method					
(HGH)人生长激素					
Beckman UniCel Dxl 600/800	μg/L	3.32	9.9%	31	825
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Tosoh AIA 900/1800	μg/L	3.8	10.5%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	μg/L	4.43	10.2%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(IgE)免疫球蛋白 IgE					
Beckman Image 800	IU/mL	85.8	8.9%	<5	<100
Dedicated \ Immunoturbidimetry (ITA)					
Tosoh AIA 900/1800	IU/mL	80.9	10.1%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS	ng/mL	207	11.6%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Siemens Immulite1000/2000	nmol/L	22.8	10.1%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
(INS)胰岛素					
Roche cobas 4000/6000/8000	mIU/L	34.3	4.4%	9	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
6 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.1	2.9%	9	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	36.9	1.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	24.4	4.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	30.8	4.8%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	25.1	2.5%	10	201
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.1	10.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.4	11.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.5	10.5%	<5	<100
(LH)促黄体生成素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	24.1	4.7%	6	101
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.3	11.3%	5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	16.7	6.3%	31	824
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	27.2	8.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	14.2	9.2%	15	381
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.7	3.8%	10	193
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.1	10.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	15.6	0.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.6	8.8%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.7	1.5%	<5	<100
(Myo / Mb)肌红蛋白					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	124	7%	31	881
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	117	6.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	120	5.8%	19	461
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	139	10.1%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	145	12%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	0.76	11.8%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	0.089	11.8%	<5	<100
(PRL)催/ 泌乳素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	21.4	5.3%	6	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	13.3	6.1%	31	842
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	15	4.1%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.1	5.7%	15	381
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.1	14.4%	11	135
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.9	9.2%	10	208
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.4	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.2	9.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.3	0.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.4	0.9%	<5	<100
OTHER Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.3	0%	<5	<100
(PROG)孕酮					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.5	4%	6	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
7 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	9.6%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.91	9.6%	32	828
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.38	5.6%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.1	8.8%	15	381
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.3	11%	12	144
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.78	3.6%	10	200
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.7	10.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	0.2%	<5	<100
(PSA)前列腺特异性抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.3	8.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.6	5.8%	31	854
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.3	3.4%	<5	<100
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.1	4.7%	21	552
Abbott ARCHITECT ci series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.8	2.7%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.6	4.5%	10	206
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	11.6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.8	3.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.7	11%	<5	<100
(PTH)甲状旁腺素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	9.27	4%	8	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	8.53	2.9%	6	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	9.07	8.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	10.2	7.8%	31	828
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.59	10.4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.56	10.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.29	1%	<5	<100
(TES)睾酮					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.475	7.2%	6	101
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.575	14.8%	5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.709	9%	46	1261
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.556	14.3%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.408	6.3%	10	205
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.487	10.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.476	11.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.516	1.3%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.366	3%	<5	<100
OTHER Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.764	0%	<5	<100
(Tg)甲状腺球蛋白					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.49	12.7%	9	228
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.72	11.4%	31	798
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.38	7.7%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
8 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.88	8.5%	24	649
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.4	18.4%	10	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.4	9.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	10.4%	<5	<100
(THEO)茶碱					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	24.7	9.5%	<5	<100
(Tnl)肌钙蛋白I					
Beckman AU series Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.426	4.5%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.393	5.7%	24	665
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.46	13%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.05	9.7%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1201	12%	<5	<100
(TSH)促甲状腺素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.3	16.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	9.55	6%	31	1653
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	8.63	6.7%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.67	4.3%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.2	5.2%	27	698
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.85	9.8%	<5	<100
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.61	6.1%	19	496
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.1	5.3%	15	213
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11	5.4%	11	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11	1.9%	10	251
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.21	5.3%	10	207
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.84	8.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	11%	<5	<100
OTHER Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.4	0%	<5	<100
(TT3)总三碘甲状原氨酸					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.5	9.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.02	9.9%	31	839
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.89	5.8%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.75	3%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.28	7.4%	27	666
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.22	7.8%	20	561
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.77	7.7%	15	192
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.74	3.9%	10	240
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.19	1.9%	10	197
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.01	8.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.18	10.1%	<5	<100
(TT4)总T4					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	121	12.4%	<5	<100

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2024-04
9 / 9								

		批号：3384F15			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	139	5.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	114	8%	31	836
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	124	13.3%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	129	4.1%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	6.9%	27	678
Abbott ARCHITECT ci series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	123	5.7%	21	550
Abbott ARCHITECT ci series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	213	1.3%	10	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	119	3.1%	10	242
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	5.6%	10	206
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	132	9.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.009	10.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	107	9.4%	<5	<100
OTHER Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	0%	<5	<100
(VB12)维生素B12					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	506	7.3%	31	842
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	535	7.2%	31	742
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	538	6.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	688	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	665	10.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	672	0%	<5	<100
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	1370	10.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1460	11%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.154	9.8%	8	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.145	10.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.121	5.8%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。