

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
1 / 8								

批号：3385042					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.78	10.4%	<5	<100
(25-OH-VD3)25羟基维生素D3					
Beckman AU series Bio-assay \ Chemiluminescence Immunoassay (CLIA)	ng/mL	29.3	0.3%	<5	<100
(25-OH-Vit D)25羟基维生素D					
Beckman AU series Bio-assay \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	3.8%	<5	<100
Beckman AU series Bio-assay \ Chemiluminescence Immunoassay (CLIA)	nmol/L	281	3.8%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23	5.5%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	nmol/L	57.5	5.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	30.9	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	77.1	10%	<5	<100
(AFP)甲胎蛋白					
YHLO iFlash Analyzer YHL \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	139	3.5%	5	102
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	135	9.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	106	10.4%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	140	2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	9.4%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	9.9%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	101	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	105	12%	<5	<100
Mindray BS series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	132	0%	<5	<100
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	135	0%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	138	0%	<5	<100
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	1411	104.5%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	167	9.6%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	72.7	8.3%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.6	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	9.86	9.9%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	35.7	10.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.4	10.3%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.4	10%	<5	<100
(ASD)雄烯二酮					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.12	11%	<5	<100
(CA 125)糖类抗原 125					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	48.6	9.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	85.7	9.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	61.4	10.3%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	70.8	12%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.9	9.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	62.3	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	53.1	10%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
2 / 8								

批号：3385042					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CA 15-3)糖类抗原 15-3					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50	9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.1	10.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	71.2	10.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.1	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	50.6	10.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.3	10.2%	<5	<100
(CA 19-9)糖类抗原 19-9					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	46.1	6.8%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	89.5	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	90.5	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	319	9.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	46.8	10.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.4	9.9%	<5	<100
(CEA)癌胚抗原					
YHLO iFlash Analyzer YHL \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.45	3.9%	5	103
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.53	7.8%	<5	<100
Abbott ARCHITECT i series Abbott \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	11.5	9.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.5	11.4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.48	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.69	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	1.7%	<5	<100
Mindray BS series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7	0%	<5	<100
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.55	0%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.1	0%	<5	<100
(CK)肌酸激酶					
TOSHIBA TBA TOSHIBA \ Rate -AMP buffer	U/L	88.8	1.2%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	89.4	10.1%	<5	<100
(CK-MB mass)肌酸激酶-MB质量					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.3	10%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
TOSHIBA TBA TOSHIBA \ Immune suppression	U/L	11.9	7.3%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	13.1	10.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	10.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.2	4.6%	<5	<100
(Cortisol)皮质醇					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	192	2.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	216	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	228	11%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	225	9.7%	<5	<100
(C-P)C-肽					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.3	3.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.16	6.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.35	9.8%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
3 / 8								

		批号：3385042			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.28	9.6%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.35	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.74	11.5%	<5	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.46	9.6%	<5	<100
(DGXN)地高辛					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.47	10%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5711	10.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6338	10.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6647	10%	<5	<100
(E2)雌二醇					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	946	4.2%	6	128
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	955	3.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1025	7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	808	10.1%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	924	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	968	10.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1042	9.9%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	282	7.1%	<5	<100
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	39.7	10%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	45.9	9.4%	<5	<100
(FA)叶酸					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.38	4%	6	137
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.73	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.37	10.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.48	10%	<5	<100
(Ferritin)铁蛋白					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	127	7.1%	5	137
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	216	5.1%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	189	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	137	10.2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	274	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	124	5.9%	<5	<100
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	210	12%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	153	10.5%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.09	7.4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.43	10.4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.64	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.74	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.92	10%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
4 / 8								

		批号：3385042			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(FSH)促卵泡刺激素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	25.2	7.1%	6	135
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	25.9	5.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	24.6	7.8%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	38.9	9.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.8	10.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.8	11.6%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	22.4	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	23.6	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.8	9.8%	<5	<100
(FT3)游离T3					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	11.7	2.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	12.2	7.5%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.5	0.8%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	pmol/L	14.8	6.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	9.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1.99	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.36	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	10.6%	<5	<100
(FT4)游离T4					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33.4	4.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	32.2	8.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	34.2	10.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.5	9.8%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.5	9.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.42	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	381	12.1%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	pmol/L	343	15%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	32.2	9.9%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	25	12%	<5	<100
(hCG)人绒毛膜促性腺激素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	135	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	137	1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	152	9.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	138	2.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	147	9.5%	<5	<100
(HCY)同型半胱氨酸					
TOSHIBA TBA TOSHIBA \ Enzymatic Cycling Assays Method	μmol/L	20.8	9.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.1	10.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	21.7	10%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μmol/L	20.5	20%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	24.3	9.5%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
5 / 8								

批号：3385042					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(HGH)人生长激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	3.97	10.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.07	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	3.53	9.6%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	82.6	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.6	10.2%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	237	10%	<5	<100
(INS)胰岛素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	34.3	4.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32	7.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.1	10.6%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.7	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	25.1	9.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.7	10.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.2	10.3%	<5	<100
(LH)促黄体生成素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	16	9.5%	6	108
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.4	3.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.1	7%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.1	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	16.6	12.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	17.9	10.6%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	28.5	9.8%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	16.3	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.3	9.8%	<5	<100
(Myo / Mb)肌红蛋白					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	126	9.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	143	9.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	135	5.7%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	100	1.2%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	139	10.1%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1269	4.7%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	149	4.7%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1184	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	139	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1246	0.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	147	0.8%	<5	<100
(PRL)催/ 泌乳素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.8	4.8%	6	127

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
6 / 8								

		批号：3385042			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.4	5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.1	7.9%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.2	9.9%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.9	8.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.1	9.9%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.8	10%	<5	<100
(PROG)孕酮					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.98	6.5%	6	133
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.7	3.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.7	8.2%	<5	<100
Siemens Immulite1000/2000 Siemens Dade Actin FS(APTT) \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.23	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	9.8%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.36	5.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.6	10.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	10.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.7	10%	<5	<100
(PSA)前列腺特异性抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14	8.5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.9	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.5	11.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.5	11.3%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.6	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.7	9.4%	<5	<100
(PTH)甲状旁腺素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	8.89	3.8%	6	136
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.78	5.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.84	10.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.46	10.3%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	2.16	8%	<5	<100
(SHBG)性激素结合球蛋白					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	78.8	12%	<5	<100
(TES)睾酮					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.638	6%	6	137
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.53	6.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.541	15.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.82	11%	<5	<100
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.59	10.2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.59	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.42	9.5%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
7 / 8								

批号：3385042					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Tg)甲状腺球蛋白					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.83	7.8%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.37	9.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.3	12.3%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	21.5	10%	<5	<100
(Tnl)肌钙蛋白I					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.37	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.47	12.8%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.336	1.6%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.889	2%	<5	<100
(TSH)促甲状腺素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.8	2.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.8	7.6%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.13	1.6%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.15	0.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.7	11.2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.15	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.76	10.3%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.78	8.8%	<5	<100
(TT3)总三碘甲状腺氨酸					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.32	3.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.49	7.1%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.64	12%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.82	9.9%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.12	1.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.15	10.2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.02	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.17	2.6%	<5	<100
(TT4)总T4					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	115	2.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	134	7.8%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	103	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	9.2%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	ng/mL	129	8.9%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	131	5.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	8.5%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	136	9.6%	<5	<100
(VB12)维生素B12					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	530	5.9%	6	100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1056	3.6%	<5	<100

REF	540-2	6×3mL	EXP	2026-09-20	LOT	3385042	更新日期：	2024-04
8 / 8								

		批号：3385042			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	771	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	318	1.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	567	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	659	10.3%	<5	<100
(β2M)β-2-微球蛋白					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1440	11.1%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	1180	10%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.166	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.136	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.139	10.1%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。