

REF

670-1

6×1mL

EXP

2023-05-12

LOT

3FE85AB

更新日期：2024-05

1 / 1

		批号：3FE85AB			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	31	2.6%	8	485
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	36.4	3.3%	5	134
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	37.5	2.4%	5	123
Sysmex CS 5100 Siemens \ Light scattering	Seconds	31.4	5.1%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	33.5	9.6%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	33.9	9.7%	<5	<100
Stago STA series Dedicated \ Magnetic Beads	Seconds	32.4	9.9%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	%	85.2	3.2%	5	128
Sysmex CS-1300 SUNBIO \ OTHER	%	88.5	3.4%	5	111
Werfen ACL series Dedicated \ Nephelometry	%	92.6	9.6%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	%	93.8	9.4%	<5	<100
(FIB)纤维蛋白原					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	g/L	3.25	5.9%	8	483
Mindray BC-3000/2800/3600/2300/30/20 series KHB \ Magnetic Beads	g/L	2.85	2.5%	8	218
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	g/L	3.34	2.9%	5	133
Sysmex CS-1300 SUNBIO \ OTHER	g/L	3.32	2.9%	5	116
Sysmex CS 5100 Siemens \ Light scattering	g/L	3.53	4.2%	<5	<100
Stago STA series Dedicated \ Magnetic Beads	g/L	3.79	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	3.56	7.9%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	3.68	10.1%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	11.7	6.9%	8	485
Mindray BC-3000/2800/3600/2300/30/20 series KHB \ Magnetic Beads	Seconds	9.82	1.1%	8	216
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	12.9	1.5%	5	123
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	12.8	1.3%	5	115
Sysmex CS 5100 Siemens \ Light scattering	Seconds	13.9	22.2%	<5	149
Stago STA series Dedicated \ Light scattering	Seconds	12.6	9.5%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	12.3	9.8%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	12.4	9.7%	<5	<100
(TT)凝血酶时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	19.2	4.6%	9	508
Mindray BC-3000/2800/3600/2300/30/20 series KHB \ Magnetic Beads	Seconds	18.2	1.8%	8	219
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	17.5	2.1%	5	113
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	16.5	3%	5	102
Sysmex CS 5100 Siemens \ Light scattering	Seconds	20.3	4.2%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	15.8	10.1%	<5	<100
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	17.3	9.8%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。