

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期: 2025-03
							1 / 13

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(CK)肌酸激酶	(CK-MB)肌酸激酶-MB (活性)	(Myo / Mb)肌红蛋白
(TnI)肌钙蛋白I	(TnT)肌钙蛋白T	(NT-proBNP)N-末端脑钠肽前体
(DGXN)地高辛	(THEO)茶碱	(TT3)总三碘甲状腺原氨酸
(FT3)游离T3	(TT4)总T4	(FT4)游离T4
(TSH)促甲状腺素	(Cortisol)皮质醇	(FSH)促卵泡刺激素
(LH)促黄体生成素	(PROG)孕酮	(PRL)催/泌乳素
(TES)睾酮	(E2)雌二醇	(C-P)C-肽
(FA)叶酸	(INS)胰岛素	(VB12)维生素B12
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(CA 125)糖类抗原 125
(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9	(β2M)β-2-微球蛋白
(Ferritin)铁蛋白	(β-hCG)人绒毛膜促性腺激素β亚基	(F-PSA)游离前列腺特异性抗原
(HCY)同型半胱氨酸	(17-OHP)17-羟孕酮	(ALD)醛固酮
(Anti-Tg)抗甲状腺球蛋白抗体	(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮
(EPO)促红细胞生成素	(HGH)人生长激素	(IGF1)胰岛素生长因子1
(Ost)骨泌素	(PTH)甲状旁腺素	(SHBG)性激素结合球蛋白
(Tg)甲状腺球蛋白	(25-OH-Vit D)25羟基维生素D	(25-OH-VD3)25羟基维生素D3
(PCT)降钙素原	(25-OH-VD2)25羟基维生素D2	(IntactPTH)全段甲状旁腺激素
(HS-CTnI) (高) 超敏肌钙蛋白I	(HS-CTnT)超敏肌钙蛋白T	(CK-MB mass)肌酸激酶-MB质量

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03 2 / 13
------------	-------	-------	------------	------------	------------	---------	-------	-------------------

批号: 33828C6			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Agilent agilent \ MS	ng/mL	0.759	0.578-0.94
(25-OH-VD2)25羟基维生素D2			
OTHER Dedicated \ OTHER	ng/mL	0.03	0.028-0.032
OTHER Dedicated \ OTHER	nmol/L	0.072	0.067-0.077
(25-OH-VD3)25羟基维生素D3			
OTHER Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.1	18.5-21.7
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20	15.2-24.8
(25-OH-Vit D)25羟基维生素D			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.9	15.8-32
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	59.8	39.6-80
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.8	16.6-27
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	54.4	41.5-67.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.7	16.6-24.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	51.6	41.5-61.7
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.9	20.6-31.2
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	64.8	51.5-78.1
(AFP)甲胎蛋白			
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.84	5.24-6.44
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.32	5.97-6.67
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.44	6.25-8.63
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.96	6.35-7.57
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.14	6.38-7.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.84	6.21-7.47
OTHER Abbott \ OTHER	ng/mL	7.09	5.84-8.34
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.41	4.87-7.95
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.6	5.94-7.26
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.34	5.59-9.09
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.26	6.82-7.7
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.58	5.72-7.44
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.61	6.64-8.58
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.51	4.95-8.07
(ALD)醛固酮			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	134	107-161
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pg/mL	70.3	53.6-87
(Anti-Tg)抗甲状腺球蛋白抗体			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	11.2	10-12.4
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	9.71	5.36-14.06
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	19.1	4.9-33.3
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	21	17.2-24.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	1.68	1.21-2.15

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								3 / 13

Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	3.55	2.73-4.37
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	17.4	17.1-17.7
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	6.3	5.42-7.18
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	9.91	8.07-11.75
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.45	1.21-1.69
OTHER Abbott \ OTHER	IU/mL	6.49	6.49-6.49

(Anti-TPO)抗甲状腺过氧化物酶抗体

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	12.8	10.7-14.9
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	7.95	6.71-9.19
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15.7	10.6-20.8
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	10.5	5.9-15.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14.6	11.6-17.6
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15.4	10-20.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	4.78	4.07-5.49
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	6.7	5.48-7.92
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	9.2	8.7-9.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	8.31	6.66-9.96
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	7.26	6.1-8.42

(CA 125)糖类抗原 125

Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.3	28.4-32.2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	31.1	24.5-37.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	36.9	33-40.8
OTHER Abbott \ OTHER	IU/mL	39.9	33.5-46.3
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	34.8	32.9-36.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	34.2	28.9-39.5
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.3	25.9-52.7
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	48.5	37-60
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	29.8	28.7-30.9
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	34.9	31.4-38.4
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	35.5	31.7-39.3
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	43.1	32.8-53.4

(CA 15-3)糖类抗原 15-3

Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.6	26.4-30.8
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.1	27.2-33
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14	11.3-16.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.6	25.8-35.4
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.4	22.4-36.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	18.7	12-25.4
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	26.6-29.8
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	26.7	24.1-29.3
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	31.3	26.4-36.2

(CA 19-9)糖类抗原 19-9

Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	26.7	22.2-31.2
---	-------	------	-----------

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								4 / 13

Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.1	23.8-32.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	24.2	22.5-25.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	23.4-33
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	20.6	15.7-25.5
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	21.5	12.3-30.7
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	27-29.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	40.8	36.1-45.5
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	33.6	30.2-37
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	27.1	17.5-36.7
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.3	34.4-56.2

(CEA)癌胚抗原

Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.8	3.47-4.13
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.32	2.73-3.91
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.22	2.96-3.48
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.71	2.87-4.55
OTHER Abbott \ OTHER	ng/mL	2.89	2.02-3.76
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.63	2.05-3.21
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.9	3.73-6.07
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	2.46-4.58
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.05	2-2.1
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.99	2.27-3.71

(CK)肌酸激酶

Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	130	104-156
--	-----	-----	---------

(CK-MB mass)肌酸激酶-MB质量

Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.91	2.37-3.45
---	-------	------	-----------

(CK-MB)肌酸激酶-MB (活性)

Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.93	2.54-3.32
Beckman AU series Dedicated \ Immune suppression	U/L	9.5	7.77-11.23
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.1	2.88-5.32
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.79	1.78-5.8
TOSHIBA TBA MAKER \ Immune suppression	U/L	8.51	7.28-9.74
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.49	3.6-5.38
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	3.87-6.29
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.25	7.42-11.08
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.09	5.93-8.25
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.87	2.31-3.43

(Cortisol)皮质醇

Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	102	91-113
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	92.4	70.4-114.4
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	115	88-142
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	134	107-161
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	94.3	78.6-110
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	99	97.6-100.4

(C-P)C-肽

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.9	0.68-1.12
---	-------	-----	-----------

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								5 / 13

YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.903	0.834-0.972
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.814	0.746-0.882
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.18	1.14-1.22
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.4	0.299-0.501
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.06	0.85-1.27
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.648	0.494-0.802
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.997	0.941-1.053
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.956	0.891-1.021
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.93	0.75-1.11

(DGXN)地高辛			
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.932	0.641-1.223
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.88	0.704-1.056
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.87	0.663-1.077

(DHEA-S)硫酸脱氢表雄酮			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3364	3115-3613
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3390	2583-4197
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3416	2740-4092
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3906	3094-4718

(E2)雌二醇			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	342	284-400
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	292	265-319
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	288	259-317
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	346	272-420
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	277	218-336
OTHER Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	343	305-381
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	335	313-357
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	342	313-371
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	331	315-347
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	263	200-326
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	246	187-305
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	312	281-343
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	298	257-339
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	275	212-338
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	268	203-333
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	365	314-416

(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	6.77	4.43-9.11
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	6.14	5.39-6.89
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	6.76	5.42-8.1

(FA)叶酸			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	1.84-5.2
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.79	2.31-3.27
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.81	2.71-2.91

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								6 / 13

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.47	3.13-3.81
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.66	2.79-4.53
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.87	2.49-3.25
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.04	2.54-3.54
(Ferritin)铁蛋白			
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	76.9	63.4-90.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	66.5	56.8-76.2
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	73.9	53.8-94
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	105	91-119
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	80.5	61.3-99.7
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	119	95-143
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	72.8	62-83.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	105	103-107
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	194	171-217
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	95.9	72.9-118.9
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	69.2	55.6-82.8
(F-PSA)游离前列腺特异性抗原			
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.97	1.81-2.13
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.75	1.56-1.94
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.16	1.95-2.37
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.74	1.57-1.91
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.06	1.57-2.55
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.48	1.3-1.66
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.93	1.68-2.18
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.81	1.38-2.24
OTHER Abbott \ OTHER	ng/mL	2.05	1.92-2.18
(FSH)促卵泡刺激素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.37	3.83-6.91
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.89	4.6-5.18
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.47	4.33-6.61
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.66	5.47-5.85
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.66	4.89-6.43
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.15	4.61-5.69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.85	5.31-6.39
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.92	4.51-7.33
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.8	3.66-5.94
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.4	4.98-5.82
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	6.34	5.41-7.27
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.3	4.97-5.63
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.77	5.22-6.32
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.34	4.34-6.34
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.51	3.43-5.59
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.07	4.37-5.77

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								7 / 13

(FT3)游离T3			
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.81	4.22-5.4
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.54	4.56-6.52
Roche cobas e411			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.92	5.94-7.9
Beckman AU series			
Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.1	4.5-5.7
Siemens ADVIA Centaur XP/XPT			
Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.13	5.71-6.55
Roche cobas 4000/6000/8000			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.66	5.35-5.97
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.07	5.06-7.08
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.36	4.6-6.12
Abbott Alinity ci-series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.43	4.81-6.05
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.5	3.43-5.57
Tosoh AIA Series			
TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.18	3.95-6.41
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	4.78	4.11-5.45
Roche cobas E601/602			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.72	6.32-7.12
Caris200			
WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.99	3.79-6.19
Mindray CL-8000i			
Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.13	5.27-6.99
OTHER			
Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.51	6.25-6.77
(FT4)游离T4			
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.8	11.1-14.5
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.7	16.1-19.3
Beckman AU series			
Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.8	10.5-15.1
Roche cobas e411			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.7	14.4-17
Siemens ADVIA Centaur XP/XPT			
Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.4	13.8-17
Roche cobas 4000/6000/8000			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.1	13.8-16.4
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.3	12.5-18.1
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.7	11.3-14.1
Abbott Alinity ci-series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.4	11-13.8
Tosoh AIA Series			
TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	14.1	10.7-17.5
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	9.4-15.2
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	11.6	10-13.2
Roche cobas E601/602			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	18.5	18.1-18.9
Caris200			
WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	9.3-15.3
OTHER			
Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	16	14.8-17.2
Mindray CL-8000i			
Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20.9	18-23.8
(hCG)人绒毛膜促性腺激素			
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	0.031	0.029-0.033
OTHER			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.62	3.34-3.9
OTHER			
ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.44	2.73-4.15
Tosoh AIA Series			
TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.22	3.22-5.22
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	0.029	0.024-0.034
Siemens Immulite1000/2000			
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.74	0.44-9.04
(HCY)同型半胱氨酸			
TOSHIBA TBA			
Bio-assay \ Enzymatic Cycling Assays Method	μmol/L	11.7	10.9-12.5
Siemens ADVIA Centaur XP/XPT			
Siemens \ Chemiluminescence Immunoassay (CLIA)	μmol/L	10.8	9.4-12.2
Beckman AU series			
Dedicated \ Enzymatic Cycling Assays Method	μmol/L	12.1	8.3-15.9

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								8 / 13

Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.9	11.1-16.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	12.1	9.7-14.5
(HGH)人生长激素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	1.51	1.15-1.87
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	1.99	1.6-2.38
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	1.41	1.25-1.57
(HS-CTnI) (高) 超敏肌钙蛋白I			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.061	0.047-0.075
(HS-CTnT)超敏肌钙蛋白T			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.019	0.015-0.023
(IgE)免疫球蛋白 IgE			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	169	133-205
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	131	105-157
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	160	122-198
(IGF1)胰岛素生长因子1			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	218	177-259
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	17.7	13.5-21.9
(INS)胰岛素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	7.34	5.5-9.18
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.25	7.39-9.11
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.39	7.79-8.99
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	6.72	5.95-7.49
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.98	9.42-10.54
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.6	7.53-9.67
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.01	7.23-10.79
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	6.95	5.3-8.6
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	8.15	8.02-8.28
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	12	10.9-13.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.83	7.06-10.6
(IntactPTH)全段甲状旁腺激素			
Abbott ARCHITECT i series Abbott \ OTHER	pg/mL	2.98	2.29-3.67
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	3.62	2.87-4.37
(LH)促黄体生成素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.51	3.03-3.99
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.38	4.99-5.77
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.62	4.48-6.76
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.81	3.4-4.22
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.46	2.85-4.07
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.41	3.08-3.74
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.47	5.29-5.65
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.62	2.76-4.48
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.79	2.89-4.69
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.92	3.44-4.4
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.43	5.14-5.72

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								9 / 13

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.76	3.13-4.39
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.17	4.61-5.73
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.37	4.88-5.86
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.93	4.51-7.35
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.57	3.93-5.21
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.66	3.58-5.74

(Myo / Mb)肌红蛋白

Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	40.7	37.9-43.5
TOSHIBA TBA OTHER \ OTHER	µg/L	42.8	29.9-55.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	28.9	25.3-32.5
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	26.9	22.1-31.7
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	29.1	23.3-34.9
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	µg/L	34.2	26.1-42.3
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	µg/L	28.5	24.5-32.5
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	µg/L	39.3	34.1-44.5
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	µg/L	41.8	37-46.6
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	µg/L	28.8	23.6-34

(NT-proBNP)N-末端脑钠肽前体

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	182	163-201
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	21.4	19.2-23.6
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	185	168-202
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	21.8	19.8-23.8
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	149	136-162
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.5	15.9-19.1
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	202	162-242
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	23.8	19.1-28.5
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pg/mL	234	221-247
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	27.5	26-29
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pg/mL	83.9	72.2-95.6
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.87	8.49-11.25
J&J Vitros ECL Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	133	133-133
J&J Vitros ECL Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.6	15.6-15.6

(Ost)骨泌素

Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.1	11.3-14.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.3	12.7-15.9

(PCT)降钙素原

Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.471	0.42-0.522
---	-------	-------	------------

(PRL)催/泌乳素

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.15	8.05-10.25
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.72	9.27-10.17
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.25	8.81-9.69
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.89	7.04-8.74
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.74	7.01-8.47
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	9.4-11.4

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								10 / 13

Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.42	6.68-8.16
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.17	5.46-8.88
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.01	6.87-11.15
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.83	7.08-10.58
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.3	5.56-7.04
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.24	8.11-8.37
Siemens Immulite1000/2000 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.64	5.95-7.33
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.04	6.91-11.17
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.37	8.08-10.66
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.12	6.93-11.31

(PROG)孕酮

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.983	0.354-1.612
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.09	0.8-1.38
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.04	1.87-2.21
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.94	0.823-1.057
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.848	0.728-0.968
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.023	0.02-0.026
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.916	0.335-1.497
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.27	1-1.54
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.896	0.599-1.193
OTHER Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.15	0.95-1.35
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.45	1.18-1.72
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.29	1.03-1.55
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.22	0.93-1.51
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.01	0.77-1.25
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.21	0.75-1.67
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.78	0.747-0.813
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.17	0.98-1.36
Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.54	1.3-1.78

(PSA)前列腺特异性抗原

YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.1	4.79-7.41
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.55	3.49-5.61
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.59	4.34-4.84
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.08	3.11-5.05
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.16	4.04-4.28
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.52	3.33-3.71
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.22	3.71-4.73
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.05	3.49-4.61
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.44	3.72-5.16
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.71	3.73-5.69
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.33	3.29-5.37
OTHER Abbott \ OTHER	ng/mL	4.16	3.59-4.73

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								11 / 13

(PTH)甲状旁腺素				
Siemens ADVIA Centaur XP/XPT		pmol/L	1.83	1.34-2.32
Siemens \ Chemiluminescence Immunoassay (CLIA)				
Beckman UniCel Dxl 600/800		pmol/L	1.91	1.65-2.17
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Roche cobas E601/602		pmol/L	2.01	1.69-2.33
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Tosoh AIA Series		pmol/L	1.94	1.48-2.4
TOSOH \ Chemiluminescence Immunoassay (CLIA)				
Roche cobas E601/602		pmol/L	2.05	1.64-2.46
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman UniCel Dxl 600/800		pmol/L	2	1.66-2.34
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)				
Roche cobas 4000/6000/8000		pmol/L	1.95	1.91-1.99
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman UniCel Dxl 600/800		pmol/L	2.26	1.82-2.7
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)				
Abbott ARCHITECT i series		pmol/L	1.96	1.58-2.34
Dedicated \ Chemiluminescence Immunoassay (CLIA)				
OTHER		pmol/L	1.88	1.59-2.17
Siemens \ Chemiluminescence Immunoassay (CLIA)				
(SHBG)性激素结合球蛋白				
Beckman UniCel Dxl 600/800		nmol/L	33.2	25.2-41.2
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)				
(TES)睾酮				
Roche cobas 4000/6000/8000		ng/mL	6.45	5.71-7.19
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Abbott ARCHITECT i series		ng/mL	6.39	5.18-7.6
Abbott \ Chemiluminescence Immunoassay (CLIA)				
Roche cobas 4000/6000/8000		ng/mL	6.71	6.39-7.03
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Roche cobas e411		ng/mL	6.02	5.55-6.49
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
OTHER		ng/mL	6.31	5.43-7.19
ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman UniCel Dxl 600/800		ng/mL	5.53	5.12-5.94
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman AU series		ng/mL	5.29	4.8-5.78
Dedicated \ Chemiluminescence Immunoassay (CLIA)				
Beckman UniCel Dxl 600/800		ng/mL	5.13	4.62-5.64
Dedicated \ Chemiluminescence Immunoassay (CLIA)				
Beckman UniCel Dxl 600/800		ng/mL	5.28	4.02-6.54
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)				
Roche cobas E601/602		ng/mL	6.17	4.91-7.43
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Tosoh AIA Series		ng/mL	6.61	5.04-8.18
TOSOH \ Chemiluminescence Immunoassay (CLIA)				
Roche cobas E601/602		ng/mL	6.13	6.02-6.24
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Siemens Immulite1000/2000		ng/mL	6.03	5.19-6.87
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)				
Mindray CL-8000i		ng/mL	9.01	7.75-10.27
Mindray \ Chemiluminescence Immunoassay (CLIA)				
(Tg)甲状腺球蛋白				
Beckman UniCel Dxl 600/800		ng/mL	9.97	9.21-10.73
Dedicated \ Chemiluminescence Immunoassay (CLIA)				
Roche cobas E601/602		ng/mL	10.7	8.7-12.7
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman UniCel Dxl 600/800		ng/mL	8.26	6.39-10.13
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)				
Abbott ARCHITECT i series		ng/mL	9.88	7.92-11.84
Abbott \ Chemiluminescence Immunoassay (CLIA)				
Beckman UniCel Dxl 600/800		ng/mL	8.04	6.4-9.68
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)				
(THEO)茶碱				
Agilent		µg/mL	8.7	6.63-10.77
agilent \ MS				
(TnI)肌钙蛋白I				
Ortho Vitros 3600/5600		ng/mL	0.079	0.066-0.092
Dedicated \ Chemiluminescence Immunoassay (CLIA)				
Beckman UniCel Dxl 600/800		ng/mL	0.066	0.019-0.113
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)				
J&J Vitros ECL		ng/mL	0.075	0.063-0.087
Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)				
Beckman UniCel Dxl 600/800		ng/mL	0.07	0.056-0.084
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)				
Mindray CL-8000i		ng/mL	0.005	0.004-0.006
Mindray \ Chemiluminescence Immunoassay (CLIA)				

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								12 / 13

(TnT)肌钙蛋白T

Roche cobas E601/602				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02		0.018-0.022
Roche cobas e411				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02		0.016-0.024
Roche cobas E601/602				
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02		0.018-0.022

(TSH)促甲状腺素

Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.76		2.25-3.27
Beckman UniCel Dxl 600/800				
Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.12		2.87-3.37
Roche cobas 4000/6000/8000				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.34		3.15-3.53
Roche cobas 4000/6000/8000				
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.34		2.99-3.69
Roche cobas e411				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.5		3.26-3.74
Siemens ADVIA Centaur XP/XPT				
Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3		2.8-3.2
OTHER				
Abbott \ OTHER	mIU/L	2.49		2.2-2.78
Roche cobas E601/602				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.48		3.15-3.81
Beckman AU series				
Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.96		2.66-3.26
Abbott Alinity ci-series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.39		2.16-2.62
Tosoh AIA Series				
TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.87		2.19-3.55
Beckman UniCel Dxl 600/800				
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.9		2.21-3.59
Beckman UniCel Dxl 600/800				
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.83		2.54-3.12
Roche cobas E601/602				
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.34		3.2-3.48
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.63		2.23-3.03
Caris200				
WanTai \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3		2.28-3.72
OTHER				
Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.86		2.61-3.11
Mindray CL-8000i				
Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.45		2.01-2.89

(TT3)总三碘甲状腺原氨酸

Beckman UniCel Dxl 600/800				
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.13		0.98-1.28
Roche cobas 4000/6000/8000				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.46		1.1-1.82
Beckman AU series				
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.17		0.99-1.35
Siemens ADVIA Centaur XP/XPT				
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.45		1.24-1.66
Roche cobas E601/602				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.27		1.03-1.51
Roche cobas e411				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.58		1.45-1.71
Abbott Alinity ci-series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.2		1.07-1.33
Beckman UniCel Dxl 600/800				
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.14		0.87-1.41
Beckman UniCel Dxl 600/800				
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.07		0.87-1.27
Abbott ARCHITECT i series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.18		0.78-1.58
Caris200				
WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.956		0.727-1.185
OTHER				
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.32		1.17-1.47

(TT4)总T4

Roche cobas 4000/6000/8000				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	77.2		68.7-85.7
Roche cobas e411				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	64.2		57.8-70.6
OTHER				
Abbott \ OTHER	ng/mL	73.7		64.3-83.1
Roche cobas E601/602				
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	93.5		65.1-121.9
Beckman AU series				
Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	75.4		64.4-86.4
Beckman UniCel Dxl 600/800				
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	69.2		52.7-85.7
Beckman UniCel Dxl 600/800				
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	76.7		63-90.4
Abbott Alinity ci-series				
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	59.7		53.1-66.3

REF	540-1	6x3mL	EXP	2025-02-16	LOT	33828C6	更新日期:	2025-03
								13 / 13

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA) Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	58.1	46.2-70
	ng/mL	55.4	42.1-68.7

(VB12)维生素B12

Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA) maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA) Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA) Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA) YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA) Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA) Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA) Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA) Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	347	278-416
	pg/mL	485	428-542
	pg/mL	610	471-749
	pg/mL	389	349-429
	pg/mL	312	177-447
	pg/mL	324	247-401
	pg/mL	509	388-630
	pg/mL	528	477-579
	pg/mL	383	333-433

(β2M)β-2-微球蛋白

Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA) Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1020	818-1222
	ng/mL	1170	892-1448

(β-hCG)人绒毛膜促性腺激素β亚基

Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA) Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA) Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA) Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA) Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA) Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA) Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA) Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA) OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA) Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA) Mindray CL-8000i Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.006	0.005-0.007
	IU/mL	0.006	0.006-0.006
	IU/mL	0.004	0.003-0.005
	IU/mL	0.003	0.002-0.004
	IU/mL	0.004	0.003-0.005
	IU/mL	0.004	0.004-0.004
	IU/mL	0.003	0.003-0.003
	IU/mL	0.003	0.003-0.003
	IU/mL	0.003	0.002-0.004
	IU/mL	0.004	0.003-0.005
	IU/mL	0.005	0.004-0.006
	IU/mL	0.005	0.004-0.006

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。