

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
1 / 6								

		批号：336C9A1			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(25-OH-Vit D)25羟基维生素D					
Roche cobas E601/602	ng/mL	35.5	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	nmol/L	88.5	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(AFP)甲胎蛋白					
maglumi system	ng/mL	129	7.9%	9	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	118	8%	8	105
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)					
Roche cobas E601/602	ng/mL	137	11.8%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	ng/mL	135	4.3%	5	125
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	ng/mL	113	7.5%	<5	<100
Snibe \ OTHER					
Maccura i3000	ng/mL	109	4.6%	<5	<100
MAKER \ OTHER					
Mindray CL-2000i/1000i/6000i series	ng/mL	126	3.5%	15	325
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	145	5.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	ng/mL	120	1.9%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	104	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	144	3.1%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS	pg/mL	3217	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	pg/mL	3171	5.3%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
(Anti-Tg)抗甲状腺球蛋白抗体					
YHLO iFlash Analyzer	IU/mL	139	5.2%	8	218
YHL \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	IU/mL	131	10%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	17.6	0.6%	<5	<100
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	20.2	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	76.3	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Dirui CM Series	IU/mL	25	10%	<5	<100
DIRUI \ Chemiluminescence Immunoassay (CLIA)					
(Anti-TPO)抗甲状腺过氧化物酶抗体					
YHLO iFlash Analyzer	IU/mL	28.4	7.3%	9	243
YHL \ Chemiluminescence Immunoassay (CLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	33.1	6.6%	10	205
Mindray \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	IU/mL	97.1	10%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	34	3.8%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	26.4	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Dirui CM Series	IU/mL	33.4	10%	<5	<100
DIRUI \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	33.1	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(CA 125)糖类抗原 125					
maglumi system	IU/mL	52.3	8.5%	9	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	49.2	11.9%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	IU/mL	51	3.6%	5	124
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	IU/mL	51.1	3.5%	<5	<100
Snibe \ OTHER					
Maccura i3000	IU/mL	86.3	9.2%	<5	<100
MAKER \ OTHER					
Mindray CL-2000i/1000i/6000i series	IU/mL	80.6	3.7%	15	314
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	51.5	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	IU/mL	61.7	8.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	51.2	2.1%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(CA 15-3)糖类抗原 15-3					
maglumi system	IU/mL	62.4	11.4%	6	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
2 / 6								

		批号：336C9A1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Mindray CL-2000i/1000i/6000i series	IU/mL	25.3	7.2%	6	130
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	48.1	9.4%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	IU/mL	76.2	4.2%	5	125
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	IU/mL	64.3	4.4%	<5	<100
Snibe \ OTHER					
Maccura i3000	IU/mL	48.6	2.5%	<5	<100
MAKER \ OTHER					
Roche cobas e 801	IU/mL	58.6	11.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	47.7	4.4%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel DxI 600/800	IU/mL	23.3	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
(CA 19-9)糖类抗原 19-9					
maglumi system	IU/mL	58.4	11.3%	7	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	45.2	10.5%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	IU/mL	56	4.8%	5	125
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	IU/mL	54.6	6.1%	<5	<100
Snibe \ OTHER					
Maccura i3000	IU/mL	33.3	2.6%	<5	<100
MAKER \ OTHER					
Mindray CL-2000i/1000i/6000i series	IU/mL	43.2	4.9%	15	333
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	IU/mL	39.3	11.4%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel DxI 600/800	IU/mL	76.3	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	43.5	2.6%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(CEA)癌胚抗原					
maglumi system	ng/mL	7.89	12.2%	9	108
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	8.37	11%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	ng/mL	7.54	5.2%	5	125
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	ng/mL	7.4	7%	<5	<100
Snibe \ OTHER					
Maccura i3000	ng/mL	9.79	3.2%	<5	<100
MAKER \ OTHER					
Mindray CL-2000i/1000i/6000i series	ng/mL	6	3.9%	15	324
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	8.39	2%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	8.91	7.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel DxI 600/800	ng/mL	8.2	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
(CK)肌酸激酶					
Roche cobas 4000/6000/8000	U/L	87.9	12%	<5	<100
Roche \ Urease-UV-rate					
(CK-MB)肌酸激酶-MB (活性)					
Roche cobas e411	ng/mL	12.2	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	U/L	15.9	12%	<5	<100
Roche \ Immunoturbidimetry (ITA)					
(Cortisol)皮质醇					
Autolumo A2000/A2000 PLUS	ng/mL	241	1.1%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	244	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	ng/mL	206	2.8%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
(C-P)C-肽					
YHLO iFlash Analyzer	ng/mL	1.97	5.1%	8	214
YHL \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	2.17	14.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Mindray CL-2000i/1000i/6000i series	ng/mL	1.67	3.4%	<5	<100
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	2.17	9%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(DGXN)地高辛					
Roche cobas e411	ng/mL	2.51	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(DHEA-S)硫酸脱氢表雄酮					
Roche cobas E601/602	ng/mL	6760	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
3 / 6								

		批号：336C9A1			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(E2)雌二醇					
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1047	4%	9	229
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	843	8.5%	8	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	972	11.6%	8	107
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	780	10.7%	5	125
maglumi system Snibe \ OTHER	pmol/L	790	6.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1039	12%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1099	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1037	2.5%	<5	<100
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	40.4	12%	<5	<100
(FA)叶酸					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.59	12%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.72	2.8%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.47	2.1%	<5	<100
(Ferritin)铁蛋白					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	208	3.1%	6	141
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	192	12%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	152	10%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.97	4.5%	6	148
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.24	10.1%	5	<100
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.64	5.4%	5	125
Maccura i3000 MAKER \ OTHER	ng/mL	5.81	4.4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.31	4.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.51	10%	<5	<100
(FSH)促卵泡刺激素					
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	26.5	3.6%	9	218
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	27.6	10.2%	8	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	25.5	9.3%	8	104
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	29.1	6.5%	5	125
maglumi system Snibe \ OTHER	IU/L	28.7	4.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	26.1	12%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/L	27.2	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.1	0.9%	<5	<100
(FT3)游离T3					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.5	10.2%	9	248
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11	9.4%	9	113
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	14.6	26.5%	6	<100
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	11.1	4.7%	5	124
maglumi system Snibe \ OTHER	pmol/L	11.7	4.7%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.6	2.5%	15	302
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.5	10%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.81	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12	9.8%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
4 / 6								

批号：336C9A1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(FT4)游离T4					
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.5	10.2%	9	126
YHLO iFlash Analyzer					
YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.9	32.6%	7	184
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	37.1	19.1%	6	<100
maglumi system					
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	30.9	3.9%	5	125
maglumi system					
Snibe \ OTHER	pmol/L	29.6	3.2%	<5	<100
Mindray CL-2000i/1000i/6000i series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	37.6	28.3%	16	360
Mindray BC-6000 series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	157	10%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	33.9	4.3%	<5	<100
Dirui CM Series					
DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	36.6	10%	<5	<100
(hCG)人绒毛膜促性腺激素					
Mindray CL-2000i/1000i/6000i series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	95.6	9%	9	177
OTHER					
Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	188	4.5%	8	212
Beckman UniCel Dxl 600/800					
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	125	12.3%	8	109
maglumi system					
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	209	3.7%	5	125
Ortho Vitros 3600/5600					
Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	188	3.6%	<5	175
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	192	10%	<5	<100
(HCY)同型半胱氨酸					
Roche cobas 4000/6000/8000					
Roche \ Enzymatic Cycling Assays Method	μmol/L	22.1	12%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.2	6.1%	<5	<100
(HGH)人生长激素					
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.14	12%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Image 800					
Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	84	12%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	241	12%	<5	<100
(INS)胰岛素					
YHLO iFlash Analyzer					
YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.3	11.6%	8	215
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.7	11.3%	5	<100
Mindray CL-2000i/1000i/6000i series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.9	2.5%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.9	2.7%	<5	<100
(LH)促黄体生成素					
OTHER					
Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	19.7	4%	9	227
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.5	13.4%	8	<100
Beckman UniCel Dxl 600/800					
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	16	8.8%	8	106
maglumi system					
Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	33.1	4.4%	5	125
maglumi system					
Snibe \ OTHER	IU/L	33.1	8.1%	<5	<100
Dirui CM Series					
DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.7	10%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	21.3	1.7%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.5	12%	<5	<100
(Myo / Mb)肌红蛋白					
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	127	12%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1204	12%	<5	<100
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	142	12%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
5 / 6								

批号：336C9A1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(PRL)催/ 泌乳素					
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.7	9.3%	8	<100
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.7	3.4%	8	216
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.8	7.9%	8	108
maglumi system Snibe \ OTHER	ng/mL	15.7	6.5%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.8	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.7	3.7%	<5	<100
(PROG)孕酮					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12	5.9%	9	180
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.28	10.7%	8	212
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.54	9.4%	8	117
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	10.2	18%	8	107
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.9	12.5%	5	<100
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11	4.5%	<5	175
maglumi system Snibe \ OTHER	ng/mL	11.1	7.8%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.2	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.1	3%	<5	<100
(PSA)前列腺特异性抗原					
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	11.9	4.7%	5	125
Maccura i3000 MAKER \ OTHER	ng/mL	12.3	2.5%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.8	1%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.4	4.4%	15	324
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.5	12%	<5	<100
(PTH)甲状旁腺素					
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	4.83	2.9%	9	226
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	7.39	12%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.07	3.9%	<5	<100
(TES)睾酮					
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.42	21.3%	8	188
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.625	15.6%	8	107
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.518	23%	6	<100
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.402	21.6%	5	125
maglumi system Snibe \ OTHER	ng/mL	0.005	19.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.502	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.51	4.8%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.444	10%	<5	<100
(Tg)甲状腺球蛋白					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.63	10%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.85	2.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.53	12%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2025-09-27	LOT	336C9A1	更新日期：	2024-04
6 / 6								

批号：336C9A1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(TnI)肌钙蛋白I					
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.376	12%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.02	12%	<5	<100
(TSH)促甲状腺素					
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.4	9.6%	9	111
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	12.3	8.4%	8	214
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	9.2%	6	<100
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	10.4	3.8%	5	125
maglumi system Snibe \ OTHER	mIU/L	10.1	3.2%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	11.1	8.4%	16	384
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.26	10%	<5	<100
Mindray BC-6000 serices Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	5.3	10%	<5	<100
(TT3)总三碘甲状腺原氨酸					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.54	10.7%	9	243
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.18	9%	9	123
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.53	10.7%	5	<100
maglumi system Snibe \ OTHER	ng/mL	2.37	3.5%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.17	22%	15	298
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.58	5.8%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.002	6.9%	<5	<100
Mindray BC-6000 serices Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.08	10%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.88	10%	<5	<100
(TT4)总T4					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	129	7.1%	9	238
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	134	6.5%	9	111
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	136	11.3%	5	<100
maglumi system Snibe \ OTHER	ng/mL	135	3.2%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	138	3.1%	16	368
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	133	10%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	139	1.1%	<5	<100
Mindray BC-6000 serices Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	126	10%	<5	<100
(VB12)维生素B12					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	778	1%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pg/mL	761	1.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	784	12%	<5	<100
(β2M)β-2-微球蛋白					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1172	7.4%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	1200	12%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.133	11.3%	5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.102	4.4%	5	115
maglumi system Snibe \ OTHER	IU/mL	0.198	5.2%	<5	<100
Dirui CM Series DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.203	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。