

REF	539-1 6x3mL	EXP 2024-08-04	LOT 336A22D	更新日期: 2025-02 1 / 5
------------	----------------	-----------------------	--------------------	------------------------

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(THEO)茶碱	(TT3)总三碘甲状腺氨酸	(FT3)游离T3
(TT4)总T4	(FT4)游离T4	(TSH)促甲状腺素
(Cortisol)皮质醇	(FSH)促卵泡刺激素	(LH)促黄体生成素
(PROG)孕酮	(PRL)催/ 泌乳素	(TES)睾酮
(E2)雌二醇	(FA)叶酸	(INS)胰岛素
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(CA 125)糖类抗原 125
(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9	(β2M)β-2-微球蛋白
(Ferritin)铁蛋白	(β-hCG)人绒毛膜促性腺激素β亚基	(F-PSA)游离前列腺特异性抗原
(HCY)同型半胱氨酸	(17-OHP)17-羟孕酮	(Cyfra 21-1)细胞角蛋白19片段
(Anti-Tg)抗甲状腺球蛋白抗体	(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮
(IGF1)胰岛素生长因子1	(PTH)甲状旁腺素	(Tg)甲状腺球蛋白

REF	539-1	6x3mL	EXP	2024-08-04	LOT	336A22D	更新日期:	2025-02
								2 / 5

批号: 336A22D			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens BEP III Dedicated \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	1.59	1.33-1.85
(AFP)甲胎蛋白			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.21	5.77-8.65
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.95	6.46-9.44
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.39	4.33-6.45
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.51	6.65-10.37
(Anti-Tg)抗甲状腺球蛋白抗体			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.89	1.47-2.31
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	2.06	1.64-2.48
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	48.1	38.5-57.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	5.89	4.74-7.04
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	26.2	21-31.4
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	74.2	59.4-89
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	7.08	5.66-8.5
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	10.1	7.9-12.3
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	7.65	6.2-9.1
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	6.56	5.31-7.81
(CA 125)糖类抗原 125			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	37.3	30.9-43.7
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	35.2	27.8-42.6
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	47.7	38.5-56.9
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	33.1	26.7-39.5
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	33.5	27.1-39.9
(CA 15-3)糖类抗原 15-3			
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	26.3	21.1-31.5
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	16.4	11.2-21.6
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	21.9	17.7-26.1
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	22.7	17.7-27.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.2	18-28.4
(CA 19-9)糖类抗原 19-9			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	34.2	27.4-41
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	22.3	17.7-26.9
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	19.6	15.8-23.4
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	25.6	20.2-31
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	21.5	17.3-25.7
(CEA)癌胚抗原			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.53	4.41-6.65

REF	539-1	6x3mL	EXP	2024-08-04	LOT	336A22D	更新日期:	2025-02
								3 / 5

Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.68	3.76-5.6
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.25	4.99-7.51
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.75	3.83-5.67
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.79	3.91-5.67
(Cortisol)皮质醇			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	111	89-133
(Cyfra 21-1)细胞角蛋白19片段			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.21	1.03-1.39
(DHEA-S)硫酸脱氢表雄酮			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2203	1886-2520
(E2)雌二醇			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	194	162-226
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	209	204-214
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	142	116-168
(FA)叶酸			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.26	3.46-5.06
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.22	1.69-2.75
(Ferritin)铁蛋白			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	34	27.1-40.9
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	47.2	37.4-57
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	47.2	37.4-57
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	42.5	41.6-43.4
(F-PSA)游离前列腺特异性抗原			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.16	0.94-1.38
(FSH)促卵泡刺激素			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.64	4.56-6.72
(FT3)游离T3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.23	3.37-5.09
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.15	3.31-4.99
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.11	3.29-4.93
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.14	3.91-6.37
(FT4)游离T4			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	14.2	11.4-17
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.8	10.6-15
Roche Modular E170 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	16	12.6-19.4
(hCG)人绒毛膜促性腺激素			
Tosoh AIA Series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.81	4.83-6.79
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.45	4.33-6.57
Siemens ADVIA Centaur XP/XPT Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.34	3.46-5.22
(HCY)同型半胱氨酸			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	13.7	11.3-16.1
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μmol/L	14.6	11.8-17.4

REF	539-1	6x3mL	EXP	2024-08-04	LOT	336A22D	更新日期:	2025-02
								4 / 5

(IgE)免疫球蛋白 IgE			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	33.8	27-40.6
Beckman Immage 800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	55.8	45-66.6
(IGF1)胰岛素生长因子1			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	152	124-180
(INS)胰岛素			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	12.6	10.4-14.8
(LH)促黄体生成素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.34	3.46-5.22
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.93	4.07-5.79
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.25	3.39-5.11
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.44	2.84-4.04
(PRL)催/ 泌乳素			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.4	6.79-10.01
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.68	6.28-9.08
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.67	7.11-10.23
(PROG)孕酮			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.836	0.672-1
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.77	0.444-1.096
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.843	0.668-1.018
(PSA)前列腺特异性抗原			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.19	3.35-5.03
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.96	3.3-4.62
Roche Modular E170 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.21	3.41-5.01
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.28	2.5-4.06
(PTH)甲状旁腺素			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.93	1.57-2.29
(TES)睾酮			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.32	5.11-7.53
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.5	5.3-7.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.41	6.03-8.79
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.47	5.25-7.69
(Tg)甲状腺球蛋白			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.3	10.7-15.9
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22.3	17.9-26.7
(THEO)茶碱			
Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	µg/mL	10.8	8.6-13
(TSH)促甲状腺素			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	15.2	12-18.4
Abbott ARCHITECT i series Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	13.5	10.7-16.3

REF	539-1	6x3mL	EXP	2024-08-04	LOT	336A22D	更新日期:	2025-02
								5 / 5

(TT3)总三碘甲状腺原氨酸

Siemens ADVIA Centaur XP/XPT Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.21	0.85-1.57
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.03	0.86-1.2
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.09	0.85-1.33

(TT4)总T4

Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	56.1	46.1-66.1
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	63.3	50.5-76.1

(β2M)β-2-微球蛋白

Roche cobas 4000/6000/8000 Dedicated \ Immunonephelometry (INA)	ng/mL	12700	10312-15088
--	-------	-------	-------------

(β-hCG)人绒毛膜促性腺激素β亚基

Roche Modular E170 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	0.004	0.003-0.005
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.004	0.003-0.005
Tosoh AIA Series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.006	0.005-0.007

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。