

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期: 2025-02
							1 / 9

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(TT3)总三碘甲状腺原氨酸	(FT3)游离T3	(TT4)总T4
(FT4 )游离T4	(TSH)促甲状腺素	(Cortisol)皮质醇
(FSH)促卵泡刺激素	(LH)促黄体生成素	(PROG)孕酮
(PRL)催/ 泌乳素	(TES)睾酮	(E2)雌二醇
(C-P)C-肽	(FA)叶酸	(INS)胰岛素
(VB12)维生素B12	(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白
(CEA)癌胚抗原	(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原
(CA 125)糖类抗原 125	(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9
(β2M)β-2-微球蛋白	(Ferritin)铁蛋白	(β-hCG)人绒毛膜促性腺激素β亚基
(F-PSA)游离前列腺特异性抗原	(HCY)同型半胱氨酸	(ALD)醛固酮
(Anti-Tg)抗甲状腺球蛋白抗体	(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮
(HGH)人生长激素	(IGF1)胰岛素生长因子1	(PTH)甲状旁腺素
(Tg)甲状腺球蛋白		

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								2 / 9

批号: 336A2F0			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
<b>(AFP)甲胎蛋白</b>			
Roche cobas e 801	ng/mL	7.06	6.45-7.67
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas E601/602	ng/mL	7.14	5.98-8.3
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
YHLO iFlash Analyzer	ng/mL	5.64	5-6.28
YHL \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	ng/mL	6.39	5.75-7.03
Mindray \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	ng/mL	5.85	4.56-7.14
MAKER \ OTHER			
Beckman UniCel Dxl 600/800	ng/mL	6.47	5.59-7.35
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)			
OTHER	ng/mL	8.46	7.83-9.09
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
maglumi system	ng/mL	7.8	5.52-10.08
Snibe \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	ng/mL	6.49	5.61-7.37
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	ng/mL	7.2	6.35-8.05
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
chemclin LiCA Series	ng/mL	7.79	6.01-9.57
CHIVD \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	ng/mL	6.9	6.36-7.44
MAKER \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	ng/mL	7.58	7.14-8.02
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	ng/mL	6.97	6.97-6.97
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)			
<b>(ALD)醛固酮</b>			
Mindray CL-2000i/1000i/6000i series	pg/mL	55.9	23.9-87.9
Mindray \ Chemiluminescence Immunoassay (CLIA )			
<b>(Anti-Tg)抗甲状腺球蛋白抗体</b>			
Beckman UniCel Dxl 600/800	IU/mL	10.1	8.1-12.1
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	IU/mL	11.9	8.9-14.9
Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)			
maglumi system	IU/mL	104	95-113
Snibe \ Chemiluminescence Immunoassay (CLIA )			
YHLO iFlash Analyzer	IU/mL	86.1	68.4-103.8
YHL \ Chemiluminescence Immunoassay (CLIA )			
OTHER	IU/mL	16.7	14.7-18.7
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	IU/mL	24.8	21.3-28.3
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	IU/mL	20.1	17-23.2
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	IU/mL	24.9	18.9-30.9
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Maccura i3000	IU/mL	21.3	19.8-22.8
MAKER \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	IU/mL	25.2	22.9-27.5
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
<b>(Anti-TPO)抗甲状腺过氧化物酶抗体</b>			
Mindray CL-2000i/1000i/6000i series	IU/mL	16.5	12.9-20.1
Mindray \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	IU/mL	18.6	12.9-24.3
Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)			
YHLO iFlash Analyzer	IU/mL	18.3	15.3-21.3
YHL \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas E601/602	IU/mL	20.4	12.9-27.9
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e411	IU/mL	18.2	11.4-25
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
OTHER	IU/mL	17.6	14.7-20.5
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/mL	10.6	9.4-11.8
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	IU/mL	23.7	17.3-30.1
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Abbott ARCHITECT i series	IU/mL	27.7	25-30.4
Abbott \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	IU/mL	24.7	23.3-26.1
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	20.1	18.6-21.6
MAKER \ Chemiluminescence Immunoassay (CLIA )			
maglumi system	IU/mL	12.4	8.4-16.4
Snibe \ Chemiluminescence Immunoassay (CLIA )			
<b>(CA 125)糖类抗原 125</b>			
Beckman UniCel Dxl 600/800	IU/mL	7.04	6.32-7.76
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas e 801	IU/mL	8.46	7.65-9.27
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								3 / 9

Roche cobas E601/602	IU/mL	8.7	7.08-10.32
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
maglumi system	IU/mL	9.72	7.76-11.68
Snibe \ Chemiluminescence Immunoassay (CLIA )			
YHLO iFlash Analyzer	IU/mL	9.73	8.66-10.8
YHL \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	IU/mL	7.74	6.46-9.02
Mindray \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	8.03	7.32-8.74
MAKER \ OTHER			
OTHER	IU/mL	8.94	8.26-9.62
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/mL	6.11	4.77-7.45
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	IU/mL	8.62	7.84-9.4
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
chemclin LiCA Series	IU/mL	10.1	7.7-12.5
CHIVD \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	IU/mL	8.95	8.29-9.61
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	8.98	8.3-9.66
MAKER \ Chemiluminescence Immunoassay (CLIA )			

**(CA 15-3)糖类抗原 15-3**

Beckman UniCel Dxl 600/800	IU/mL	10.6	9.6-11.6
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas E601/602	IU/mL	12.9	10.2-15.6
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas e 801	IU/mL	13.5	12.6-14.4
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
maglumi system	IU/mL	11.6	9.9-13.3
Snibe \ Chemiluminescence Immunoassay (CLIA )			
YHLO iFlash Analyzer	IU/mL	12.4	11.4-13.4
YHL \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	7.9	7.3-8.5
MAKER \ OTHER			
Mindray CL-2000i/1000i/6000i series	IU/mL	10.7	9.5-11.9
Mindray \ Chemiluminescence Immunoassay (CLIA )			
OTHER	IU/mL	16.5	15.1-17.9
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/mL	9.7	7.95-11.45
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	IU/mL	13.1	11.4-14.8
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
chemclin LiCA Series	IU/mL	12.1	10.1-14.1
CHIVD \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	10.9	9.7-12.1
MAKER \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	IU/mL	11.6	10.6-12.6
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	IU/mL	10.9	10.1-11.7
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(CA 19-9)糖类抗原 19-9**

Beckman UniCel Dxl 600/800	IU/mL	2.8	2.31-3.29
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)			
maglumi system	IU/mL	5.46	2.11-8.81
Snibe \ Chemiluminescence Immunoassay (CLIA )			
YHLO iFlash Analyzer	IU/mL	5.17	4.38-5.96
YHL \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	IU/mL	4.83	3.33-6.33
Mindray \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	3.87	3.36-4.38
MAKER \ OTHER			
Roche cobas E601/602	IU/mL	4.94	4.08-5.8
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
OTHER	IU/mL	5.59	5.17-6.01
DIRUI \ Chemiluminescence Immunoassay (CLIA )			
Beckman UniCel Dxl 600/800	IU/mL	3.04	2.37-3.71
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas 4000/6000/8000	IU/mL	4.73	3.98-5.48
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
chemclin LiCA Series	IU/mL	4.81	3.89-5.73
CHIVD \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	IU/mL	4.04	3.82-4.26
MAKER \ Chemiluminescence Immunoassay (CLIA )			
Lifotronic ecl Series	IU/mL	5.07	3.9-6.24
Lifotronic \ Chemiluminescence Immunoassay (CLIA )			
Roche cobas e411	IU/mL	4.25	3.38-5.12
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			

**(CEA)癌胚抗原**

Beckman UniCel Dxl 600/800	ng/mL	2.62	2.44-2.8
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)			
Roche cobas e 801	ng/mL	3.26	3.03-3.49
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
Roche cobas E601/602	ng/mL	3.18	2.61-3.75
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)			
YHLO iFlash Analyzer	ng/mL	3.16	2.93-3.39
YHL \ Chemiluminescence Immunoassay (CLIA )			
Mindray CL-2000i/1000i/6000i series	ng/mL	2.5	1.99-3.01
Mindray \ Chemiluminescence Immunoassay (CLIA )			
Maccura i3000	ng/mL	3.61	3.4-3.82
MAKER \ OTHER			

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								4 / 9

OTHER								
DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL			3.25				3-3.5
maglumi system								
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL			3.35				2.73-3.97
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL			2.86				2.51-3.21
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			3.07				2.81-3.33
chemclin LiCA Series								
CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL			4.05				3.56-4.54
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL			2.86				2.74-2.98
Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL			2.62				2.22-3.02

**(Cortisol)皮质醇**

Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL			114				99-129
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			119				90-148
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL			102				78-126

**(C-P)C-肽**

maglumi system								
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL			1.52				1.39-1.65
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL			1.07				0.95-1.19
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			1.27				1.21-1.33
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			1.26				1.02-1.5
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			1.24				0.94-1.54
Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL			1.31				1.24-1.38
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL			1.22				1.07-1.37
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL			0.93				0.874-0.986

**(DHEA-S)硫酸脱氢表雄酮**

Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			3780				2873-4687

**(E2)雌二醇**

Siemens ADVIA Chemistry XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L			260				232-288
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L			351				310-392
maglumi system								
Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L			212				180-244
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA )	pmol/L			301				149-453
Roche cobas e411								
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L			280				255-305
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L			221				118-324
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L			253				218-288
OTHER								
Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L			247				189-305
OTHER								
DIRUI \ Chemiluminescence Immunoassay (CLIA )	pmol/L			405				370-440
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L			231				58-404
Siemens ADVIA Centaur XP/XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L			161				143-179
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L			266				202-330
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L			230				214-246

**(FA)叶酸**

Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL			2.9				2.67-3.13
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL			3.15				2.39-3.91
Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL			2.75				2.38-3.12
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL			3.83				2.76-4.9

**(Ferritin)铁蛋白**

Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL			118				102-134

**(F-PSA)游离前列腺特异性抗原**

Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL			0.089				0.074-0.104
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL			0.101				0.085-0.117
YHLO iFlash Analyzer								
YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL			0.124				0.105-0.143

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								5 / 9

Maccura i3000 MAKER \ OTHER	ng/mL	0.085	0.076-0.094
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.097	0.085-0.109
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.134	0.11-0.158
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.143	0.108-0.178
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.178	0.12-0.236
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.084	0.074-0.094
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.113	0.086-0.14

**(FSH)促卵泡刺激素**

Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	6.15	5.45-6.85
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.82	5.34-6.3
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.45	4.83-6.07
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.28	4.69-5.87
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.38	5.15-5.61
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.48	4.22-6.74
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.2	4.54-5.86
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.91	4.38-5.44
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.88	3.55-4.21
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.81	4.27-5.35
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.07	3.85-6.29
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	7.5	6.84-8.16

**(FT3)游离T3**

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	6	4.99-7.01
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	6.72	6.09-7.35
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.03	4.37-5.69
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.79	5.45-6.13
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.23	4.32-6.14
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.86	5.19-6.53
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.77	5.32-6.22
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.28	5.87-6.69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.97	4.16-5.78
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.63	4.35-4.91
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.86	4.7-5.02
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.56	4.42-4.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.16	4.68-7.64

**(FT4 )游离T4**

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	11.8	10.7-12.9
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	17.7	16.1-19.3
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	16	13.8-18.2
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	21.2	19.6-22.8
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	pmol/L	14.3	6.8-21.8
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.1	14.7-19.5
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	pmol/L	21.4	19.2-23.6
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.1	15.6-18.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	11.9	10.3-13.5
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	pmol/L	15.6	14.1-17.1
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	12.5	11.8-13.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	13.7	12.1-15.3

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								6 / 9

Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.3	12.4-20.2
<b>(hCG)人绒毛膜促性腺激素</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.68	2.94-4.42
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.66	1.71-5.61
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/L	4.87	0.35-9.39
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.35	3.44-5.26
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.24	2.41-4.07
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	2.95	2.82-3.08
<b>(HCY)同型半胱氨酸</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	μmol/L	12.5	12.2-12.8
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.4	10.2-16.6
<b>(HGH)人生长激素</b>			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	1.82	1.38-2.26
<b>(IgE)免疫球蛋白 IgE</b>			
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	28.4	27.9-28.9
Beckman Image 800 Roche \ Immunoturbidimetry (ITA)	IU/mL	32	24.3-39.7
<b>(IGF1)胰岛素生长因子1</b>			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	223	169-277
<b>(INS)胰岛素</b>			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	mIU/L	13	12.2-13.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.3	9.3-13.3
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	mIU/L	11.5	9.9-13.1
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	9.4-12.4
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA )	mIU/L	10.1	8.9-11.3
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	mIU/L	10.8	10.3-11.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	mIU/L	12.2	7.1-17.3
<b>(LH)促黄体生成素</b>			
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.63	3.39-3.87
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.11	4.76-5.46
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/L	4.86	4.46-5.26
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	IU/L	4.56	4.2-4.92
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	6.18	5.91-6.45
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.1	2.48-3.72
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.56	3.1-4.02
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.42	3.87-4.97
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	IU/L	4.89	4.53-5.25
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.05	4.73-5.37
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.91	3.73-6.09
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	1.38	1.31-1.45
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.19	3.07-3.31
<b>(PRL)催/泌乳素</b>			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.9	6.62-9.18
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.17	5.64-6.7
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	10.3	9.1-11.5
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	9.09	8.56-9.62
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.58	8.55-8.61
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.75	4.83-8.67



<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								7 / 9

OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.65	7.7-9.6
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.39	7.85-8.93
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.46	6.62-8.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.49	6.45-10.53
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.02	6.19-7.85
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.94	4.79-5.09

**(PROG)孕酮**

Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.29	1.12-1.46
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.012	0.01-0.014
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.775	0.653-0.897
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.03	0.81-1.25
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.33	1.23-1.43
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.19	1-1.38
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.08	0.67-1.49
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.867	0.758-0.976
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.5	0.81-2.19
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.37	1.06-1.68
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.37	1.04-1.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.414	0.399-0.429
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.654	0.564-0.744

**(PSA)前列腺特异性抗原**

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.568	0.52-0.616
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.009	0-0.029
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.198	0-0.758
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.811	0.727-0.895
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.542	0.487-0.597
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.459	0.305-0.613
Maccura i3000 MAKER \ OTHER	ng/mL	0.638	0.601-0.675
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.636	0.563-0.709
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.592	0.476-0.708
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.589	0.452-0.726
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.922	0.821-1.023
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.486	0.361-0.611
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.49	0.444-0.536
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.47	0.47-0.47

**(PTH)甲状旁腺素**

Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.16	1.12-1.2
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	0.619	0.425-0.813
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.22	0.84-1.6

**(TES)睾酮**

Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.072	0.059-0.085
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.52	8.04-9
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.4	4.98-5.82
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.24	5.09-7.39
Roche cobas e411 Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.16	7.65-8.67
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.43	4.95-9.91

<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								8 / 9

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.22	3.92-8.52
OTHER Lifotronic \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.59	5.44-7.74
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.17	6.8-7.54
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.03	6.51-7.55
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.279	0.237-0.321
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.65	5.05-8.25
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.71	1.6-1.82

**(Tg)甲状腺球蛋白**

YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.07	3.67-4.47
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.7	1.7-37.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	7.85	7.25-8.45
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	IU/mL	2.76	0.91-4.61
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.784	0.654-0.914
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.09	2.32-3.86

**(TSH)促甲状腺素**

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.96	2.78-3.14
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	mIU/L	3.04	2.82-3.26
Mindray CL-2000i/1000i/6000i series Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.2	2.95-3.45
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	mIU/L	3.42	3.03-3.81
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	mIU/L	3.01	2.34-3.68
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.12	2.86-3.38
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	mIU/L	3.08	2.84-3.32
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.01	2.85-3.17
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	2.54	2.27-2.81
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	mIU/L	3.05	2.61-3.49
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	2.41	2.22-2.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.23	2.45-4.01
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	0.434	0.414-0.454
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	2.71	2.71-2.71

**(TT3)总三碘甲状腺原氨酸**

Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.68	1.4-1.96
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.25	1.08-1.42
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.41	0.85-1.97
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.36	1.18-1.54
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.43	1.22-1.64
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.17	1.05-1.29
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.5	1.26-1.74
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.18	0.77-1.59
chemclin LiCA Series CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.37	1.23-1.51
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.96	0.891-1.029
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.43	1.09-1.77
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.05	0.92-1.18

**(TT4)总T4**

maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	71.2	63.4-79
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	65.1	61.6-68.6
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	66.9	56.9-76.9
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	68.5	61.9-75.1



<b>REF</b>	539-1	6x3mL	<b>EXP</b>	2026-04-10	<b>LOT</b>	336A2F0	更新日期:	2025-02
								9 / 9

OTHER								
DIRUI \ Chemiluminescence Immunoassay (CLIA )	ng/mL		66.4				62.4-70.4	
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		66.3				59.5-73.1	
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL		69.2				46.8-91.6	
chemclin LiCA Series								
CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL		58.5				52.9-64.1	
Siemens ADVIA Centaur XP/XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL		57.5				55-60	
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL		62.1				47.2-77	
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL		65.6				61.3-69.9	
Siemens ADVIA Chemistry XPT								
Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL		89.9				89.9-89.9	

**(VB12)维生素B12**

Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL		366				330-402	
Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	pg/mL		441				422-460	
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL		436				331-541	
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	pg/mL		345				330-360	
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	pg/mL		435				399-471	

**(β2M)β-2-微球蛋白**

Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	ng/mL		728				639-817	
Roche cobas 4000/6000/8000								
Roche \ Immunoturbidimetry (ITA)	ng/mL		890				676-1104	
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL		1286				1222-1350	
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL		816				731-901	

**(β-hCG)人绒毛膜促性腺激素β亚基**

Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL		0.004				0.003-0.005	
Roche cobas e411								
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL		0.003				0.003-0.003	
Mindray CL-2000i/1000i/6000i series								
Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL		0.002				0.001-0.003	
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL		0.003				0.003-0.003	
Lifotronic ecl Series								
Lifotronic \ Chemiluminescence Immunoassay (CLIA )	IU/mL		0.003				0.003-0.003	
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL		0.003				0.002-0.004	

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。