

REF	600-2	30×1.7mL	EXP	2026-07-04	LOT	393DDBF	更新日期：	2024-05
1 / 2								

批号：393DDBF					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Ca)钙					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	1.09	0.8%	11	261
Roche cobas b 121/123/221 Dedicated \ ISE (indirect)	mmol/L	1.13	12%	<5	<100
Radiometer ABL Series Dedicated \ ISE (indirect)	mmol/L	1.1	12%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	1.15	12%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (direct)	mmol/L	1.07	0%	<5	<100
(Cl)氯					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	90.4	0.6%	12	307
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	101	12%	<5	<100
IL Gem Premier 3000 series Dedicated \ ISE (direct)	mmol/L	95.2	12%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	93.2	12%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (direct)	mmol/L	91	0%	<5	<100
(K)钾					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	4.1	0%	5	<100
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	4.15	12%	<5	<100
IL Gem Premier 3000 series Dedicated \ ISE (direct)	mmol/L	4.07	12%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	4.13	12%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (direct)	mmol/L	4	0%	<5	<100
(Na)钠					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	134	0.6%	12	306
Radiometer AQTSeries Dedicated \ ISE (direct)	mmol/L	140	12%	<5	<100
IL Gem Premier 3000 series Dedicated \ ICP-MS	mmol/L	146	12%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	142	12%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ ISE (indirect)	mmol/L	143	12%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (direct)	mmol/L	133	0%	<5	<100
(PCO2)二氧化碳分压					
Roche cobas b 121/123/221 Dedicated \ Electrode	mmHg	40	11.8%	<5	<100
Roche cobas b 121/123/221 Roche \ Electrode	mmHg	40.7	3.1%	11	291
Radiometer ABL Series Radiometer \ Electrode	mmHg	39.1	1.5%	11	251
Nova Stat Profile pHox series Dedicated \ Electrode	mmHg	32.9	12%	<5	<100
Radiometer ABL Series Dedicated \ Electrode	mmHg	41	12%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ Glucose oxidase (oxygen electrode)	mmHg	32.1	12%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	36.3	12%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	mmHg	43	12%	<5	<100
(pH)酸碱度					
Roche cobas b 121/123/221 Dedicated \ Electrode	No Units	7.5	11.8%	<5	<100
Roche cobas b 121/123/221 Roche \ Electrode	No Units	7.38	0.5%	11	294
Radiometer ABL Series Radiometer \ Electrode	No Units	7.4	0.1%	11	261
IL Gem Premier 3000 series Dedicated \ Electrode	No Units	7.53	12%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ Electrode	No Units	7.49	12%	<5	<100
Radiometer ABL Series Dedicated \ Electrode	No Units	7.56	12%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	No Units	7.45	12%	<5	<100
(PO2)氧分压					
Roche cobas b 121/123/221 Dedicated \ Electrode	mmHg	109	11.8%	<5	<100
Radiometer ABL Series Radiometer \ Electrode	mmHg	97.7	2.7%	12	253
Roche cobas b 121/123/221 Roche \ Electrode	mmHg	103	3.7%	11	292
Radiometer AQTSeries Dedicated \ Electrode	mmHg	95.5	12%	<5	<100

REF	600-2	30×1.7mL	EXP	2026-07-04	LOT	393DDBF	更新日期：	2024-05
2 / 2								

		批号：393DDBF			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Radiometer ABL Series Dedicated \ Electrode	mmHg	120	12%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	97.2	12%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ Electrode	mmHg	123	12%	<5	<100
Nova Stat Profile pHox series Dedicated \ Electrode	mmHg	124	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。