

REF	600-2	30×1.7mL	EXP	2025-12-17	LOT	393DD5F	更新日期：	2024-05
1 / 1								

批号：393DD5F					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Ca)钙					
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	1.14	10.5%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (indirect)	mmol/L	1.12	9.8%	<5	<100
Radiometer ABL Series Dedicated \ ISE (indirect)	mmol/L	1.09	10.1%	<5	<100
(Cl)氯					
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	99.6	9.1%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	92.3	9.9%	<5	<100
IL Gem Premier 3000 series Dedicated \ ISE (direct)	mmol/L	94.3	9.7%	<5	<100
(K)钾					
IL Gem Premier 3000 series Dedicated \ ISE (direct)	mmol/L	4.03	9.7%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	4.11	9.5%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	4.09	9.3%	<5	<100
(Na)钠					
Radiometer AQTSeries Dedicated \ ISE (direct)	mmol/L	139	8.6%	<5	<100
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	141	8.5%	<5	<100
IL Gem Premier 3000 series Dedicated \ ICP-MS	mmol/L	145	8.3%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ ISE (indirect)	mmol/L	142	7.7%	<5	<100
(PCO2)二氧化碳分压					
Radiometer ABL Series Dedicated \ Electrode	mmHg	40.6	9.6%	<5	<100
Nova Stat Profile pHx series Dedicated \ Electrode	mmHg	32.6	9%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	mmHg	42.6	9.2%	<5	<100
Roche cobas b 121/123/221 Dedicated \ Electrode	mmHg	39.6	9.3%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ Glucose oxidase (oxygen electrode)	mmHg	31.8	10.2%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	35.9	8.2%	<5	<100
(pH)酸碱度					
Siemens RAPIDPoint 400/500 Siemens \ Electrode	No Units	7.42	8.6%	<5	<100
Roche cobas b 121/123/221 Dedicated \ Electrode	No Units	7.43	9.6%	<5	<100
Radiometer ABL Series Dedicated \ Electrode	No Units	7.49	9.1%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	No Units	7.38	9.3%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	No Units	7.46	9.2%	<5	<100
(PO2)氧分压					
Nova Stat Profile pHx series Dedicated \ Electrode	mmHg	123	8%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	mmHg	94.6	9.4%	<5	<100
Radiometer ABL Series Dedicated \ Electrode	mmHg	119	9.2%	<5	<100
Roche cobas b 121/123/221 Dedicated \ Electrode	mmHg	108	8.3%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ Electrode	mmHg	122	8.6%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	96.2	8.9%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。