

质控品名： 特殊生化与蛋白非定值质控品

|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 583-1 | 6×3mL | EXP | 2026-02-19 | LOT | 379C66F | 更新日期： | 2024-05 |
| 1 / 2 |       |       |     |            |     |         |       |         |

| 批号：379C66F                                                           |        |         |       |          |      |
|----------------------------------------------------------------------|--------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                          | 单位     | 均值 \ CV |       | Lab \ 点数 |      |
| (5'NT)5`核苷酸酶                                                         |        |         |       |          |      |
| Beckman AU series<br>Beckman Coulter \ Enzymatic                     | U/L    | 16      | 7.8%  | <5       | <100 |
| Beckman AU series<br>BSBE \ Enzymatic                                | U/L    | 59.5    | 7.3%  | <5       | <100 |
| MedicalSystem MS-i3080<br>MedicalSystem \ Enzymatic                  | U/L    | 32.1    | 2.1%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>chongqiong zybio \ Enzymatic  | U/L    | 37.2    | 4.4%  | <5       | <100 |
| (ADA)腺苷脱氨酶                                                           |        |         |       |          |      |
| Roche Hitachi/Hitachi Systems<br>MAKER \ Trinder reaction            | U/L    | 20.7    | 5.6%  | 9        | 220  |
| Roche Hitachi/Hitachi Systems<br>MAKER \ Peroxidase                  | U/L    | 21.2    | 2.7%  | 7        | 204  |
| Beckman AU series<br>Whitman \ Peroxidase                            | U/L    | 23.7    | 2.3%  | 6        | <100 |
| Beckman AU series<br>BSBE \ Peroxidase                               | U/L    | 20.1    | 9.9%  | 5        | 130  |
| Beckman AU series<br>Beckman Coulter \ Peroxidase                    | U/L    | 22.8    | 9.8%  | <5       | <100 |
| Roche Hitachi/Hitachi Systems<br>MAKER \ OTHER                       | U/L    | 21.4    | 9.6%  | <5       | <100 |
| Siemens ADVIA Chemistry XPT<br>KHB \ Trinder reaction                | U/L    | 28.4    | 2.5%  | <5       | <100 |
| Mindray BC-5100/BC-5300 Series<br>LEADMAN \ Peroxidase               | U/L    | 24      | 3%    | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Peroxidase                 | U/L    | 20.4    | 0.7%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>chongqiong zybio \ Peroxidase | U/L    | 21.4    | 3.9%  | <5       | <100 |
| MedicalSystem MS-i3080<br>MedicalSystem \ Peroxidase                 | U/L    | 21      | 4.2%  | <5       | <100 |
| Mindray BS series<br>OTHER \ Peroxidase                              | U/L    | 20.2    | 1.6%  | <5       | <100 |
| (AFU)α-L-岩藻糖苷酶                                                       |        |         |       |          |      |
| Beckman AU series<br>BSBE \ CNPF                                     | U/L    | 9.97    | 5.5%  | 6        | 155  |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ CNPF                       | U/L    | 11.4    | 0.7%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>chongqiong zybio \ CNPF       | U/L    | 4.92    | 5.9%  | <5       | <100 |
| Roche cobas c 501/502/701/702<br>Roche \ CNPF                        | U/L    | 12.9    | 2.7%  | <5       | <100 |
| MedicalSystem MS-i3080<br>MedicalSystem \ CNPF                       | U/L    | 12      | 6.6%  | <5       | <100 |
| Mindray BC-5100/BC-5300 Series<br>LEADMAN \ CNPF                     | U/L    | 6.73    | 3.9%  | <5       | <100 |
| Beckman AU series<br>Beckman Coulter \ CNPF                          | U/L    | 13.6    | 9.6%  | <5       | <100 |
| Mindray BS series<br>OTHER \ OTHER                                   | U/L    | 11.6    | 2.6%  | <5       | <100 |
| (anti CCP)抗环瓜氨酸抗体                                                    |        |         |       |          |      |
| Beckman AU series<br>OTHER \ immunoturbidimetric                     | U/mL   | 14      | 18.8% | <5       | <100 |
| Beckman AU series<br>MedicalSystem \ immunoturbidimetric             | U/mL   | 59.9    | 1.8%  | <5       | <100 |
| Roche cobas e411<br>Roche \ Chemiluminescence Immunoassay (CLIA )    | U/mL   | 14.5    | 9.7%  | <5       | <100 |
| Beckman AU series<br>Beckman Coulter \ immunoturbidimetric           | U/mL   | 13.2    | 2.8%  | <5       | <100 |
| (ASO)抗链球菌溶血素 O 抗体                                                    |        |         |       |          |      |
| Beckman AU series<br>Beckman Coulter \ Immunonephelometry (INA)      | IU/mL  | 68.6    | 0.7%  | <5       | <100 |
| (CO2 (HCO3-) )二氧化碳（碳酸氢根）                                             |        |         |       |          |      |
| Roche Hitachi/Hitachi Systems<br>MAKER \ OTHER                       | mmol/L | 8.9     | 9.4%  | 7        | 209  |
| Roche Hitachi/Hitachi Systems<br>MAKER \ Colorimetry                 | mmol/L | 8.17    | 7.8%  | 11       | 238  |
| Siemens ADVIA Chemistry XPT<br>DiaSys \ Electrode                    | mmol/L | 7.7     | 23.8% | <5       | <100 |
| Beckman AU series<br>Beckman Coulter \ PEPC                          | mmol/L | 14.3    | 9.8%  | <5       | <100 |
| Roche cobas c 501/502/701/702<br>Roche \ PEPC                        | mmol/L | 8.78    | 19.3% | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ PEPC                       | mmol/L | 9.7     | 5%    | <5       | <100 |
| Mindray BS series<br>MedicalSystem \ OTHER                           | mmol/L | 8.27    | 4.4%  | <5       | <100 |
| (Cys C)半胱氨酸蛋白酶抑制剂C/ 胱抑素C                                             |        |         |       |          |      |
| Beckman AU series<br>OTHER \ OTHER                                   | mg/L   | 0.305   | 4.9%  | 6        | <100 |
| Beckman AU series<br>chongqiong zybio \ Immunoturbidimetry (ITA)     | mg/L   | 0.506   | 7.1%  | 6        | 153  |
| Siemens ADVIA Chemistry XPT<br>Dedicated \ Immunonephelometry (INA)  | mg/L   | 0.519   | 5.1%  | <5       | <100 |
| Roche cobas c 501/502/701/702<br>Roche \ Immunonephelometry (INA)    | mg/L   | 0.568   | 1.9%  | <5       | <100 |

质控品名：  特殊生化与蛋白非定值质控品

|                |       |       |                |            |                |         |       |         |
|----------------|-------|-------|----------------|------------|----------------|---------|-------|---------|
| <div>REF</div> | 583-1 | 6×3mL | <div>EXP</div> | 2026-02-19 | <div>LOT</div> | 379C66F | 更新日期： | 2024-05 |
| 2 / 2          |       |       |                |            |                |         |       |         |

|                                                         |        |            |      |          |      |
|---------------------------------------------------------|--------|------------|------|----------|------|
|                                                         |        | 批号：379C66F |      |          |      |
| 项目\仪器\试剂\方法                                             | 单位     | 均值 \ CV    |      | Lab \ 点数 |      |
| MedicalSystem MS-i3080                                  | mg/L   | 0.448      | 2.6% | <5       | <100 |
| MedicalSystem \ Chemiluminescence Immunoassay (CLIA )   |        |            |      |          |      |
| Mindray BC-5100/BC-5300 Series                          | mg/L   | 0.602      | 7.4% | <5       | <100 |
| LEADMAN \ Immunoturbidimetry (ITA)                      |        |            |      |          |      |
| Beckman AU series                                       | mg/L   | 0.54       | 9.3% | <5       | <100 |
| Beckman Coulter \ Chemiluminescence Immunoassay (CLIA ) |        |            |      |          |      |
| Mindray BS series                                       | mg/L   | 0.496      | 2.6% | <5       | <100 |
| MedicalSystem \ OTHER                                   |        |            |      |          |      |
| Siemens ADVIA Centaur XP/XPT                            | mg/L   | 0.384      | 4.5% | <5       | <100 |
| Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)  |        |            |      |          |      |
| Roche Hitachi/Hitachi Systems                           | mg/L   | 0.436      | 4.4% | <5       | <100 |
| DiaSys \ Immunonephelometry (INA)                       |        |            |      |          |      |
| (FRUC (GSP) )果糖胺（糖化血清蛋白）                                |        |            |      |          |      |
| Roche cobas e411                                        | mmol/L | 0.21       | 10%  | <5       | <100 |
| Roche \ NBT                                             |        |            |      |          |      |
| Beckman AU series                                       | mmol/L | 1.52       | 1%   | <5       | <100 |
| Beckman Coulter \ NBT                                   |        |            |      |          |      |
| Siemens ADVIA Centaur XP/XPT                            | mmol/L | 0.113      | 4.7% | <5       | <100 |
| Siemens \ NBT                                           |        |            |      |          |      |
| Roche Hitachi/Hitachi Systems                           | mmol/L | 1.33       | 1.9% | <5       | <100 |
| MedicalSystem \ NBT                                     |        |            |      |          |      |
| (HCY)同型半胱氨酸                                             |        |            |      |          |      |
| Beckman AU series                                       | μmol/L | 4.48       | 2.2% | <5       | <100 |
| Beckman Coulter \ Immunonephelometry (INA)              |        |            |      |          |      |
| (Lp(a))脂蛋白a                                             |        |            |      |          |      |
| Beckman AU series                                       | mg/L   | 132        | 9.8% | <5       | <100 |
| Beckman Coulter \ Immunoturbidimetry (ITA)              |        |            |      |          |      |
| Beckman AU series                                       | mg/L   | 137        | 1.2% | <5       | <100 |
| Beckman Coulter \ Immunonephelometry (INA)              |        |            |      |          |      |
| (NEFA)游离脂肪酸                                             |        |            |      |          |      |
| Beckman AU series                                       | mmol/L | 1.58       | 1.2% | <5       | <100 |
| Beckman Coulter \ ACS - ACOD                            |        |            |      |          |      |
| (RF)类风湿因子                                               |        |            |      |          |      |
| Beckman AU series                                       | IU/mL  | 7.5        | 2.8% | <5       | <100 |
| Beckman Coulter \ Immunonephelometry (INA)              |        |            |      |          |      |
| (β-Hb)Beta-羟丁酸                                          |        |            |      |          |      |
| Siemens ADVIA Chemistry XPT                             | mmol/L | 0.216      | 4.7% | <5       | <100 |
| KHB \ Dehydrogenase                                     |        |            |      |          |      |
| Beckman AU series                                       | mmol/L | 0.26       | 7.7% | <5       | <100 |
| Beckman Coulter \ Dehydrogenase                         |        |            |      |          |      |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。