

REF	600-1	30×1.7mL	EXP	2026-07-12	LOT	393B6AF	更新日期：	2024-05
1 / 1								

批号：393B6AF					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Ca)钙					
OTHER OTHER \ ISE (direct)	mmol/L	1.58	2.2%	13	310
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	1.52	0.9%	12	268
Radiometer ABL Series Dedicated \ ISE (indirect)	mmol/L	1.56	12%	<5	<100
(Cl)氯					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	76	0.1%	12	309
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	76.8	12%	<5	<100
(K)钾					
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	2.3	0%	6	105
OTHER OTHER \ ISE (direct)	mmol/L	2.4	1%	12	262
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	2.28	12%	<5	<100
(Na)钠					
OTHER OTHER \ ISE (direct)	mmol/L	113	0.6%	13	320
Radiometer ABL Series Radiometer \ ISE (direct)	mmol/L	109	0.4%	12	308
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	109	12%	<5	<100
(PCO2)二氧化碳分压					
OTHER OTHER \ Electrode	mmHg	77.5	6.7%	13	306
Radiometer ABL Series Radiometer \ Electrode	mmHg	63.6	2.1%	12	259
Radiometer ABL Series Dedicated \ Electrode	mmHg	67.7	12%	<5	<100
(pH)酸碱度					
OTHER OTHER \ Electrode	No Units	7.14	0.2%	13	320
Radiometer ABL Series Radiometer \ Electrode	No Units	7.18	0.1%	12	267
Radiometer ABL Series Dedicated \ Electrode	No Units	7.22	12%	<5	<100
(PO2)氧分压					
OTHER OTHER \ Electrode	mmHg	87.1	3.8%	13	274
Radiometer ABL Series Radiometer \ Electrode	mmHg	60.5	4.4%	12	257
Radiometer ABL Series Dedicated \ Electrode	mmHg	65.9	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。