

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期: 2025-03
							1 / 14

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(CK)肌酸激酶	(CK-MB)肌酸激酶-MB (活性)	(Myo / Mb)肌红蛋白
(TnI)肌钙蛋白I	(TnT)肌钙蛋白T	(BNP)B型利钠肽
(NT-proBNP)N-末端脑钠肽前体	(DGXN)地高辛	(TT3)总三碘甲状腺原氨酸
(FT3)游离T3	(TT4)总T4	(FT4 )游离T4
(TSH)促甲状腺素	(Cortisol)皮质醇	(FSH)促卵泡刺激素
(LH)促黄体生成素	(PROG)孕酮	(PRL)催/泌乳素
(TES)睾酮	(E2)雌二醇	(C-P)C-肽
(FA)叶酸	(INS)胰岛素	(VB12)维生素B12
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(CA 125)糖类抗原 125
(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9	(β2M)β-2-微球蛋白
(Ferritin)铁蛋白	(β-hCG)人绒毛膜促性腺激素β亚基	(F-PSA)游离前列腺特异性抗原
(HCY)同型半胱氨酸	(NSE)神经元特异性烯醇化酶	(Cyfra 21-1)细胞角蛋白19片段
(SCC)鳞状细胞癌抗原	(CA 50)糖类抗原 50	(CA 72-4)糖类抗原 72-4
(ASD)雄烯二酮	(ALD)醛固酮	(Anti-Tg)抗甲状腺球蛋白抗体
(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮	(EPO)促红细胞生成素
(HGH)人生长激素	(IGF1)胰岛素生长因子1	(PTH)甲状旁腺素
(SHBG)性激素结合球蛋白	(Tg)甲状腺球蛋白	(25-OH-Vit D)25羟基维生素D
(ACTH)促肾上腺皮质激素	(Renin)肾素	(CT)降钙素
(IL6)白介素6	(HE4)人附睾蛋白4	(AMH)抗缪勒氏管激素
(25-OH-VD2)25羟基维生素D2	(IntactPTH)全段甲状旁腺激素	(Digitoxin)洋地黄毒苷
(HS-CTnI) (高) 超敏肌钙蛋白I	(CK-MB mass)肌酸激酶-MB质量	

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								2 / 14

批号: 3382931			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
<b>(25-OH-VD2)25羟基维生素D2</b>			
Autolumo A6000 Series	ng/mL	23.8	20.6-27
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	23.8	20.6-27
Autolumo A6000 Series	nmol/L	57.6	49.9-65.3
Autobio \ Chemiluminescence Immunoassay (CLIA )	nmol/L	57.6	49.9-65.3
<b>(25-OH-Vit D)25羟基维生素D</b>			
Roche cobas 4000/6000/8000	ng/mL	26	22.2-29.8
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	26	22.2-29.8
Roche cobas 4000/6000/8000	nmol/L	64.9	55.4-74.4
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	64.9	55.4-74.4
OTHER	ng/mL	25.7	15.6-35.8
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	25.7	15.6-35.8
OTHER	nmol/L	64.2	39-89.4
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	64.2	39-89.4
OTHER	ng/mL	39.6	34.7-44.5
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	39.6	34.7-44.5
OTHER	nmol/L	99	86.7-111.3
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	99	86.7-111.3
maglumi system	ng/mL	26.8	14.9-38.7
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	26.8	14.9-38.7
maglumi system	nmol/L	66.9	37.2-96.6
Snibe \ Chemiluminescence Immunoassay (CLIA )	nmol/L	66.9	37.2-96.6
Maccura i3000	ng/mL	40.2	33.8-46.6
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	40.2	33.8-46.6
Maccura i3000	nmol/L	101	85-117
MAKER \ Chemiluminescence Immunoassay (CLIA )	nmol/L	101	85-117
Abbott ARCHITECT i series	ng/mL	24.6	21.5-27.7
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	24.6	21.5-27.7
Abbott ARCHITECT i series	nmol/L	61.5	53.8-69.2
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	61.5	53.8-69.2
Roche cobas E601/602	ng/mL	24.3	19.4-29.2
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.3	19.4-29.2
Roche cobas E601/602	nmol/L	60.6	48.5-72.7
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	60.6	48.5-72.7
AB Sciex API Series	ng/mL	40	39.7-40.3
ABX \ Chemiluminescence Immunoassay (CLIA )	ng/mL	40	39.7-40.3
AB Sciex API Series	nmol/L	99.9	99.1-100.7
ABX \ Chemiluminescence Immunoassay (CLIA )	nmol/L	99.9	99.1-100.7
<b>(ACTH)促肾上腺皮质激素</b>			
Mindray CL-2000i/1000i/6000i series	pg/mL	1	1-1
Mindray \ Chemiluminescence Immunoassay (CLIA )	pg/mL	1	1-1
<b>(AFP)甲胎蛋白</b>			
Autolumo A6000 Series	ng/mL	8.38	7.46-9.3
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.38	7.46-9.3
Autolumo A2000/A2000 PLUS	ng/mL	8.24	5.32-11.16
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	8.24	5.32-11.16
OTHER	ng/mL	6.25	4.16-8.34
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.25	4.16-8.34
OTHER	ng/mL	7.5	6.24-8.76
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.5	6.24-8.76
maglumi system	ng/mL	6.04	3-9.08
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.04	3-9.08
maglumi system	ng/mL	9.43	8.3-10.56
Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.43	8.3-10.56
Abbott ARCHITECT i series	ng/mL	6.02	5.62-6.42
Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.02	5.62-6.42
Roche cobas E601/602	ng/mL	6.88	5.6-8.16
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.88	5.6-8.16
Beckman UniCel Dxl 600/800	ng/mL	6.75	5.79-7.71
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.75	5.79-7.71
Abbott ARCHITECT i series	ng/mL	6.59	6.22-6.96
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.59	6.22-6.96
Beckman UniCel Dxl 600/800	ng/mL	6.31	5.51-7.11
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.31	5.51-7.11
Maccura i3000	ng/mL	7.19	6.21-8.17
MAKER \ OTHER	ng/mL	7.19	6.21-8.17
Roche cobas 4000/6000/8000	ng/mL	7.07	4.93-9.21
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.07	4.93-9.21
Abbott ARCHITECT ci series	ng/mL	6.13	6.06-6.2
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.13	6.06-6.2
Mindray CL-2000i/1000i/6000i series	ng/mL	5.95	5.13-6.77
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.95	5.13-6.77
Maccura i3000	ng/mL	7.1	5.68-8.52
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.1	5.68-8.52
Mindray BS series	ng/mL	7.46	7.46-7.46
Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.46	7.46-7.46
Roche cobas e411	ng/mL	6.68	6.68-6.68
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.68	6.68-6.68
Systemx HISCL-5000	ng/mL	6.9	6.9-6.9
systemx \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.9	6.9-6.9

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								3 / 14

<b>(ALD)醛固酮</b>			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	pg/mL	100	0-233.6
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pg/mL	46.5	29.8-63.2
<b>(AMH)抗缪勒氏管激素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.34	4.89-5.79
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.25	4.62-5.88
<b>(Anti-Tg)抗甲状腺球蛋白抗体</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	0.791	0.696-0.886
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	19.4	16-22.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	1.31	1.14-1.48
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	8.14	6.89-9.39
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	9.54	5.93-13.15
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	16.2	11.8-20.6
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	IU/mL	1.32	1.16-1.48
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.5	11.8-45.2
Maccura i3000 MAKER \ OTHER	IU/mL	11.7	6.5-16.9
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	4.71	3.53-5.89
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14.3	14.3-14.3
<b>(Anti-TPO)抗甲状腺过氧化物酶抗体</b>			
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	IU/mL	11.6	9.9-13.3
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15	10-20
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	7.2	6.02-8.38
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	8.49	7.4-9.58
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	12.2	7.8-16.6
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	5.79	4.82-6.76
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/mL	19.6	16.5-22.7
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	29.1	6.3-51.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14.5	11.8-17.2
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	IU/mL	7.22	6.38-8.06
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	5.08	4.17-5.99
Maccura i3000 MAKER \ OTHER	IU/mL	6.09	5.15-7.03
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	11.9	11.9-11.9
<b>(ASD)雄烯二酮</b>			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.92	1.72-2.12
<b>(BNP)B型利钠肽</b>			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	291	221-361
<b>(CA 125)糖类抗原 125</b>			
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	36.3	34.6-38
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	30.9	21.3-40.5
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	27.8	20.7-34.9
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	32.6	29.5-35.7
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	38	34.4-41.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	23-33.4
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/mL	33.5	30.3-36.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	30.3	27.4-33.2

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								4 / 14

Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	39.3	36.4-42.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	42.8	33.7-51.9
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	25.6	24.1-27.1
Maccura i3000 MAKER \ OTHER	IU/mL	43.1	40.5-45.7
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	34.9	34.2-35.6
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	43.7	39-48.4
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	44.1	35.3-52.9
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA )	IU/mL	7.47	7.47-7.47

**(CA 15-3)糖类抗原 15-3**

OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	43.1	16.8-69.4
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	29.1	26.5-31.7
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	29.3	24.2-34.4
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	24.5	18.4-30.6
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/mL	28.7	25.4-32
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	25.9	21.8-30
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/mL	29.6	26.4-32.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	15.3	13.7-16.9
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	29.1	26-32.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	16.3	14.8-17.8
Maccura i3000 MAKER \ OTHER	IU/mL	23.7	21.6-25.8
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	29.2	23.5-34.9
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	25.8	25.3-26.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	16.2	12.4-20
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	24.4	19.5-29.3
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA )	IU/mL	2.44	2.44-2.44

**(CA 19-9)糖类抗原 19-9**

Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	26.1	24.7-27.5
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	27.3	20.8-33.8
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	22.3	18.1-26.5
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	22.4	18.7-26.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	25.2	22-28.4
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/mL	29.6	27.4-31.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	42.8	37.7-47.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	22.9	10.8-35
Maccura i3000 MAKER \ OTHER	IU/mL	17	15.5-18.5
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	25.1	20.1-30.1
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	21.4	18.7-24.1
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	IU/mL	17.3	13.8-20.8
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA )	IU/mL	12.5	12.5-12.5

**(CA 50)糖类抗原 50**

Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/mL	18.4	15.7-21.1
--	-------	------	-----------

**(CA 72-4)糖类抗原 72-4**

Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	5.44	4.46-6.42
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	5.44	4.94-5.94

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								5 / 14

(CEA)癌胚抗原			
Autolumo A6000 Series			
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.09	4.69-5.49
Mindray CL-2000i/1000i/6000i series			
Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.99	1.7-2.28
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.92	2.88-4.96
OTHER			
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.46	1.48-3.44
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.66	2.42-2.9
Mindray CL-2000i/1000i/6000i series			
Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.14	1.88-2.4
maglumi system			
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.25	2.49-4.01
maglumi system			
Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.66	2.12-3.2
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.89	2.2-3.58
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.51	1.93-3.09
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.3	2.05-2.55
Abbott ARCHITECT i series			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.88	2.57-3.19
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.87	1.51-2.23
Maccura i3000			
MAKER \ OTHER	ng/mL	2.69	2.37-3.01
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.48	1.52-3.44
Abbott ARCHITECT ci series			
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.93	2.84-3.02
Maccura i3000			
MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.69	2.15-3.23
Mindray BS series			
Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.45	2.45-2.45
Roche cobas e411			
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.38	2.38-2.38
(CK)肌酸激酶			
Roche cobas 4000/6000/8000			
Roche \ Urease-UV-rate	U/L	126	101-151
(CK-MB mass)肌酸激酶-MB质量			
Beckman UniCel DxC 600/800			
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.94	3.66-4.22
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.69	1.26-2.12
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.19	3.51-4.87
(CK-MB)肌酸激酶-MB (活性)			
Ortho Vitros 3600/5600			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.97	2.5-3.44
Roche cobas E601/602			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.48	4-4.96
Roche cobas 4000/6000/8000			
Roche \ Immunoturbidimetry (ITA)	U/L	6.9	5.52-8.28
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.13	3.79-4.47
Roche cobas e411			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.44	3.55-5.33
(Cortisol)皮质醇			
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	87.9	77.7-98.1
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	108	95-121
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	109	100-118
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	116	88-144
Autolumo A6000 Series			
Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	107	95-119
(C-P)C-肽			
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.04	0.97-1.11
Roche cobas 4000/6000/8000			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.03	0.95-1.11
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.04	0.89-1.19
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.971	0.651-1.291
OTHER			
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.24	0.8-1.68

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								6 / 14

Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA) maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.87	0.807-0.933
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.01	0.92-1.1
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.585	0.495-0.675
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.841	0.755-0.927
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.643	0.566-0.72
<b>(CT)降钙素</b>			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	1.71	1.06-2.36
<b>(Cyfra 21-1)细胞角蛋白19片段</b>			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.23	2.01-2.45
Maccura i3000 MAKER \ OTHER	ng/mL	2.18	1.83-2.53
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.22	1.78-2.66
<b>(DGXN)地高辛</b>			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.892	0.812-0.972
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.09	0.87-1.31
<b>(DHEA-S)硫酸脱氢表雄酮</b>			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2536	2323-2749
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3651	2921-4381
<b>(Digitoxin)洋地黄毒苷</b>			
Roche cobas 4000/6000/8000 Dedicated \ Chemiluminescence Immunoassay (CLIA )	nmol/L	1.17	1.09-1.25
<b>(E2)雌二醇</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	254	235-273
Autolumo A6000 Series Autobio \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	345	285-405
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	264	231-297
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	pmol/L	442	21-863
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	228	204-252
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	292	246-338
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	253	224-282
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	258	212-304
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	214	166-262
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	241	210-272
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	74.2	60.4-88
<b>(EPO)促红细胞生成素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	6.65	5.87-7.43
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/mL	5.73	4.58-6.88
<b>(FA)叶酸</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.32	3.74-4.9
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.2	2.59-3.81
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.06	3.52-4.6
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.09	3.39-4.79
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.89	3.03-4.75
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.3	3.6-5
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.85	3.08-4.62
<b>(Ferritin)铁蛋白</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	58.5	46.6-70.4



<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								7 / 14

Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	104	87-121
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	62	50.5-73.5
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	77.6	67.8-87.4
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	91.9	73.3-110.5
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	100	37-163
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	78.3	60.9-95.7
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	93.1	53.6-132.6
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	86.3	71.6-101
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	90.6	66.1-115.1
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	97.3	81.3-113.3
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	158	126-190
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	120	96-144
<b>(F-PSA)游离前列腺特异性抗原</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.69	1.57-1.81
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.41	0-3.32
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.45	2.18-2.72
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.91	1.26-2.56
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.22	2.04-2.4
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.7	1.4-2
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.67	1.52-1.82
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.25	2.05-2.45
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.07	1.91-2.23
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.16	2-2.32
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.19	2.12-2.26
<b>(FSH)促卵泡刺激素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.46	4.77-6.15
Autolumo A6000 Series Autobio \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	6.77	5.97-7.57
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/L	7	5.71-8.29
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.64	5.02-6.26
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.55	3.15-7.95
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.8	4.53-5.07
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.12	4.67-5.57
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.13	4.66-5.6
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.21	4.73-5.69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.3	4.49-6.11
<b>(FT3)游离T3</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.84	5.36-6.32
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.61	5.19-6.03
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	4.8	4.14-5.46
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.49	4.84-6.14
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	6.32	2.34-10.3
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.99	5.13-6.85
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.41	4.07-4.75
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.7	4.14-5.26

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								8 / 14

Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.93	4.17-5.69
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	pmol/L	4.87	4.51-5.23
Mindray BS series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.58	4.93-6.23
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.63	4.85-6.41
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	4.88	4.12-5.64
Maccura i3000 MAKER \ OTHER	pmol/L	6.28	5.46-7.1
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	5.55	5.43-5.67
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	5.57	4.83-6.31
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA )	pmol/L	9.11	8.45-9.77
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	pmol/L	6	4.8-7.2

**(FT4)游离T4**

Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.9	15.8-18
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	pmol/L	12.4	10.2-14.6
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.7	15.3-18.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	11.7	10.1-13.3
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	18	16.7-19.3
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	73.1	0-219.9
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	53.9	47.6-60.2
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	pmol/L	16.6	13.3-19.9
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	12.5	11.4-13.6
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	13	11.7-14.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	12	10.4-13.6
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	pmol/L	15.5	14.8-16.2
Mindray BS series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	22	20.2-23.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.1	14-18.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	12.6	11-14.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	13.3	11.8-14.8
Maccura i3000 MAKER \ OTHER	pmol/L	17.4	15-19.8
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	13	12.6-13.4
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	21.4	19.3-23.5

**(hCG)人绒毛膜促性腺激素**

OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.33	3.13-5.53
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	6.73	5.21-8.25
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.77	2.93-4.61
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.81	2.42-7.2
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/L	5.76	4.61-6.91
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.7	3.37-4.03
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	IU/L	4.34	3.94-4.74

**(HCY)同型半胱氨酸**

Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	μmol/L	12.3	11.6-13
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.5	10.8-16.2
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA )	μmol/L	1.65	1.56-1.74

**(HE4)人附睾蛋白4**

YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	pmol/L	30.8	26.1-35.5
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	43.3	38.2-48.4



<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								9 / 14

(HGH)人生长激素			
maglumi system			
Snibe \ Chemiluminescence Immunoassay (CLIA )	µg/L	2.57	2.34-2.8
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	1.86	1.49-2.23
(HS-CTnI) (高) 超敏肌钙蛋白I			
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.071	0.064-0.078
(IgE)免疫球蛋白 IgE			
Beckman Image 800			
Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	142	114-170
(IGF1)胰岛素生长因子1			
maglumi system			
Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	178	162-194
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	221	177-265
(IL6)白介素6			
OTHER			
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	9.34	5.83-12.85
OTHER			
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	6.35	0-14.07
(INS)胰岛素			
Roche cobas 4000/6000/8000			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	7.94	7.3-8.58
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8	7.34-8.66
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	6.57	6.14-7
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	mIU/L	10.6	8.8-12.4
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	7.68	4.67-10.69
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.47	7.05-9.89
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	9.9	8.91-10.89
OTHER			
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.7	5.59-13.81
Chemclin LICA 800			
CHIVD \ Chemiluminescence Immunoassay (CLIA )	mIU/L	1.16	1.12-1.2
Mindray CL-2000i/1000i/6000i series			
Mindray \ Chemiluminescence Immunoassay (CLIA )	mIU/L	8.49	8.05-8.93
Maccura i3000			
MAKER \ Chemiluminescence Immunoassay (CLIA )	mIU/L	7.53	6.72-8.34
maglumi system			
Snibe \ Chemiluminescence Immunoassay (CLIA )	mIU/L	0.909	0.827-0.991
(IntactPTH)全段甲状旁腺激素			
Abbott ARCHITECT i series			
Abbott \ OTHER	pg/mL	3	2.54-3.46
(LH)促黄体生成素			
Beckman UniCel Dxl 600/800			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.66	3.33-3.99
Autolumo A6000 Series			
Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.97	4.67-5.27
Autolumo A2000/A2000 PLUS			
Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/L	6.12	4.64-7.6
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.81	3.47-4.15
Tosoh AIA Series			
OTHER \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.75	3.5-4
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.3	4.49-6.11
Abbott ARCHITECT i series			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.03	2.75-3.31
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	3.09	2.86-3.32
Siemens Immulite1000/2000			
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.81	4.2-5.42
Roche cobas 4000/6000/8000			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.39	4.74-6.04
Beckman UniCel Dxl 600/800			
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.62	2.96-4.28
(Myo / Mb)肌红蛋白			
Ortho Vitros 3600/5600			
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	µg/L	47.1	38.9-55.3
Roche cobas E601/602			
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	30.6	27.2-34
Autolumo A6000 Series			
Autobio \ Chemiluminescence Immunoassay (CLIA )	µg/L	23.2	20-26.4

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								10 / 14

OTHER								
OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)		µg/L		41.3				29.4-53.2
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )		µg/L		30.4				26.8-34
OTHER								
Norman \ Chemiluminescence Immunoassay (CLIA )		µg/L		38.7				32.3-45.1
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		µg/L		28.2				22.6-33.8
Ortho Vitros 3600/5600								
Ortho \ Chemiluminescence Immunoassay (CLIA )		µg/L		22.6				21.4-23.8
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		µg/L		34				31.6-36.4
<b>(NSE)神经元特异性烯醇化酶</b>								
Maccura i3000								
MAKER \ OTHER		ng/mL		4.68				3.37-5.99
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )		ng/mL		5.13				4.11-6.15
<b>(NT-proBNP)N-末端脑钠肽前体</b>								
Roche cobas e411								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)		pg/mL		190				175-205
Roche cobas e411								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)		pmol/L		22.3				20.5-24.1
Roche cobas E601/602								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)		pg/mL		183				171-195
Roche cobas E601/602								
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)		pmol/L		21.5				20.1-22.9
OTHER								
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)		pg/mL		185				127-243
OTHER								
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)		pmol/L		21.8				14.9-28.7
Zybio ZS 400								
chongqiong zybio \ Immunofluorescence (IFA)		pg/mL		235				186-284
Zybio ZS 400								
chongqiong zybio \ Immunofluorescence (IFA)		pmol/L		27.6				21.9-33.3
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		pg/mL		183				173-193
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		pmol/L		21.6				20.4-22.8
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		pg/mL		173				138-208
Roche cobas e411								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		pmol/L		20.4				16.3-24.5
<b>(PRL)催乳素</b>								
Beckman UniCel Dxl 600/800								
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)		ng/mL		6.81				6.27-7.35
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		7.54				5.23-9.85
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		8.35				7.77-8.93
Autolumo A2000/A2000 PLUS								
Autobio \ Chemiluminescence Immunoassay (CLIA )		ng/mL		9.44				7.95-10.93
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		9.48				7.32-11.64
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )		ng/mL		6.87				5.84-7.9
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)		ng/mL		8.48				7.22-9.74
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA )		ng/mL		8.34				7.86-8.82
Siemens Immulite1000/2000								
Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		6.55				5.84-7.26
<b>(PROG)孕酮</b>								
Beckman UniCel Dxl 600/800								
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)		ng/mL		1.44				1.13-1.75
Beckman UniCel Dxl 600/800								
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )		ng/mL		1.36				1.17-1.55
Autolumo A6000 Series								
Autobio \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		1.29				1.04-1.54
Autolumo A2000/A2000 PLUS								
Autobio \ Chemiluminescence Immunoassay (CLIA )		ng/mL		1.05				0.5-1.6
Roche cobas E601/602								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		1.13				0.92-1.34
Maccura i3000								
MAKER \ Chemiluminescence Immunoassay (CLIA )		ng/mL		1.23				1.1-1.36
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)		ng/mL		1.05				0.85-1.25
Abbott ARCHITECT i series								
Abbott \ Chemiluminescence Immunoassay (CLIA )		ng/mL		1.11				0.99-1.23
Siemens Immulite1000/2000								
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)		ng/mL		1.44				1.24-1.64
Roche cobas 4000/6000/8000								
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)		ng/mL		0.89				0-1.801

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								11 / 14

Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA) maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.31	0.82-1.8
ng/mL	1.35	1.09-1.61	
<b>(PSA)前列腺特异性抗原</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.38	4.08-4.68
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.07	3.59-4.55
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.45	4.33-4.57
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.96	4.2-5.72
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.36	5.05-5.67
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.2	3.6-4.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.65	3.96-5.34
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.08	3.71-4.45
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.01	4.57-5.45
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.97	3.61-4.33
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4	3.66-4.34
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.09	4.02-4.16
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4	3.2-4.8
<b>(PTH)甲状旁腺素</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	2.08	1.81-2.35
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	2.4	2-2.8
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	2.32	1.24-3.4
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pmol/L	1.5	0.01-2.99
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	2.96	2.72-3.2
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	pmol/L	1.91	1.59-2.23
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.79	1.41-2.17
Autolumo A6000 Series Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	2.15	1.82-2.48
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	2.18	1.91-2.45
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.92	1.72-2.12
<b>(Renin)肾素</b>			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	μU/mL	15.8	8.9-22.7
<b>(SCC)鳞状细胞癌抗原</b>			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.04	0.8-1.28
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.17	0.91-1.43
<b>(SHBG)性激素结合球蛋白</b>			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	nmol/L	35.3	28.1-42.5
<b>(TES)睾酮</b>			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.12	5.47-6.77
Autolumo A6000 Series Autobio \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.68	6.41-8.95
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.07	0-12.83
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.22	5.69-6.75
Tosoh AIA Series OTHER \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.72	7.43-8.01
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.67	6.41-8.93
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.66	6.57-8.75
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	7.35	6.59-8.11
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.44	6.23-8.65
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.86	0-10.69

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								12 / 14

Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.81	4.94-6.68
<b>(Tg)甲状腺球蛋白</b>			
Autolumo A6000 Series Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.45	3.2-3.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.56	7.98-9.14
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	9.4-11.4
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.88	0.31-9.45
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.1	8.3-11.9
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.14	3.9-4.38
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	3.8	3.04-4.56
<b>(TnI)肌钙蛋白I</b>			
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.058	0.047-0.069
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.01	0.01-0.01
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02	0.02-0.02
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.057	0.05-0.064
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.06	0.048-0.072
<b>(TnT)肌钙蛋白T</b>			
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.019	0.016-0.022
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.019	0.018-0.02
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02	0.013-0.027
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.018	0.014-0.022
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.019	0.018-0.02
<b>(TSH)促甲状腺素</b>			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.41	3.23-3.59
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.22	3.04-3.4
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.8	3.44-4.16
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	2.8	2.49-3.11
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.08	2.55-3.61
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	2.81	2.01-3.61
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.01	2.43-3.59
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.61	2.43-2.79
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.7	2.37-3.03
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.84	2.57-3.11
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.62	3.33-3.91
Mindray BS series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.55	3.34-3.76
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.71	2.5-2.92
Maccura i3000 MAKER \ OTHER	mIU/L	3.06	2.79-3.33
Abbott ARCHITECT ci series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.68	2.6-2.76
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.32	2.86-3.78
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.44	3.12-3.76
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	4.15	3.32-4.98
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.6	2.4-2.8
<b>(TT3)总三碘甲状腺原氨酸</b>			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.25	1.1-1.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.28	1.17-1.39

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								13 / 14

Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.1	0.88-1.32
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2	1.76-2.24
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.19	1.04-1.34
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.11	0.97-1.25
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.06	1.01-1.11
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.34	1.17-1.51
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.961	0.899-1.023
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.95	0.87-1.03
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.992	0.867-1.117
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.26	1.23-1.29
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.46	1.17-1.75

**(TT4)总T4**

Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	62.1	55.8-68.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	64.3	59.4-69.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	71	58.9-83.1
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	61.9	42.2-81.6
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	72.6	65.3-79.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	73.5	63.8-83.2
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA )	ng/mL	54.4	52-56.8
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	62	55.8-68.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	60.6	55.5-65.7
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	66.1	53.3-78.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	81.8	54.2-109.4
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA )	ng/mL	67.4	66.9-67.9

**(VB12)维生素B12**

Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	281	251-311
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	432	393-471
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	pg/mL	455	403-507
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA )	pg/mL	412	357-467
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	348	287-409
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	304	267-341
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	pg/mL	192	182-202
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	477	382-572

**(β2M)β-2-微球蛋白**

Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	ng/mL	963	840-1086
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	890	712-1068

**(β-hCG)人绒毛膜促性腺激素β亚基**

Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	0.003-0.005
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	0.004	0.004-0.004
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	0.004-0.004
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA )	IU/mL	0.003	0.003-0.003
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.46	4.16-4.76
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/mL	0.004	0.004-0.004

<b>REF</b>	540-1	6x3mL	<b>EXP</b>	2026-09-20	<b>LOT</b>	3382931	更新日期:	2025-03
								14 / 14

Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.005	0.004-0.006
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	0.003-0.005

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。