

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
1 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(a-HBDH)α-羟丁酸脱氢酶					
Beckman AU series chongqiong zybïo \ OTHER	U/L	131	2.3%	9	268
Siemens ADVIA 1800/2400 Szybio \ DGKC	U/L	110	6.8%	9	103
Beckman AU series OTHER \ OTHER	U/L	137	9.5%	8	189
Beckman AU series Szybio \ DGKC	U/L	133	2.6%	8	188
Roche Hitachi/Hitachi Systems MedicalSystem \ Enzymatic Cycling Assays Method	U/L	118	4%	7	845
KHB ZY series KeFang(Guangzhou) \ OTHER	U/L	129	3.5%	7	161
Roche Hitachi/Hitachi Systems OTHER \ DGKC	U/L	117	5.4%	6	129
TOSHIBA TBA chongqiong zybïo \ Enzymatic Cycling Assays Method	U/L	116	2.5%	5	116
TOSHIBA TBA OTHER \ OTHER	U/L	118	2.3%	5	103
Mindray BC-2000 series Mindray \ DGKC	U/L	119	4%	5	101
Beckman UniCel DxÇ 600/800 DiaSys \ OTHER	U/L	97.2	7.4%	47	964
Beckman AU series DENUO \ Enzymatic Cycling Assays Method	U/L	128	1.9%	<5	<100
Beckman AU series BSBE \ OTHER	U/L	108	5.4%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	U/L	120	5.3%	33	939
Mindray BS series MAKER \ Enzymatic Cycling Assays Method	U/L	117	4.7%	<5	<100
Roche Hitachi/Hitachi Systems MAKER \ OTHER	U/L	129	2.2%	<5	<100
Roche Hitachi/Hitachi Systems BSBE \ Enzymatic Cycling Assays Method	U/L	115	2.7%	<5	<100
Mindray BS series Mindray \ OTHER	U/L	123	5.2%	<5	<100
Beckman AU series Dedicated \ Enzymatic Cycling Assays Method	U/L	116	5.9%	29	1381
(Alb)白蛋白					
Beckman AU series Dedicated \ Bromocresol green	g/L	42.6	1.5%	91	3796
DiruiCS series DIRUI \ Bromocresol green	g/L	46.7	4.5%	9	<100
Beckman AU series Szybio \ Bromocresol green	g/L	43.9	1.7%	9	197
Mindray BS series Mindray \ Bromocresol green	g/L	44.6	2.2%	9	158
Beckman AU series OTHER \ Bromocresol green	g/L	45.3	2.2%	8	197
Roche Hitachi/Hitachi Systems MedicalSystem \ Bromocresol green	g/L	44.5	2.2%	7	977
Beckman AU series KeFang(Guangzhou) \ Bromocresol green	g/L	42.8	1%	7	154
Roche Hitachi/Hitachi Systems OTHER \ Bromocresol green	g/L	44.2	1.6%	7	104
Roche cobas c 501/502/701/702 Roche \ Bromocresol green	g/L	44.4	1.7%	6	173
Roche Modular P series G-CELL \ Bromocresol green	g/L	45	2.2%	5	<100
Abbott ARCHITECT c series MAKER \ Bromocresol green	g/L	43.7	2.5%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Bromocresol green	g/L	44.7	0.7%	5	<100
KHB ZY series KHB \ Bromocresol green	g/L	43.3	3.1%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Bromocresol green	g/L	44.1	1.3%	5	<100
Roche Hitachi/Hitachi Systems Wako \ Bromocresol green	g/L	43.7	0.9%	5	138
TOSHIBA TBA chongqiong zybïo \ Bromocresol green	g/L	45.1	1%	5	115
TOSHIBA TBA OTHER \ Bromocresol green	g/L	43.9	0.9%	5	109
Mindray BC-2000 series Mindray \ Bromocresol green	g/L	44	1.8%	5	103
Beckman UniCel DxÇ 600/800 Beckman Coulter \ Bromocresol purple	g/L	42.4	2%	49	1014
(ALP)碱性磷酸酶					
Beckman AU series Dedicated \ Rate -AMP buffer	U/L	99.5	5%	91	3778
Beckman AU series Beckman Coulter \ OTHER	U/L	109	2.2%	9	265
Roche Hitachi/Hitachi Systems MAKER \ Rate -AMP buffer	U/L	77.2	12.7%	9	209
Beckman AU series OTHER \ OTHER	U/L	95.2	13.8%	8	187
Beckman AU series Szybio \ OTHER	U/L	98.6	3.9%	8	185
Beckman AU series Beckman Coulter \ Rate -AMP buffer	U/L	93.8	7.9%	7	<100

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
2 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche Hitachi/Hitachi Systems MedicalSystem \ Rate -AMP buffer	U/L	96.6	1.9%	7	911
Beckman AU series KeFang(Guangzhou) \ Rate -AMP buffer	U/L	93.6	2.8%	7	155
Roche Hitachi/Hitachi Systems OTHER \ Rate -AMP buffer	U/L	89.1	6.3%	7	105
KHB ZY series KHB \ Rate -AMP buffer	U/L	92.9	6.5%	6	<100
Mindray BS series Mindray \ dry chemistry	U/L	81	6.8%	5	<100
Roche Modular P series G-CELL \ Rate -AMP buffer	U/L	78.3	5.7%	5	<100
Abbott ARCHITECT c series MAKER \ Rate -AMP buffer	U/L	81.6	3.9%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Rate -AMP buffer	U/L	94.5	1.8%	5	<100
TOSHIBA TBA chongqiong zybio \ Rate -AMP buffer	U/L	86.8	2.6%	5	118
Mindray BC-2000 series Mindray \ Rate -AMP buffer	U/L	90.7	2%	5	103
TOSHIBA TBA OTHER \ Rate -AMP buffer	U/L	85.9	2.9%	5	102
Beckman UniCel DxC 600/800 Beckman Coulter \ Rate -AMP buffer	U/L	96.3	5.2%	49	1023
Roche cobas 4000/6000/8000 Roche \ Rate -AMP buffer	U/L	95.6	14.8%	43	1060
(ALT/ SGPT)丙氨酸氨基转移酶					
Beckman AU series Dedicated \ UV-rate without P5P	U/L	40.3	3.6%	91	3728
Beckman AU series OTHER \ OTHER	U/L	42.8	4.7%	8	184
Beckman AU series Szybjo \ OTHER	U/L	46	3.7%	8	174
Roche Hitachi/Hitachi Systems Wako \ OTHER	U/L	42.2	4.6%	8	171
Roche Hitachi/Hitachi Systems MedicalSystem \ UV-rate with P5P	U/L	44.7	3.6%	7	977
Roche Hitachi/Hitachi Systems OTHER \ OTHER	U/L	42.2	3.9%	7	<100
Beckman AU series KeFang(Guangzhou) \ OTHER	U/L	43.1	1.8%	7	150
KHB ZY series KHB \ UV-rate with P5P	U/L	43.4	6%	6	<100
TOSHIBA TBA chongqiong zybjo \ UV-rate without P5P	U/L	50.1	2.5%	6	112
Roche Modular P series G-CELL \ UV-rate without P5P	U/L	42.3	4%	5	<100
Mindray BS series Mindray \ dry chemistry	U/L	47.2	3.4%	5	<100
Abbott ARCHITECT c series MAKER \ UV-rate without P5P	U/L	42.1	3.8%	5	<100
Siemens ADVIA 1800/2400 Siemens \ UV-rate without P5P	U/L	41.5	2.7%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ UV-rate without P5P	U/L	45.9	3.7%	5	<100
Mindray BS series Dedicated \ UV-rate without P5P	U/L	43.8	3.4%	5	294
TOSHIBA TBA OTHER \ OTHER	U/L	42	4.3%	5	243
Roche cobas c 501/502/701/702 Roche \ UV-rate without P5P	U/L	37.8	3.9%	5	149
Roche Hitachi/Hitachi Systems Wako \ UV-rate without P5P	U/L	42.7	2.2%	5	138
Mindray BC-2000 series Mindray \ OTHER	U/L	46.6	4.5%	5	103
(AMY)淀粉酶					
Abbott ARCHITECT ci series BIOSINO \ PNP-G7	U/L	95.9	1.9%	9	242
Roche Hitachi/Hitachi Systems WanTai DRD \ PNP-G7	U/L	103	11.8%	8	177
Beckman AU series KeFang(Guangzhou) \ CNP_G3	U/L	84.8	2.9%	7	112
TOSHIBA TBA WanTai \ PNP-G7	U/L	87.9	1.6%	5	<100
Roche Modular P series G-CELL \ PNP-G7	U/L	94.9	4.1%	5	<100
Abbott ARCHITECT c series MAKER \ PNP-G7	U/L	100	2.6%	5	<100
Roche Hitachi/Hitachi Systems WanTai \ PNP-G7	U/L	85.4	1.8%	5	<100
KHB ZY series KHB \ CNP_G3	U/L	98.1	4%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ CNP_G3	U/L	94.4	3%	5	<100
Mindray BS series Dedicated \ Iodine colorimetry	U/L	117	2.4%	5	295
Roche Hitachi/Hitachi Systems MAKER \ PNP-G5	U/L	86.9	0.6%	5	123
TOSHIBA TBA chongqiong zybjo \ PNP-G7	U/L	95.6	2%	5	119
Beckman AU series DENUO \ CNP_G3	U/L	84.3	2.7%	<5	<100

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
3 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
OTHER	U/L	96.4	2.7%	<5	<100
3V-Diagnostic \ PNP-G7	U/L	103	0.8%	<5	<100
ABII Real-Time PCR	U/L	96.7	2.2%	<5	107
MedicalSystem \ PNP-G7	U/L	92.3	4.7%	34	688
Beckman AU series	U/L	97.1	1.9%	32	866
Szybio \ Iodine colorimetry	U/L	99.3	1.1%	<5	990
Beckman UniCel DxC 600/800	U/L				
KHB \ Iodine colorimetry	U/L				
Roche cobas 4000/6000/8000	U/L				
Roche \ Iodine colorimetry	U/L				
Roche Hitachi/Hitachi Systems	U/L				
MedicalSystem \ PNP-G7	U/L				

(Apo-A I)载脂蛋白A I

Beckman AU series	g/L	1.69	6.7%	9	263
KINGSBIO \ OTHER	g/L	1.55	3.5%	7	139
Beckman AU series	g/L	1.65	4.8%	62	2865
KeFang(Guangzhou) \ immunoturbidimetric	g/L	1.44	5.6%	6	<100
Beckman AU series	g/L	1.93	3.4%	5	<100
DongOu \ immunoturbidimetric	g/L	1.55	3.3%	5	<100
KHB ZY series	g/L	1.59	2.9%	5	<100
KHB \ immunoturbidimetric	g/L	1.62	1.1%	5	<100
Roche cobas 4000/6000/8000	g/L	1.62	1.1%	5	<100
Dedicated \ immunoturbidimetric	g/L	1.65	5.9%	5	168
TOSHIBA TBA	g/L	1.61	2.4%	5	147
OTHER \ OTHER	g/L	1.72	1.4%	5	100
Roche Hitachi/Hitachi Systems	g/L	1.78	4%	49	1020
OTHER \ OTHER	g/L	1.43	4.9%	<5	<100
Roche Hitachi/Hitachi Systems	g/L	1.67	3.2%	<5	<100
BSBE \ immunoturbidimetric	g/L	1.54	7.7%	<5	<100
Roche Hitachi/Hitachi Systems	g/L	1.76	6.6%	32	877
MAKER \ immunoturbidimetric	g/L	1.76	5%	31	1471
Roche cobas c 501/502/701/702	g/L	1.79	3.3%	22	264
Roche \ immunoturbidimetric	g/L	1.67	4%	21	547
Mindray BC-2000 series	g/L				
Mindray \ immunoturbidimetric	g/L				
Beckman UniCel Dx C 600/800	g/L				
Beckman Coulter \ immunoturbidimetric	g/L				
Beckman AU series	g/L				
KUAKE \ immunoturbidimetric	g/L				
Beckman AU series	g/L				
Szybio \ immunoturbidimetric	g/L				
Beckman AU series	g/L				
BSBE \ immunoturbidimetric	g/L				
Roche cobas 4000/6000/8000	g/L				
Roche \ immunoturbidimetric	g/L				
Siemens ADVIA 1800/2400	g/L				
BIOSINO \ immunoturbidimetric	g/L				
Beckman AU series	g/L				
HuaChenBio \ OTHER	g/L				
Roche Hitachi/Hitachi Systems	g/L				
MAKER \ OTHER	g/L				

(Apo-B)载脂蛋白B

Beckman AU series	g/L	0.649	3.9%	9	252
KINGSBIO \ OTHER	g/L	0.773	5%	7	116
Beckman AU series	g/L	0.684	6.2%	62	2884
KeFang(Guangzhou) \ immunoturbidimetric	g/L	0.673	5.2%	6	<100
Beckman AU series	g/L	0.7	3.9%	6	127
DongOu \ immunoturbidimetric	g/L	0.752	1.4%	5	<100
KHB ZY series	g/L	0.612	2.4%	5	<100
KHB \ immunoturbidimetric	g/L	0.754	4.1%	5	<100
Roche cobas 4000/6000/8000	g/L	0.781	3.2%	5	160
Dedicated \ immunoturbidimetric	g/L	0.687	3.3%	5	147
Mindray BC-2000 series	g/L	0.571	1.8%	5	109
Mindray \ immunoturbidimetric	g/L	0.734	6.8%	49	1016
Roche Hitachi/Hitachi Systems	g/L	0.717	5.3%	<5	<100
OTHER \ OTHER	g/L	0.732	6.4%	<5	<100
Roche Hitachi/Hitachi Systems	g/L	0.837	6.1%	32	891
BSBE \ immunoturbidimetric	g/L	0.763	4.2%	31	1472
Roche cobas 4000/6000/8000	g/L	0.759	5.3%	26	373
Roche \ immunoturbidimetric	g/L	0.791	4.3%	21	561
Siemens ADVIA 1800/2400	g/L				
BIOSINO \ immunoturbidimetric	g/L				
Beckman AU series	g/L				
HuaChenBio \ OTHER	g/L				
Roche Hitachi/Hitachi Systems	g/L				
MAKER \ OTHER	g/L				
Beckman AU series	g/L				
Szybio \ immunoturbidimetric	g/L				

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
4 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(AST/ SGOT)天门冬氨酸氨基转移酶					
Beckman AU series Dedicated \ UV-rate without P5P	U/L	38.1	5.2%	92	3773
Abbott ARCHITECT ci series BIOSINO \ UV-rate with P5P	U/L	35.2	9.5%	9	245
Beckman AU series Szybio \ OTHER	U/L	36.5	3.5%	9	191
Beckman AU series OTHER \ OTHER	U/L	32.4	4.5%	8	194
Roche Hitachi/Hitachi Systems Wako \ OTHER	U/L	38.6	5%	8	162
Beckman AU series Beckman Coulter \ UV-rate without P5P	U/L	38.4	7.7%	8	108
Roche Hitachi/Hitachi Systems MedicalSystem \ UV-rate with P5P	U/L	39.8	3.5%	7	999
Roche Hitachi/Hitachi Systems OTHER \ OTHER	U/L	35.2	5.3%	7	<100
Roche Hitachi/Hitachi Systems BSBE \ UV-rate without P5P	U/L	38.4	6.5%	7	<100
Beckman AU series KeFang(Guangzhou) \ OTHER	U/L	40.5	2.8%	6	<100
KHB ZY series KHB \ UV-rate with P5P	U/L	37.3	5.8%	6	<100
TOSHIBA TBA chongqiong zybjo \ UV-rate without P5P	U/L	40.1	2.1%	6	116
Roche Modular P series G-CELL \ UV-rate without P5P	U/L	33.3	3.5%	5	<100
Mindray BS series Mindray \ dry chemistry	U/L	36.1	5.1%	5	<100
Abbott ARCHITECT c series MAKER \ UV-rate without P5P	U/L	36.8	4.3%	5	<100
Mindray BS series Dedicated \ UV-rate without P5P	U/L	37.2	4%	5	294
TOSHIBA TBA OTHER \ OTHER	U/L	36	3.8%	5	225
Roche cobas c 501/502/701/702 Roche \ UV-rate without P5P	U/L	34.1	5.5%	5	147
Roche Hitachi/Hitachi Systems Wako \ UV-rate without P5P	U/L	37.5	1.9%	5	138
(Ca)钙					
Beckman AU series Dedicated \ Arsenazo III	mmol/L	2.31	1.4%	91	3787
Abbott ARCHITECT ci series BIOSINO \ Arsenazo III	mmol/L	2.21	1.9%	9	246
Beckman AU series Szybio \ Arsenazo III	mmol/L	2.32	2.2%	8	157
Roche Modular P series G-CELL \ Arsenazo III	mmol/L	2.28	5.1%	5	<100
Mindray BC-2000 series Mindray \ Arsenazo III	mmol/L	2.35	1.8%	5	<100
Abbott ARCHITECT c series MAKER \ Arsenazo III	mmol/L	2.27	3.1%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Arsenazo III	mmol/L	2.31	1.4%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Arsenazo III	mmol/L	2.25	2.4%	5	<100
Roche cobas c 501/502/701/702 Roche \ BAPTA	mmol/L	2.29	1.2%	5	146
Roche Hitachi/Hitachi Systems Wako \ MTB colorimetry	mmol/L	2.23	1.5%	5	138
TOSHIBA TBA MedicalSystem \ Arsenazo III	mmol/L	2.29	1.8%	5	125
Beckman UniCel DxC 600/800 Beckman Coulter \ ISE (indirect)	mmol/L	2.26	1.5%	49	1023
Beckman AU series G-CELL \ Arsenazo III	mmol/L	2.26	2.5%	<5	<100
OTHER 3V-Diagnostic \ Arsenazo III	mmol/L	2.18	3%	<5	<100
PSD OTHER \ ISE (direct)	mmol/L	2.42	1.3%	<5	<100
Beckman AU series BSBE \ Arsenazo III	mmol/L	2.21	1.8%	<5	<100
Roche cobas 4000/6000/8000 Roche \ ISE (indirect)	mmol/L	2.28	4%	34	895
Roche Hitachi/Hitachi Systems MAKER \ Arsenazo III	mmol/L	2.3	3.6%	31	780
Mindray BS series MAKER \ Arsenazo III	mmol/L	2.16	1.8%	<5	<100
(CHE)胆碱酯酶					
Abbott ARCHITECT ci series BIOSINO \ DTNB	U/L	8002	2.9%	9	243
Beckman AU series Autobio \ Butyrylthiocholine Method	U/L	7887	7.6%	9	176
Beckman AU series Beckman Coulter \ OTHER	U/L	7770	4.9%	9	150
Roche Hitachi/Hitachi Systems MedicalSystem \ Butyrylthiocholine Method	U/L	7897	8%	9	1089
Beckman AU series Szybio \ Butyrylthiocholine Method	U/L	7858	2.8%	8	187
Beckman AU series KeFang(Guangzhou) \ OTHER	U/L	8163	1.2%	7	149

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
5 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA	U/L	7175	8.8%	7	122
OTHER \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	6932	9.2%	7	106
OTHER \ OTHER					
Roche cobas 4000/6000/8000	U/L	7803	8.1%	54	1317
Roche \ Butyrylthiocholine Method					
Beckman UniCel Dx C 600/800	U/L	7774	9.1%	52	1144
KHB \ Butyrylthiocholine Method					
Roche Hitachi/Hitachi Systems	U/L	7208	17%	5	<100
chongqiong zybio \ Butyrylthiocholine Method					
Roche cobas 4000/6000/8000	U/L	7185	2.1%	5	119
Roche \ DTNB					
TOSHIBA TBA	U/L	7638	1.4%	5	118
chongqiong zybio \ Butyrylthiocholine Method					
Mindray BC-2000 series	U/L	8048	1.4%	5	102
Mindray \ Butyrylthiocholine Method					
Roche cobas c 501/502/701/702	U/L	7572	28.5%	<5	<100
Roche \ Butyrylthiocholine Method					
TOSHIBA TBA	U/L	5015	8.2%	<5	<100
MedicalSystem \ Butyrylthiocholine Method					
Roche Hitachi/Hitachi Systems	U/L	8526	1.9%	<5	<100
OTHER \ DTNB					
Beckman AU series	U/L	7830	1.3%	<5	<100
LEADMAN \ Butyrylthiocholine Method					
Beckman AU series	U/L	9359	2.1%	<5	<100
DongOu \ DTNB					
(CHOL/ TC)总胆固醇					
Mindray BS series	mmol/L	3.91	2.2%	9	160
Mindray \ Cholesterol oxidase					
Roche Hitachi/Hitachi Systems	mmol/L	3.88	1.2%	7	999
MedicalSystem \ Cholesterol oxidase					
Roche Hitachi/Hitachi Systems	mmol/L	3.77	2.7%	7	<100
BSBE \ Cholesterol oxidase					
Beckman AU series	mmol/L	3.98	1.2%	7	160
KeFang(Guangzhou) \ Cholesterol oxidase					
Roche Hitachi/Hitachi Systems	mmol/L	3.92	1.6%	7	103
OTHER \ Cholesterol oxidase					
Beckman AU series	mmol/L	3.85	2%	62	2889
Dedicated \ Cholesterol oxidase					
Roche Modular P series	mmol/L	3.81	2.4%	5	<100
G-CELL \ Cholesterol oxidase					
Mindray BC-2000 series	mmol/L	3.81	1.4%	5	<100
Mindray \ Cholesterol oxidase					
Abbott ARCHITECT c series	mmol/L	3.87	3.3%	5	<100
MAKER \ Cholesterol oxidase					
Siemens ADVIA 1800/2400	mmol/L	3.86	1.6%	5	<100
LEADMAN \ Cholesterol oxidase					
KHB ZY series	mmol/L	3.78	5.1%	5	<100
KHB \ Cholesterol oxidase					
Roche cobas c 501/502/701/702	mmol/L	3.75	2.2%	5	149
Roche \ Cholesterol oxidase					
Roche Hitachi/Hitachi Systems	mmol/L	3.73	0.9%	5	138
Wako \ Cholesterol oxidase					
TOSHIBA TBA	mmol/L	3.82	1.8%	5	113
chongqiong zybio \ Cholesterol oxidase					
TOSHIBA TBA	mmol/L	3.85	0.9%	5	107
OTHER \ Cholesterol oxidase					
Beckman UniCel Dx C 600/800	mmol/L	3.92	3.1%	49	1008
Beckman Coulter \ OTHER					
Beckman AU series	mmol/L	3.99	1.9%	<5	<100
BSBE \ Cholesterol oxidase					
Beckman AU series	mmol/L	3.96	1.5%	<5	<100
Szybio \ OTHER					
OTHER	mmol/L	4	1.7%	<5	<100
3V-Diagnostic \ Cholesterol oxidase					
(CK)肌酸激酶					
Abbott ARCHITECT ci series	U/L	92.4	6.1%	9	246
BIOSINO \ Urease-UV-rate					
Beckman AU series	U/L	81.1	5.5%	9	208
Beckman Coulter \ OTHER					
Beckman AU series	U/L	104	15.8%	8	181
OTHER \ OTHER					
Beckman AU series	U/L	97.2	3.8%	8	178
Szybio \ OTHER					
Beckman AU series	U/L	91.4	5%	62	2786
WEGO \ Rate -AMP buffer					
Beckman AU series	U/L	72.5	12.7%	6	<100
KeFang(Guangzhou) \ Rate -other buffers					
Mindray BC-2000 series	U/L	85.5	12.9%	5	<100
Mindray \ OTHER					
Roche Modular P series	U/L	98.7	5%	5	<100
G-CELL \ Rate -other buffers					
Mindray BS series	U/L	109	6.8%	5	<100
Mindray \ dry chemistry					
Abbott ARCHITECT c series	U/L	71.8	6.4%	5	<100
MAKER \ Rate -other buffers					
Roche Hitachi/Hitachi Systems	U/L	93.6	4.4%	5	<100
BSBE \ Urease-UV-rate					
Mindray BS series	U/L	97.4	4.2%	5	292
Dedicated \ Rate -other buffers					
Roche Hitachi/Hitachi Systems	U/L	75.9	5.7%	5	137
MAKER \ Rate -AMP buffer					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
6 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA	U/L	89.7	5.2%	5	123
OTHER \ Rate -other buffers					
TOSHIBA TBA	U/L	89.2	5.7%	5	116
chongqiong zybio \ Rate -AMP buffer					
Beckman UniCel DxC 600/800	U/L	97.8	5.8%	49	1021
DiaSys \ OTHER					
Roche cobas 4000/6000/8000	U/L	93.8	10.2%	42	983
Roche \ Rate -AMP buffer					
Beckman AU series	U/L	113	4.5%	<5	<100
DENUO \ Urease-UV-rate					
Roche Hitachi/Hitachi Systems	U/L	98.4	5%	<5	209
BIOSINO \ OTHER					
(CK-MB)肌酸激酶-MB (活性)					
Abbott ARCHITECT ci series	U/L	15	12.1%	9	243
BIOSINO \ Immune suppression					
Beckman AU series	U/L	18.5	5.5%	9	241
chongqiong zybio \ OTHER					
Siemens ADVIA 1800/2400	U/L	14.5	4.5%	9	130
Szybio \ Immune suppression					
Roche Hitachi/Hitachi Systems	U/L	14	9.8%	8	172
MAKER \ OTHER					
Beckman AU series	U/L	12.2	10%	7	114
KeFang(Guangzhou) \ Immune suppression					
KHB ZY series	U/L	15.4	9.5%	6	<100
KHB \ Immune suppression					
Beckman AU series	U/L	20.2	6.7%	6	200
WEGO \ Immune suppression					
Roche Hitachi/Hitachi Systems	U/L	14.8	5.8%	6	129
OTHER \ DGKC					
Mindray BC-2000 series	U/L	11.2	26.1%	5	<100
Mindray \ Immune suppression					
Beckman AU series	U/L	18.1	10.5%	5	<100
Szybio \ Immune suppression					
Roche Hitachi/Hitachi Systems	U/L	14.1	12.9%	5	<100
BSBE \ DGKC					
Mindray BS series	U/L	16.3	7%	5	293
BIOSINO \ Immune suppression					
Beckman UniCel DxC 600/800	U/L	13.6	11.7%	49	1008
DiaSys \ OTHER					
Beckman AU series	U/L	19.4	6.7%	<5	<100
DENUO \ DGKC					
Roche cobas 4000/6000/8000	U/L	15.9	11.2%	32	895
Roche \ DGKC					
Mindray BS series	U/L	14.1	9.4%	<5	<100
MAKER \ Enzyme immunoassay (EIA)					
Beckman AU series	U/L	14.1	8.2%	<5	<100
BSBE \ Immune suppression					
Roche Hitachi/Hitachi Systems	U/L	13.6	10.3%	27	697
MAKER \ Immune suppression					
Beckman AU series	U/L	11.1	18.2%	22	932
Dedicated \ Immune suppression					
(Cl)氯					
Beckman AU series	mmol/L	117	1.1%	92	3842
Dedicated \ ISE (indirect)					
Abbott ARCHITECT ci series	mmol/L	120	1.1%	9	245
BIOSINO \ ISE (indirect)					
Beckman AU series	mmol/L	118	3.5%	9	197
OTHER \ ISE (direct)					
Beckman AU series	mmol/L	115	1.8%	8	176
Szybio \ ISE (indirect)					
MZKL	mmol/L	114	4.1%	7	132
OTHER \ ISE (direct)					
TOSHIBA TBA	mmol/L	116	0.9%	6	118
chongqiong zybio \ ISE (indirect)					
Beckman UniCel DxC 600/800	mmol/L	115	1.5%	50	1018
Beckman Coulter \ ISE (indirect)					
Mindray BC-2000 series	mmol/L	119	1%	5	<100
Mindray \ ISE (indirect)					
Roche Modular P series	mmol/L	111	2.3%	5	<100
Roche \ ISE (direct)					
Abbott ARCHITECT c series	mmol/L	114	1.9%	5	<100
Abbott \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	113	1.6%	5	<100
BSBE \ ISE (direct)					
Mindray BS series	mmol/L	119	1.3%	5	295
Dedicated \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	105	0.6%	5	138
Wako \ ISE (indirect)					
Roche cobas c 501/502/701/702	mmol/L	116	1%	5	137
Roche \ ISE (direct)					
TOSHIBA TBA	mmol/L	117	0.7%	5	109
OTHER \ ISE (indirect)					
OTHER	mmol/L	99	1.8%	<5	<100
3V-Diagnostic \ ISE (indirect)					
PSD	mmol/L	121	1.1%	<5	<100
OTHER \ ISE (direct)					
ABI Real-Time PCR	mmol/L	114	1.3%	<5	<100
Dedicated \ ISE (direct)					
Roche Hitachi/Hitachi Systems	mmol/L	112	1.6%	<5	115
Hitachi \ ISE (direct)					

质控品名： 生化与激素非定值质控品

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
7 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Cr)肌酐					
DiruiCS series	μmol/L	82	8.6%	9	<100
DIRUI \ Enzymatic					
Beckman UniCel DxC 600/800	μmol/L	76.8	7.1%	9	<100
BGH \ OTHER					
Abbott ARCHITECT ci series	μmol/L	73.9	3.8%	9	244
BIOSINO \ Picric acid-Blank Compensation					
Roche Hitachi/Hitachi Systems	μmol/L	76.9	3.2%	8	199
Whitman \ Enzymatic					
Beckman AU series	μmol/L	75.5	2.2%	8	182
Szybio \ OTHER					
Beckman AU series	μmol/L	75.6	1.4%	7	136
KeFang(Guangzhou) \ Enzymatic					
Beckman AU series	μmol/L	72	4%	60	2326
Dedicated \ Picric acid-No blank Compensation					
Roche cobas 4000/6000/8000	μmol/L	86.9	2.3%	6	<100
Dedicated \ Enzymatic					
KHB ZY series	μmol/L	79.6	4.3%	6	<100
KHB \ Enzymatic					
Mindray BS series	μmol/L	79	1.5%	6	121
MAKER \ Enzymatic					
Roche Modular P series	μmol/L	76.5	4.1%	5	<100
WanTai DRD \ Enzymatic					
Roche Hitachi/Hitachi Systems	μmol/L	76.5	2%	5	<100
OTHER \ Enzymatic					
Abbott ARCHITECT c series	μmol/L	71	17.3%	5	<100
MAKER \ Enzymatic					
Siemens ADVIA 1800/2400	μmol/L	78.2	2%	5	<100
OTHER \ Enzymatic					
Roche Hitachi/Hitachi Systems	μmol/L	78.1	4.8%	5	<100
BSBE \ Enzymatic					
Mindray BS series	μmol/L	80.1	2.2%	5	294
Dedicated \ Enzymatic					
TOSHIBA TBA	μmol/L	76.2	2.1%	5	227
OTHER \ Enzymatic					
Roche cobas c 501/502/701/702	μmol/L	83.3	3%	5	149
Roche \ Enzymatic					
Roche Hitachi/Hitachi Systems	μmol/L	71	1.5%	5	133
Wako \ Enzymatic					
(CRP)C反应蛋白					
Mindray BS series	mg/dL	1.03	26.9%	6	331
SHSN \ Immunonephelometry (INA)					
Beckman AU series	mg/dL	0.982	4%	<5	<100
G-CELL \ OTHER					
Beckman AU series	mg/dL	0.753	3.7%	<5	<100
OTHER \ Immunonephelometry (INA)					
Roche Hitachi/Hitachi Systems	mg/dL	0.63	71%	13	531
DiaSys \ Immunoturbidimetry (ITA)					
Beckman AU series	mg/dL	0.616	4.2%	10	168
SekiSui \ Immunonephelometry (INA)					
Beckman Immage 800	mg/dL	0.589	9.8%	<5	<100
Beckman Coulter \ Immunoturbidimetry (ITA)					
(Cu)铜					
Beckman AU series	μmol/L	17.2	7.7%	6	<100
BSBE \ PAESA					
OTHER	μmol/L	15.6	4.5%	<5	<100
3V-Diagnostic \ Colorimetry					
Beckman AU series	μmol/L	18.8	6.8%	36	1655
BSBE \ OTHER					
Roche Hitachi/Hitachi Systems	μmol/L	17.3	6.5%	10	153
BSBE \ PAESA					
Beckman AU series	μmol/L	20.4	10%	<5	<100
Beckman Coulter \ ICP-MS					
(DBIL)直接胆红素					
Beckman AU series	μmol/L	3.94	4.4%	92	3759
Dedicated \ Diazotization -Blank Compensation					
Beckman AU series	μmol/L	4.1	4%	9	268
Beckman Coulter \ OTHER					
Abbott ARCHITECT ci series	μmol/L	4.56	5.8%	9	243
BIOSINO \ Vanadate					
Beckman AU series	μmol/L	5.93	6.7%	8	120
Szybio \ OTHER					
Beckman AU series	μmol/L	4.54	11.1%	7	134
KeFang(Guangzhou) \ Vanadate					
KHB ZY series	μmol/L	6.15	14.2%	6	<100
KHB \ Vanadate					
Beckman AU series	μmol/L	6.98	4.6%	6	143
OTHER \ OTHER					
Roche Modular P series	μmol/L	6.39	11.2%	5	<100
G-CELL \ Vanadate					
Mindray BS series	μmol/L	7.03	3.9%	5	<100
Mindray \ Diazotization					
Abbott ARCHITECT c series	μmol/L	6.2	5.1%	5	<100
Abbott \ Vanadate					
Roche Hitachi/Hitachi Systems	μmol/L	6.6	3.5%	5	<100
Szybio \ Vanadate					
Roche Hitachi/Hitachi Systems	μmol/L	7.27	6.2%	5	<100
BSBE \ Vanadate					
Mindray BS series	μmol/L	4.59	4.8%	5	294
Dedicated \ Vanadate					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
8 / 17								

		批号：376B805			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA Szybio \ Vanadate	μmol/L	3.38	7.3%	5	237
Roche cobas c 501/502/701/702 Roche \ Diazotization -No blank Compensation	μmol/L	5.98	5.1%	5	148
TOSHIBA TBA chongqiong zzybio \ Vanadate	μmol/L	6.67	2.6%	5	114
Mindray BC-2000 series Mindray \ Vanadate	μmol/L	5.88	4.3%	5	103
Mindray BS series Mindray \ Vanadate	μmol/L	6.17	4.8%	<5	<100
Beckman AU series Wako \ Vanadate	μmol/L	8.06	6.1%	<5	<100
(DGXN)地高辛					
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.75	10%	<5	<100
(E2)雌二醇					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	129	10.1%	<5	<100
(FA)叶酸					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.804	9.3%	<5	<100
(Fe)铁					
Roche Hitachi/Hitachi Systems BSBE \ Ferrozine colorimetry	μmol/L	18	2.1%	5	<100
Mindray BC-2000 series Mindray \ Ferrozine colorimetry	μmol/L	18.5	1.4%	5	103
OTHER 3V-Diagnostic \ Ferrozine colorimetry	μmol/L	19.6	5.5%	<5	<100
Beckman AU series Szybio \ Ferrozine colorimetry	μmol/L	18	1.4%	<5	<100
Beckman AU series Autobio \ Ferrozine colorimetry	μmol/L	19.6	4%	<5	<100
Beckman AU series BSBE \ OTHER	μmol/L	18.1	3.7%	<5	<100
Siemens ADVIA 1800/2400 BIOSINO \ Ferrozine colorimetry	μmol/L	18.2	3.6%	31	1469
Beckman AU series Dedicated \ Ferrozine colorimetry	μmol/L	18.3	2.6%	31	1457
Roche cobas 4000/6000/8000 Roche \ Ferrozine colorimetry	μmol/L	19.5	11.7%	<5	<100
Beckman AU series Beckman Coulter \ Ferrozine colorimetry	μmol/L	20.9	10%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Ferrozine colorimetry	μmol/L	19.4	2.6%	<5	<100
(Ferritin)铁蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	58.6	9.7%	<5	<100
(FT3)游离T3					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.94	9.8%	<5	<100
(FT4)游离T4					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.3	9.2%	<5	<100
(GGT)L-γ-谷氨酰基转移酶					
Beckman AU series Dedicated \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	38.3	2.2%	91	3804
Beckman AU series Beckman Coulter \ OTHER	U/L	38.6	1.7%	9	268
Beckman AU series OTHER \ OTHER	U/L	39.5	9.3%	8	181
Beckman AU series Szybio \ OTHER	U/L	35.3	4.2%	8	180
Roche Hitachi/Hitachi Systems Wako \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	37.7	2%	8	177
Roche Hitachi/Hitachi Systems MAKER \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	38.7	2.6%	8	150
Beckman AU series BSBE \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	39	2.7%	8	114
Beckman AU series KeFang(Guangzhou) \ OTHER	U/L	41.2	1.8%	7	140
Roche Hitachi/Hitachi Systems OTHER \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	36.1	6.4%	7	103
Roche cobas 4000/6000/8000 Dedicated \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	36.7	3.2%	6	<100
Roche Modular P series G-CELL \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	37.6	2.9%	5	<100
Mindray BS series Mindray \ dry chemistry	U/L	38.6	3.1%	5	<100
Abbott ARCHITECT c series MAKER \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	38	6.5%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	33.8	1.8%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	39.2	3.3%	5	<100

质控品名： 生化与激素非定值质控品

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
9 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Mindray BS series	U/L	38.1	2.4%	5	294
Dedicated \ γ-glutamyl -3-carboxy-4-nitroaniline					
TOSHIBA TBA	U/L	35.3	2.4%	5	240
OTHER \ γ-glutamyl -3-carboxy-4-nitroaniline					
Roche cobas c 501/502/701/702	U/L	37.8	2.8%	5	147
Roche \ γ-glutamyl -3-carboxy-4-nitroaniline					
TOSHIBA TBA	U/L	37.2	1.9%	5	112
chongqiong zybio \ γ-glutamyl -3-carboxy-4-nitroaniline					
(GLDH)谷氨酸脱氢酶					
Beckman AU series	U/L	4.93	8.5%	60	2259
Roche \ Direct Measure					
Roche Hitachi/Hitachi Systems	U/L	4.96	7.2%	6	137
OTHER \ DGKC					
Roche Hitachi/Hitachi Systems	U/L	5.59	5.5%	5	<100
BSBE \ Direct Measure					
Beckman AU series	U/L	4.95	7.5%	31	1443
Roche \ OTHER					
Beckman AU series	U/L	4.53	9.5%	<5	<100
Beckman Coulter \ DGKC					
Roche cobas 4000/6000/8000	U/L	4.78	9.8%	<5	<100
DiaSys \ Colorimetry					
(GLU)葡萄糖					
Beckman AU series	mmol/L	3.98	2.9%	91	3778
Dedicated \ Hexokinase					
Roche Hitachi/Hitachi Systems	mmol/L	4.29	3.7%	8	198
WanTai DRD \ Hexokinase					
Beckman AU series	mmol/L	4.14	1.9%	8	190
Szybio \ Hexokinase					
Beckman AU series	mmol/L	4.25	1.5%	7	163
KeFang(Guangzhou) \ Glucose oxidase (GOD)					
Roche Hitachi/Hitachi Systems	mmol/L	3.98	1.8%	7	131
OTHER \ Hexokinase					
Roche cobas 4000/6000/8000	mmol/L	4.22	2.6%	6	<100
Dedicated \ Hexokinase					
KHB ZY series	mmol/L	4.14	4%	6	<100
KHB \ Glucose oxidase (GOD)					
TOSHIBA TBA	mmol/L	4.03	2.2%	6	116
chongqiong zybio \ Hexokinase					
Roche Modular P series	mmol/L	4.1	2.8%	5	<100
G-CELL \ Glucose oxidase (GOD)					
Mindray BS series	mmol/L	4.11	2.2%	5	<100
Mindray \ Dehydrogenase					
Abbott ARCHITECT c series	mmol/L	4.14	3.8%	5	<100
MAKER \ Hexokinase					
Siemens ADVIA 1800/2400	mmol/L	4.01	1.2%	5	<100
LEADMAN \ Hexokinase					
Beckman AU series	mmol/L	3.61	2.2%	5	<100
HuaChenBio \ Glucose oxidase (GOD)					
Roche Hitachi/Hitachi Systems	mmol/L	4.24	3.5%	5	<100
BSBE \ Glucose oxidase (GOD)					
DiruiCS series	mmol/L	4.21	3.4%	5	<100
DIRUI \ Dehydrogenase					
Roche cobas c 501/502/701/702	mmol/L	3.99	2.1%	5	147
Roche \ Hexokinase					
Roche Hitachi/Hitachi Systems	mmol/L	3.88	0.9%	5	136
Wako \ Hexokinase					
Mindray BC-2000 series	mmol/L	3.93	1.7%	5	103
Mindray \ Hexokinase					
Beckman UniCel DxС 600/800	mmol/L	3.92	1.7%	49	1016
Beckman Coulter \ OTHER					
(HCY)同型半胱氨酸					
Beckman AU series	μmol/L	11.1	3.3%	9	246
OTHER \ OTHER					
Siemens ADVIA 1800/2400	μmol/L	13.5	3.6%	9	141
Szybio \ Immunonephelometry (INA)					
Beckman AU series	μmol/L	13.6	4.7%	62	2873
DongOu \ OTHER					
Mindray BC-2000 series	μmol/L	14.8	3.1%	5	<100
Mindray \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	12	9.3%	5	<100
OTHER \ Enzymatic Cycling Assays Method					
Roche Hitachi/Hitachi Systems	μmol/L	14.1	3.8%	5	<100
WanTai \ OTHER					
Roche Hitachi/Hitachi Systems	μmol/L	14.4	3.8%	5	<100
BSBE \ Enzymatic Cycling Assays Method					
TOSHIBA TBA	μmol/L	14.1	5.8%	5	114
WanTai \ OTHER					
Beckman UniCel DxС 600/800	μmol/L	14.8	8.3%	49	975
REEGEN \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	13.6	7.8%	<5	<100
KUAKE \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	14.3	10.3%	<5	<100
Szybio \ Enzymatic Cycling Assays Method					
Roche Hitachi/Hitachi Systems	μmol/L	14.8	2.9%	<5	118
MAKER \ Enzymatic Cycling Assays Method					
Roche cobas 4000/6000/8000	μmol/L	11.2	27.7%	37	940
Roche \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	13.5	12.7%	22	893
REEGEN \ Enzymatic Cycling Assays Method					
Siemens ADVIA 1800/2400	μmol/L	13.1	5.1%	21	981
BIOSINO \ Immunoturbidimetry (ITA)					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
								10 / 17

		批号：376B805			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche Hitachi/Hitachi Systems MAKER \ OTHER	μmol/L	14.7	7.9%	21	561
Beckman AU series BSBE \ Enzymatic Cycling Assays Method	μmol/L	12.8	3.5%	<5	<100
Roche Hitachi/Hitachi Systems OTHER \ Enzymatic Cycling Assays Method	μmol/L	17.7	10.4%	<5	<100
Mindray BS series Mindray \ Enzymatic Cycling Assays Method	μmol/L	11.3	18.4%	<5	<100
(HDL-C)高密度脂蛋白胆固醇					
Roche Hitachi/Hitachi Systems BSBE \ selected suppression	mmol/L	1.11	4.9%	9	<100
Beckman AU series Beckman Coulter \ Direct Measure	mmol/L	1.22	1.7%	9	268
Abbott ARCHITECT ci series BIOSINO \ selected suppression	mmol/L	1.13	2.5%	9	243
Mindray BS series Mindray \ Direct Measure	mmol/L	1.09	4.9%	9	160
Beckman AU series OTHER \ selected suppression	mmol/L	0.99	5.7%	8	191
Roche cobas 4000/6000/8000 Dedicated \ Direct Measure	mmol/L	1.18	11.8%	8	114
Roche Hitachi/Hitachi Systems OTHER \ selected suppression	mmol/L	1.07	3.3%	7	<100
Beckman AU series KeFang(Guangzhou) \ Direct Measure	mmol/L	0.998	1.6%	7	158
Beckman AU series WEGO \ Direct Measure	mmol/L	1.08	3.8%	62	2894
KHB ZY series KHB \ Direct Measure	mmol/L	1.24	11.5%	6	<100
Abbott ARCHITECT c series MAKER \ selected suppression	mmol/L	1.13	7.6%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ selected suppression	mmol/L	1.12	1.6%	5	<100
Roche cobas c 501/502/701/702 Roche \ Direct Measure	mmol/L	0.945	2.1%	5	146
Roche Hitachi/Hitachi Systems Wako \ selected suppression	mmol/L	1.09	1.8%	5	137
TOSHIBA TBA OTHER \ selected suppression	mmol/L	1.01	1.2%	5	111
Mindray BC-2000 series Mindray \ OTHER	mmol/L	1.2	1.9%	5	101
Beckman UniCel DxC 600/800 LEADMAN \ selected suppression	mmol/L	1.14	3.7%	49	1021
Beckman AU series SekiSui \ Enzymatic	mmol/L	1.32	2.6%	<5	<100
Beckman AU series Szybio \ selected suppression	mmol/L	1.07	1.5%	<5	<100
(IgA)免疫球蛋白 IgA					
Roche Hitachi/Hitachi Systems DiaSys \ Immunoturbidimetry (ITA)	mg/dL	185	2.5%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunonephelometry (INA)	mg/dL	220	2.7%	5	<100
Roche Hitachi/Hitachi Systems MAKER \ OTHER	mg/dL	209	4.4%	5	173
TOSHIBA TBA DiaSys \ Immunoturbidimetry (ITA)	mg/dL	203	1.9%	5	113
Beckman AU series OTHER \ OTHER	mg/dL	176	5.3%	<5	105
Roche Hitachi/Hitachi Systems OTHER \ Immunoturbidimetry (ITA)	mg/dL	229	3.6%	<5	<100
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	197	9.6%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Immunonephelometry (INA)	mg/dL	205	1.6%	<5	<100
Beckman AU series OTHER \ Immunonephelometry (INA)	mg/dL	172	8.3%	<5	<100
(IgG)免疫球蛋白 IgG					
TOSHIBA TBA DiaSys \ Immunoturbidimetry (ITA)	mg/dL	976	3.5%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunonephelometry (INA)	mg/dL	1087	4.2%	5	<100
Roche Hitachi/Hitachi Systems MAKER \ OTHER	mg/dL	1067	7.7%	5	169
Beckman AU series OTHER \ OTHER	mg/dL	961	5.5%	<5	105
Roche Hitachi/Hitachi Systems OTHER \ Immunoturbidimetry (ITA)	mg/dL	1284	6.5%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Immunonephelometry (INA)	mg/dL	1018	1.6%	<5	<100
Beckman AU series OTHER \ Immunonephelometry (INA)	mg/dL	933	10.9%	<5	<100
Beckman Immage 800 Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	1130	9.7%	<5	<100
(IgM)免疫球蛋白 IgM					
Roche Hitachi/Hitachi Systems DiaSys \ Immunoturbidimetry (ITA)	mg/dL	94.9	2.3%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunonephelometry (INA)	mg/dL	102	7.1%	5	<100
Roche Hitachi/Hitachi Systems MAKER \ OTHER	mg/dL	92.9	9.4%	5	175

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
11 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA	mg/dL	93	1.6%	5	120
DiaSys \ Immunoturbidimetry (ITA)					
Beckman AU series	mg/dL	84	6.1%	<5	105
OTHER \ OTHER					
Roche Hitachi/Hitachi Systems	mg/dL	84.1	2.4%	<5	<100
OTHER \ Immunoturbidimetry (ITA)					
Beckman Image 800	mg/dL	92	9.8%	<5	<100
Beckman Coulter \ Immunoturbidimetry (ITA)					
Roche cobas c 501/502/701/702	mg/dL	87.2	4.2%	<5	<100
Roche \ Immunonephelometry (INA)					
Beckman AU series	mg/dL	88.5	9%	<5	<100
OTHER \ Immunonephelometry (INA)					

(K)钾					
Beckman AU series	mmol/L	3.64	1%	92	3867
Dedicated \ ISE (indirect)					
Abbott ARCHITECT ci series	mmol/L	3.62	1.7%	9	244
BIOSINO \ ISE (indirect)					
Beckman AU series	mmol/L	3.7	2%	9	191
OTHER \ ISE (direct)					
Beckman AU series	mmol/L	3.76	2%	8	171
Szybio \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	3.65	1.4%	7	955
Dedicated \ ISE (direct)					
MZKL	mmol/L	3.9	2%	7	132
OTHER \ ISE (direct)					
Roche cobas c 501/502/701/702	mmol/L	3.72	1.1%	6	171
Roche \ ISE (indirect)					
TOSHIBA TBA	mmol/L	3.63	0.8%	6	114
chongqiong zyblio \ ISE (indirect)					
Roche Modular P series	mmol/L	3.77	2.2%	5	<100
Roche \ ISE (direct)					
Abbott ARCHITECT c series	mmol/L	3.65	2.1%	5	<100
Abbott \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	3.57	1.9%	5	<100
BSBE \ ISE (direct)					
Mindray BS series	mmol/L	3.79	1.6%	5	297
Dedicated \ ISE (indirect)					
TOSHIBA TBA	mmol/L	3.67	1.2%	5	116
OTHER \ ISE (indirect)					
Mindray BC-2000 series	mmol/L	3.73	1%	5	101
Mindray \ ISE (indirect)					
Beckman UniCel DxC 600/800	mmol/L	3.63	1.6%	49	1017
Beckman Coulter \ ISE (indirect)					
OTHER	mmol/L	3.98	1.9%	<5	<100
3V-Diagnostic \ ISE (indirect)					
PSD	mmol/L	3.83	1.2%	<5	<100
OTHER \ ISE (direct)					
Roche Hitachi/Hitachi Systems	mmol/L	3.79	1.7%	<5	208
Dedicated \ ICP-MS					
Roche Hitachi/Hitachi Systems	mmol/L	3.63	1.1%	<5	112
Hitachi \ ISE (direct)					

(LA)乳酸					
Roche Hitachi/Hitachi Systems	mmol/L	1.88	2.1%	8	190
LEADMAN \ OTHER					
Roche Hitachi/Hitachi Systems	mmol/L	1.98	0.9%	<5	<100
ABX \ dry chemistry					
Beckman AU series	mmol/L	1.96	10.1%	<5	<100
Beckman Coulter \ Lactate to pyruvate (L-P)					
Roche cobas 4000/6000/8000	mmol/L	1.82	10.4%	<5	<100
Roche \ Lactate to pyruvate (L-P)					

(LDH)乳酸脱氢酶					
Beckman AU series	U/L	133	4.2%	92	3734
Dedicated \ Lactate to pyruvate (L-P)					
Beckman AU series	U/L	143	2.8%	9	268
Beckman Coulter \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	131	2.8%	8	192
Whitman \ Lactate to pyruvate (L-P)					
Beckman AU series	U/L	138	8.8%	8	184
OTHER \ OTHER					
Beckman AU series	U/L	133	1.9%	8	173
Szybio \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	117	5.6%	8	149
MAKER \ Lactate to pyruvate (L-P)					
Roche Hitachi/Hitachi Systems	U/L	135	1.4%	7	933
MedicalSystem \ Lactate to pyruvate (L-P)					
Roche cobas 4000/6000/8000	U/L	139	3.8%	7	<100
Dedicated \ Lactate to pyruvate (L-P)					
Beckman AU series	U/L	138	2.4%	7	156
KeFang(Guangzhou) \ OTHER					
KHB ZY series	U/L	134	3.6%	6	<100
KHB \ Lactate to pyruvate (L-P)					
Mindray BC-2000 series	U/L	139	1.6%	5	<100
Mindray \ OTHER					
Mindray BS series	U/L	140	3.3%	5	<100
Mindray \ dry chemistry					
Abbott ARCHITECT c series	U/L	130	3.7%	5	<100
MAKER \ Lactate to pyruvate (L-P)					
Roche Modular P series	U/L	136	2.8%	5	<100
G-CELL \ Lactate to pyruvate (L-P)					
Roche Hitachi/Hitachi Systems	U/L	134	2.5%	5	<100
BSBE \ Lactate to pyruvate (L-P)					

质控品名： 生化与激素非定值质控品

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
12 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Mindray BS series	U/L	142	3.3%	5	294
Dedicated \ Lactate to pyruvate (L-P)					
TOSHIBA TBA chongqiong zybio \ Lactate to pyruvate (L-P)	U/L	131	2.3%	5	114
TOSHIBA TBA OTHER \ OTHER	U/L	128	1.9%	5	113
Roche Hitachi/Hitachi Systems OTHER \ OTHER	U/L	132	2.8%	5	109
(LDL-C)低密度脂蛋白胆固醇					
Beckman AU series	mmol/L	2.63	1.4%	9	268
Beckman Coulter \ Direct Measure					
Abbott ARCHITECT ci series	mmol/L	1.92	1.8%	9	243
BIOSINO \ selected suppression					
Beckman AU series	mmol/L	1.87	3.2%	7	127
KeFang(Guangzhou) \ Direct Measure					
Beckman AU series	mmol/L	1.96	3.1%	62	2896
WEGO \ Direct Measure					
Roche Hitachi/Hitachi Systems	mmol/L	2.21	14%	6	<100
OTHER \ SUR					
KHB ZY series	mmol/L	1.83	16.9%	6	<100
KHB \ Direct Measure					
Roche cobas 4000/6000/8000	mmol/L	2.38	4.3%	6	107
Dedicated \ Direct Measure					
Roche Modular P series	mmol/L	2.11	3.5%	5	<100
MAKER \ Enzymatic					
Mindray BC-2000 series	mmol/L	2.38	2.2%	5	<100
Mindray \ OTHER					
Mindray BS series	mmol/L	2.44	3.2%	5	<100
Mindray \ SUR					
Abbott ARCHITECT c series	mmol/L	1.96	3.6%	5	<100
MAKER \ SUR					
Siemens ADVIA 1800/2400	mmol/L	2.02	1.7%	5	<100
LEADMAN \ SUR					
Roche cobas c 501/502/701/702	mmol/L	2.41	2.6%	5	146
Roche \ Direct Measure					
Roche Hitachi/Hitachi Systems	mmol/L	2.25	1.7%	5	137
Wako \ selected suppression					
TOSHIBA TBA	mmol/L	2.73	5%	5	114
chongqiong zybio \ Enzymatic					
Beckman UniCel DxC 600/800	mmol/L	2.11	5%	49	1019
LEADMAN \ SUR					
TOSHIBA TBA	mmol/L	2.01	2.6%	<5	<100
OTHER \ SUR					
Mindray BS series	mmol/L	2.4	1.8%	<5	<100
Mindray \ Direct Measure					
Beckman AU series	mmol/L	2.3	2.2%	<5	<100
SekiSui \ Enzymatic					
(Li)锂					
Beckman AU series	mmol/L	0.79	8.9%	<5	<100
Beckman Coulter \ Colorimetry					
(LIP)脂肪酶					
Beckman AU series	U/L	46.9	8.1%	9	253
KINGSBIO \ OTHER					
Siemens ADVIA 1800/2400	U/L	52.9	4.1%	9	147
HOMA \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	78.2	5%	8	182
BSBE \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	53.7	5.3%	6	140
OTHER \ Glucose oxidase (GOD)					
Beckman AU series	U/L	64.1	20.7%	5	100
OTHER \ OTHER					
Roche cobas 4000/6000/8000	U/L	52.9	11.6%	31	819
Roche \ Enzymatic Cycling Assays Method					
Beckman AU series	U/L	50.3	14.3%	27	301
Autobio \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	4.31	4.6%	<5	<100
ABX \ dry chemistry					
Siemens ADVIA Chemistry XPT	U/L	75.9	8.9%	16	733
KHB \ Enzymatic Cycling Assays Method					
Beckman AU series	U/L	55.6	4.7%	11	542
Beckman Coulter \ OTHER					
Roche Hitachi/Hitachi Systems	U/L	43.4	11%	<5	<100
MAKER \ OTHER					
Roche cobas 4000/6000/8000	U/L	44.2	9.5%	<5	<100
DiaSys \ Enzymatic Cycling Assays Method					
Beckman AU series	U/L	48.5	9.9%	<5	<100
Beckman Coulter \ Glucose oxidase (GOD)					
Roche cobas c 501/502/701/702	U/L	52.8	1.1%	<5	<100
Roche \ OTHER					
(Mg)镁					
Beckman AU series	mmol/L	0.883	1.8%	92	3800
Dedicated \ Calmagite colorimetry					
Abbott ARCHITECT ci series	mmol/L	0.918	5.4%	9	244
BIOSINO \ Calmagite colorimetry					
Roche Hitachi/Hitachi Systems	mmol/L	0.901	5.7%	8	1081
MedicalSystem \ Arsenazo □					
Beckman AU series	mmol/L	0.791	3.8%	7	<100
HuaChenBio \ Arsenazo □					
Roche cobas c 501/502/701/702	mmol/L	0.874	2.6%	6	176
Roche \ Calmagite colorimetry					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
13 / 17								

		批号：376B805			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche Modular P series	mmol/L	0.919	10.2%	5	<100
G-CELL \ Xylidyl Blue					
Abbott ARCHITECT c series	mmol/L	0.848	3.5%	5	<100
MAKER \ Xylidyl Blue					
Roche Hitachi/Hitachi Systems	mmol/L	0.959	4.5%	5	<100
BSBE \ Xylidyl Blue					
TOSHIBA TBA	mmol/L	0.946	3.9%	5	120
MedicalSystem \ Arsenazo □					
TOSHIBA TBA	mmol/L	0.942	4.5%	5	115
chongqiong zybio \ Xylidyl Blue					
Mindray BC-2000 series	mmol/L	0.934	2.1%	5	103
Mindray \ Xylidyl Blue					
Beckman UniCel DxC 600/800	mmol/L	0.892	3.3%	49	1019
Beckman Coulter \ OTHER					
Beckman AU series	mmol/L	0.881	7.8%	<5	<100
Beckman Coulter \ Xylidyl Blue					
Beckman AU series	mmol/L	0.884	2.7%	<5	<100
Szybio \ OTHER					
OTHER	mmol/L	0.848	3.8%	<5	<100
3V-Diagnostic \ Xylidyl Blue					
Beckman AU series	mmol/L	0.877	11%	<5	<100
BSBE \ Xylidyl Blue					
Roche Hitachi/Hitachi Systems	mmol/L	0.922	4.7%	<5	<100
MedicalSystem \ Calmagite colorimetry					
Roche Hitachi/Hitachi Systems	mmol/L	0.861	2.9%	<5	110
MAKER \ Calmagite colorimetry					
Roche cobas 4000/6000/8000	mmol/L	0.902	6.3%	32	892
Roche \ ICP-MS					
(Na)钠					
Beckman AU series	mmol/L	140	0.9%	92	3860
Dedicated \ ISE (indirect)					
Abbott ARCHITECT ci series	mmol/L	140	1.1%	9	245
BIOSINO \ ISE (indirect)					
Beckman AU series	mmol/L	142	1.1%	9	187
OTHER \ ISE (direct)					
Beckman AU series	mmol/L	142	0.8%	8	172
Szybio \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	143	0.7%	7	955
Dedicated \ ISE (direct)					
MZKL	mmol/L	146	1.4%	7	132
OTHER \ ISE (direct)					
Roche cobas c 501/502/701/702	mmol/L	143	0.8%	6	170
Roche \ ISE (indirect)					
TOSHIBA TBA	mmol/L	142	0.8%	6	112
chongqiong zybio \ ISE (indirect)					
Roche Modular P series	mmol/L	142	1.8%	5	<100
Roche \ ISE (direct)					
Mindray BC-2000 series	mmol/L	142	0.6%	5	<100
Mindray \ ISE (indirect)					
Abbott ARCHITECT c series	mmol/L	143	1.3%	5	<100
Abbott \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	142	1.2%	5	<100
BSBE \ ISE (direct)					
Mindray BS series	mmol/L	144	1.3%	5	298
Dedicated \ ISE (indirect)					
Roche Hitachi/Hitachi Systems	mmol/L	144	0.5%	5	137
Wako \ ISE (indirect)					
TOSHIBA TBA	mmol/L	141	0.8%	5	106
OTHER \ ISE (indirect)					
Beckman UniCel DxC 600/800	mmol/L	141	1.2%	49	1022
Beckman Coulter \ ISE (indirect)					
OTHER	mmol/L	143	1.2%	<5	<100
3V-Diagnostic \ ISE (indirect)					
PSD	mmol/L	147	1.4%	<5	<100
OTHER \ ISE (direct)					
Roche Hitachi/Hitachi Systems	mmol/L	143	1.3%	<5	114
Hitachi \ ISE (direct)					
(NEFA)游离脂肪酸					
Beckman AU series	mmol/L	0.515	7.1%	9	260
KINGSBIO \ OTHER					
Beckman AU series	mmol/L	0.475	3.5%	62	2904
DongOu \ ACS - ACOD					
Roche Hitachi/Hitachi Systems	mmol/L	0.475	2.2%	6	135
OTHER \ ACS - ACOD					
Siemens ADVIA 1800/2400	mmol/L	0.5	6.6%	31	1473
OTHER \ ACS - ACOD					
Roche cobas 4000/6000/8000	mmol/L	0.42	11.9%	<5	<100
DiaSys \ ACS - ACOD					
(P)磷					
Beckman AU series	mmol/L	1.04	2.2%	92	3805
Dedicated \ Phosphomolybdate UV-rate					
Abbott ARCHITECT ci series	mmol/L	1.04	4.9%	9	244
BIOSINO \ Phosphomolybdate UV-rate					
TOSHIBA TBA	mmol/L	0.87	12.2%	9	191
MedicalSystem \ Phosphomolybdate-UV-end point					
Beckman AU series	mmol/L	0.999	1.6%	8	184
Szybio \ Phosphomolybdate-UV-end point					
Roche cobas 4000/6000/8000	mmol/L	0.997	2.9%	7	<100
Dedicated \ Phosphomolybdate-UV-end point					
Roche Modular P series	mmol/L	0.994	4.3%	5	<100
G-CELL \ Phosphomolybdate-UV-end point					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
14 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT c series MAKER \ Phosphomolybdate-UV-end point	mmol/L	1.01	3%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Phosphomolybdate UV-rate	mmol/L	1.03	3.2%	5	<100
Roche cobas c 501/502/701/702 Roche \ Phosphomolybdate-UV-end point	mmol/L	1.01	2.2%	5	145
Roche Hitachi/Hitachi Systems Wako \ Phosphomolybdate-UV-end point	mmol/L	0.996	0.6%	5	137
TOSHIBA TBA chongqiong zybio \ Phosphomolybdate-UV-end point	mmol/L	1.08	1%	5	117
Mindray BC-2000 series Mindray \ OTHER	mmol/L	0.979	2.5%	5	100
Beckman UniCel DxC 600/800 Beckman Coulter \ OTHER	mmol/L	1.07	3.3%	49	1018
Beckman AU series BSBE \ OTHER	mmol/L	1.06	2.7%	<5	<100
Roche Hitachi/Hitachi Systems MAKER \ Phosphomolybdate-UV-end point	mmol/L	1.1	1.7%	<5	116
Roche cobas 4000/6000/8000 Roche \ Phosphomolybdate-UV-end point	mmol/L	1.13	6%	32	839
Mindray BS series Mindray \ OTHER	mmol/L	0.899	2%	<5	<100
Beckman AU series HuaChenBio \ Phosphomolybdate UV-rate	mmol/L	0.985	9.4%	28	476
Roche Hitachi/Hitachi Systems MAKER \ OTHER	mmol/L	1.08	3.1%	26	729
(PAMY)胰淀粉酶					
Siemens ADVIA Chemistry XPT Dedicated \ OTHER	U/L	49.4	11.8%	19	910
Beckman AU series Reebio \ OTHER	U/L	44.1	2.7%	11	505
Beckman AU series Beckman Coulter \ Enzymatic	U/L	46.7	9.9%	<5	<100
(PRO/ TP)(总) 蛋白					
Beckman AU series Dedicated \ Biuret-endpoint	g/L	68.9	1.6%	91	3808
DiruiCS series DIRUI \ Biuret-endpoint	g/L	74.8	3.8%	9	<100
Abbott ARCHITECT ci series BIOSINO \ Biuret-rate	g/L	69.1	1.7%	9	243
Beckman AU series Szybio \ Biuret-endpoint	g/L	69.8	1.5%	9	186
Mindray BS series Mindray \ Biuret-endpoint	g/L	69.2	8.9%	9	160
Roche Hitachi/Hitachi Systems MAKER \ Biuret-endpoint	g/L	68.3	1.4%	8	193
Beckman AU series OTHER \ OTHER	g/L	65.9	3.6%	8	187
Roche Hitachi/Hitachi Systems MedicalSystem \ Biuret-endpoint	g/L	67.7	0.9%	7	944
Beckman AU series KeFang(Guangzhou) \ Biuret-rate	g/L	69.3	1.5%	7	156
Roche cobas c 501/502/701/702 Roche \ Biuret-endpoint	g/L	67.9	1.7%	6	176
Roche Modular P series G-CELL \ Biuret-endpoint	g/L	69.5	2%	5	<100
Abbott ARCHITECT c series MAKER \ Biuret-endpoint	g/L	69	1.9%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Biuret-endpoint	g/L	68.8	0.7%	5	<100
Mindray BS series Dedicated \ Biuret-endpoint	g/L	69.8	1.6%	5	295
TOSHIBA TBA BSBE \ Biuret-endpoint	g/L	69.2	2.2%	5	152
Roche Hitachi/Hitachi Systems Wako \ Biuret-endpoint	g/L	70.9	1%	5	137
TOSHIBA TBA chongqiong zybio \ Biuret-endpoint	g/L	68.4	1.4%	5	117
Mindray BC-2000 series Mindray \ Biuret-endpoint	g/L	71.6	1.4%	5	103
Beckman UniCel DxC 600/800 Beckman Coulter \ Biuret-rate	g/L	68.5	1.8%	49	1016
(PROG)孕酮					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.384	9.8%	<5	<100
(TBA)总胆汁酸					
Beckman AU series DongOu \ Enzymatic Cycling Assays Method	μmol/L	13.3	3.3%	92	3786
Abbott ARCHITECT ci series BIOSINO \ Enzymatic Cycling Assays Method	μmol/L	11.4	4.3%	9	243
Siemens ADVIA 1800/2400 Szybio \ Enzymatic Cycling Assays Method	μmol/L	11.1	5.8%	9	136
Roche Hitachi/Hitachi Systems OTHER \ Enzymatic Cycling Assays Method	μmol/L	13.4	9.8%	8	164
Beckman AU series KeFang(Guangzhou) \ Enzymatic Cycling Assays Method	μmol/L	12.9	6.6%	7	113
Roche Hitachi/Hitachi Systems MedicalSystem \ Enzymatic Cycling Assays Method	μmol/L	10.4	9.2%	6	988
TOSHIBA TBA Szybio \ Enzymatic Cycling Assays Method	μmol/L	14.3	5.8%	5	<100

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
15 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche Hitachi/Hitachi Systems	μmol/L	11.1	7%	5	<100
Szybio \ Enzymatic Cycling Assays Method					
TOSHIBA TBA	μmol/L	14.4	4.3%	5	114
chongqiong zybjo \ Enzymatic Cycling Assays Method					
Beckman UniCel DxC 600/800	μmol/L	11.4	6.3%	49	940
DiaSys \ Enzymatic Cycling Assays Method					
TOSHIBA TBA	μmol/L	10.4	3.7%	<5	<100
MedicalSystem \ Enzymatic Cycling Assays Method					
Mindray BS series	μmol/L	14.7	3.5%	<5	<100
Mindray \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	14.3	3.8%	<5	<100
OTHER \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	14.6	6.8%	<5	<100
BSBE \ Enzymatic Cycling Assays Method					
Siemens ADVIA 1800/2400	μmol/L	14.5	3.9%	31	1469
BIOSINO \ Enzymatic Cycling Assays Method					
Beckman AU series	μmol/L	8.51	6.1%	<5	<100
HuaChenBio \ Enzymatic Cycling Assays Method					
Roche cobas 4000/6000/8000	μmol/L	15.4	18.7%	<5	<100
Dedicated \ Enzymatic Cycling Assays Method					
Mindray BS series	μmol/L	14.4	6.7%	<5	<100
Mindray \ OTHER					
Siemens ADVIA Chemistry XPT	μmol/L	13.8	6.1%	16	732
KHB \ Enzymatic Cycling Assays Method					

(TBIL)总胆红素

Beckman AU series	μmol/L	17.2	1.9%	9	267
Beckman Coulter \ OTHER					
Abbott ARCHITECT ci series	μmol/L	19.2	3.9%	9	244
BIOSINO \ Vanadate					
Beckman AU series	μmol/L	15.6	3.9%	9	186
Szybio \ OTHER					
Beckman AU series	μmol/L	19.7	10.5%	8	183
OTHER \ OTHER					
Beckman AU series	μmol/L	20.8	2%	7	141
KeFang(Guangzhou) \ Vanadate					
Beckman AU series	μmol/L	18.3	2.9%	60	2325
Dedicated \ Diazotization					
Roche Hitachi/Hitachi Systems	μmol/L	15.7	6%	6	<100
MAKER \ Vanadate					
KHB ZY series	μmol/L	18.4	1.9%	6	<100
KHB \ Vanadate					
Roche Hitachi/Hitachi Systems	μmol/L	17.1	2.5%	6	<100
BSBE \ Vanadate					
Roche Modular P series	μmol/L	17.2	3.6%	5	<100
G-CELL \ Vanadate					
Mindray BS series	μmol/L	15.5	5.8%	5	<100
Mindray \ Diazotization					
Roche Hitachi/Hitachi Systems	μmol/L	17.5	2.5%	5	<100
Szybio \ Vanadate					
Abbott ARCHITECT c series	μmol/L	15.4	2.7%	5	<100
MAKER \ Vanadate					
Beckman AU series	μmol/L	18	4.9%	5	<100
Beckman Coulter \ Diazotization					
TOSHIBA TBA	μmol/L	16.8	2.6%	5	220
Szybio \ Vanadate					
Roche cobas c 501/502/701/702	μmol/L	15.3	2.1%	5	148
Roche \ Diazotization					
TOSHIBA TBA	μmol/L	20.4	1.9%	5	115
chongqiong zybjo \ Vanadate					
Mindray BC-2000 series	μmol/L	16.2	3%	5	103
Mindray \ Vanadate					
Beckman UniCel DxC 600/800	μmol/L	21.1	8.3%	49	1014
Beckman Coulter \ Diazotization					

(TES)睾酮

Roche cobas E601/602	ng/mL	3.63	8.7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

(TG)甘油三酯

Roche cobas 4000/6000/8000	mmol/L	1.14	2.7%	9	<100
Dedicated \ GPO-POD					
Abbott ARCHITECT ci series	mmol/L	0.986	3.9%	9	242
BIOSINO \ GK-PK-LDH					
Beckman AU series	mmol/L	1.07	1.6%	7	149
KeFang(Guangzhou) \ GPO-POD					
KHB ZY series	mmol/L	1.04	2.8%	6	<100
KHB \ GPO-POD					
Beckman AU series	mmol/L	1.03	2.6%	58	2769
Dedicated \ Dehydrogenase					
Roche Hitachi/Hitachi Systems	mmol/L	1.05	1.9%	5	992
MedicalSystem \ GPO-POD					
Roche Modular P series	mmol/L	1.07	3.1%	5	<100
G-CELL \ GPO-POD					
Mindray BS series	mmol/L	0.995	3.1%	5	<100
Mindray \ Acetyl - acetone photometry					
Abbott ARCHITECT c series	mmol/L	1.11	4.6%	5	<100
MAKER \ GPO-POD					
Roche Hitachi/Hitachi Systems	mmol/L	1.14	4.2%	5	<100
OTHER \ OTHER					
Roche Hitachi/Hitachi Systems	mmol/L	1.06	2.4%	5	<100
BSBE \ GPO-POD					
TOSHIBA TBA	mmol/L	1.15	4%	5	119
OTHER \ OTHER					

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
16 / 17								

批号：376B805					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA chongqiong zybjo \ GPO-POD	mmol/L	0.932	2.9%	5	116
Mindray BC-2000 series Mindray \ OTHER	mmol/L	1.04	1.6%	5	102
Beckman UniCel DxC 600/800 Beckman Coulter \ OTHER	mmol/L	1.17	2.8%	49	1019
Mindray BS series Mindray \ GPO-POD	mmol/L	1.04	2.5%	<5	<100
Beckman AU series BSBE \ GPO-POD	mmol/L	1.12	3%	<5	<100
Beckman AU series Szybio \ OTHER	mmol/L	1.05	1.8%	<5	<100
DiruiCS series DIRUI \ GPO-POD	mmol/L	1.04	4.6%	<5	<100
(THEO)茶碱					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/mL	12.2	9.9%	<5	<100
(TRF)转铁蛋白					
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	258	3.4%	<5	<100
Siemens ADVIA 1800/2400 DiaSys \ Immunoturbidimetry (ITA)	mg/dL	275	2.8%	<5	125
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	280	2.5%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/dL	289	9.7%	<5	<100
(TSH)促甲状腺素					
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	mIU/L	1.47	9.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	1.79	9.5%	<5	<100
(TT3)总三碘甲状原氨酸					
Mindray BS series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.613	13%	9	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.14	10%	8	144
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.378	12.1%	<5	<100
(TT4)总T4					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	56.3	10.1%	<5	<100
(UA)尿酸					
Beckman AU series Dedicated \ Uricase colorimetry	µmol/L	309	1.6%	91	3796
Abbott ARCHITECT ci series BIOSINO \ Phosphotungstate colorimetry	µmol/L	316	2.2%	9	243
Roche Hitachi/Hitachi Systems Whitman \ Uricase colorimetry	µmol/L	313	4%	8	201
Beckman AU series Szybio \ OTHER	µmol/L	309	1.6%	8	184
Roche cobas 4000/6000/8000 Dedicated \ Uricase colorimetry	µmol/L	278	2.6%	8	101
Roche Hitachi/Hitachi Systems MedicalSystem \ Uricase colorimetry	µmol/L	313	3.6%	7	801
Beckman AU series KeFang(Guangzhou) \ Uricase-UV	µmol/L	309	2.4%	7	125
Mindray BS series Mindray \ Uricase colorimetry	µmol/L	295	2%	7	104
Roche Hitachi/Hitachi Systems OTHER \ Uricase colorimetry	µmol/L	298	1.4%	7	101
Roche cobas c 501/502/701/702 Roche \ Uricase colorimetry	µmol/L	295	1.5%	6	170
Roche Modular P series G-CELL \ Uricase colorimetry	µmol/L	298	2%	5	<100
Mindray BS series Mindray \ dry chemistry	µmol/L	309	2.9%	5	<100
Abbott ARCHITECT c series MAKER \ Uricase-UV	µmol/L	297	2.2%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Uricase colorimetry	µmol/L	315	1.4%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Uricase colorimetry	µmol/L	305	1%	5	<100
Mindray BS series Dedicated \ Uricase colorimetry	µmol/L	304	2.8%	5	294
Roche Hitachi/Hitachi Systems Wako \ Uricase colorimetry	µmol/L	298	0.6%	5	137
TOSHIBA TBA chongqiong zybjo \ Uricase colorimetry	µmol/L	293	1.9%	5	115
TOSHIBA TBA OTHER \ Uricase colorimetry	µmol/L	300	1%	5	106
(UIBC)不饱和铁结合力					
Roche cobas c 501/502/701/702 Dedicated \ Ferrozine colorimetry	µmol/L	42.6	3.2%	6	119
Beckman AU series Dedicated \ Ferrozine colorimetry	µmol/L	43.6	2.6%	31	1449

REF	581-1	6×5mL	EXP	2023-04-30	LOT	376B805	更新日期：	2024-04
17 / 17								

		批号：376B805			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Siemens ADVIA 1800/2400 DiaSys \ Ferrozine colorimetry	μmol/L	47.1	2.3%	10	475
Beckman UniCel DxI 600/800 Beckman Coulter \ Ferrozine colorimetry	μmol/L	41.8	9.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Ferrozine colorimetry	μmol/L	41.7	9.8%	<5	<100
(Urea)尿素					
Beckman AU series Beckman Coulter \ Urease colorimetry	mmol/L	4.57	3%	9	268
Abbott ARCHITECT ci series BIOSINO \ Diacetyl monoxime	mmol/L	3.83	3.5%	9	245
Beckman AU series OTHER \ OTHER	mmol/L	4.2	8.4%	9	134
Beckman AU series Szybio \ OTHER	mmol/L	4.2	3.5%	8	188
Roche cobas 4000/6000/8000 Dedicated \ Urease-UV-rate	mmol/L	4.46	5.2%	7	<100
Beckman AU series KeFang(Guangzhou) \ Urease-UV-rate	mmol/L	4.48	3.7%	7	140
KHB ZY series KHB \ Urease-UV-rate	mmol/L	4.43	6%	6	<100
Roche Modular P series G-CELL \ Urease-UV-rate	mmol/L	4.12	3.4%	5	<100
Abbott ARCHITECT c series MAKER \ Urease-UV-rate	mmol/L	4.18	2.5%	5	<100
Siemens ADVIA 1800/2400 LEADMAN \ Urease-UV-rate	mmol/L	4.43	1.9%	5	<100
Roche Hitachi/Hitachi Systems BSBE \ Urease colorimetry	mmol/L	4.24	5.7%	5	<100
DiruiCS series DIRUI \ Diacetyl monoxime	mmol/L	4.01	9%	5	<100
TOSHIBA TBA OTHER \ Urease-UV-rate	mmol/L	4.28	2.9%	5	231
Roche cobas c 501/502/701/702 Roche \ Urease-UV-rate	mmol/L	4.22	2.1%	5	143
Roche Hitachi/Hitachi Systems Wako \ Urease-UV-rate	mmol/L	4.52	3.5%	5	138
Roche Hitachi/Hitachi Systems OTHER \ Urease-UV-rate	mmol/L	4.34	4.8%	5	111
Mindray BC-2000 series Mindray \ OTHER	mmol/L	4.31	2.9%	5	100
Mindray BS series Mindray \ OTHER	mmol/L	4.38	1.4%	<5	<100
Beckman AU series G-CELL \ Urease-UV-rate	mmol/L	4.03	4.4%	<5	<100
(VB12)维生素B12					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	421	10%	<5	<100
(Zn)锌					
OTHER 3V-Diagnostic \ Colorimetry	μmol/L	20.4	4.7%	<5	<100
Beckman AU series Szybio \ OTHER	μmol/L	26.9	5%	<5	<100
Beckman AU series BSBE \ OTHER	μmol/L	22.7	6.4%	39	1649
Roche Hitachi/Hitachi Systems BSBE \ Colorimetry	μmol/L	23.3	4.6%	10	137
Beckman AU series Beckman Coulter \ ICP-MS	μmol/L	27.5	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。