

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
1 / 9								

		批号：3387751			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.85	11.5%	<5	<100
(25-OH-VD3)25羟基维生素D3					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	112	2%	<5	<100
(25-OH-Vit D)25羟基维生素D					
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	49.9	17.3%	5	<100
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	125	17.3%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	72.5	8.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	181	8.4%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	51.3	5.6%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	nmol/L	128	5.6%	<5	<100
AB Sciex API Series ABX \ Chemiluminescence Immunoassay (CLIA)	ng/mL	104	0.7%	<5	<100
AB Sciex API Series ABX \ Chemiluminescence Immunoassay (CLIA)	nmol/L	261	0.7%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	104	2%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	nmol/L	259	2%	<5	<100
(AFP)甲胎蛋白					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	265	11.5%	9	187
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	223	5.7%	9	185
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	229	7.5%	7	207
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	223	6.9%	5	<100
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	337	7.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	161	61.1%	<5	139
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	263	2.5%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	232	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	235	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	243	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	193	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	204	3.1%	<5	<100
Maccura i3000 MAKER \ OTHER	ng/mL	266	2.5%	<5	<100
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	5414	12%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	19.6	7.3%	9	204
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	34	12.8%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	25	5.8%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	5.55	2.4%	<5	<100
Maccura i3000 MAKER \ OTHER	IU/mL	30.2	8.1%	<5	<100
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.7	0%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	24	12.2%	9	207
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.7	2.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	21.7	32.8%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.9	5.4%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.3	3.6%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	23.3	17.6%	<5	<100

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
2 / 9								

		批号：3387751			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	34.6	10%	<5	<100
Maccura i3000 MAKER \ OTHER	IU/mL	18.1	3.2%	<5	<100
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	8.3	0%	<5	<100
(CA 125)糖类抗原 125					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	83.4	3.5%	9	209
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	126	7.3%	8	163
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	147	4.2%	7	208
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	87.2	8.4%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	95.6	8.5%	<5	147
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	168	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	184	1.2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	146	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	132	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	107	4.7%	<5	<100
Maccura i3000 MAKER \ OTHER	IU/mL	170	0.9%	<5	<100
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA)	IU/mL	54.6	0%	<5	<100
(CA 15-3)糖类抗原 15-3					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	81.8	9.9%	9	206
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	106	6.4%	8	162
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	41.1	14.1%	7	209
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	76.9	13.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	79	6.3%	<5	136
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	104	2%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.6	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	81.9	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	31.6	7.2%	<5	<100
Maccura i3000 MAKER \ OTHER	IU/mL	106	2.1%	<5	<100
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.8	0%	<5	<100
(CA 19-9)糖类抗原 19-9					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	100	8.1%	9	208
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	90.1	7.3%	7	209
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	95.7	11.8%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	105	7.5%	<5	153
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	79.5	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	195	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	116	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	69.5	1.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	208	5.6%	<5	<100
Maccura i3000 MAKER \ OTHER	IU/mL	69.9	2%	<5	<100
Roche cobas e 801 Roche \ Chemiluminescence Immunoassay (CLIA)	IU/mL	106	0%	<5	<100
(CEA)癌胚抗原					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.7	6.6%	9	212
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28.5	13.1%	9	157

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
3 / 9								

		批号：3387751			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA) OTHER	ng/mL	16.8	8.5%	7	207
Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA) maglumi system	ng/mL	19.2	8.4%	5	<100
Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.4	3.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.8	7.3%	<5	153
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	30.9	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	27.8	1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	43.4	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.9	1.6%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.6	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.3	2%	<5	<100
Maccura i3000 MAKER \ OTHER	ng/mL	27.5	0.6%	<5	<100
(CK)肌酸激酶					
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	118	12%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	32.1	6%	6	151
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	27.9	3.6%	<5	139
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23	1.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	47.3	12%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	27.9	12%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	29.8	12%	<5	<100
(Cortisol)皮质醇					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	375	7.9%	8	177
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	398	13.6%	<5	152
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	324	12%	<5	<100
Autolumo A6000 Series Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	365	4%	<5	<100
(C-P)C-肽					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.72	3.8%	6	151
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.25	19.7%	5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.63	8.2%	5	111
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.53	5.4%	<5	<100
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.37	5.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.37	31.2%	<5	175
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.67	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.01	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.21	12%	<5	<100
(DGXN)地高辛					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.02	12%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9973	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9378	12%	<5	<100
(E2)雌二醇					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1472	7.3%	8	189
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1762	2.8%	6	152
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1871	11.3%	<5	<100

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
4 / 9								

		批号：3387751			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1476	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1525	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	436	2.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1490	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1770	1.5%	<5	<100
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	150	12%	<5	<100
(FA)叶酸					
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.6	7.8%	8	157
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.2	6.7%	6	150
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.6	7.6%	<5	<100
(Ferritin)铁蛋白					
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	567	9.9%	8	163
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	578	4.6%	7	175
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	497	4.3%	6	150
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	515	6.9%	5	<100
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	459	1.9%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	403	5.1%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	653	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	424	12%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	452	12%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.28	3.5%	6	152
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.87	20.8%	<5	<100
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.1	3.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.6	7.6%	<5	168
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.9	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.93	12%	<5	<100
(FSH)促卵泡刺激素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	56.7	4.6%	9	201
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	67	5.7%	8	211
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	57.6	8.5%	<5	177
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	39.6	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	46	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	58.4	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	49.5	10%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	66.3	12%	<5	<100
(FT3)游离T3					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	21.8	8.5%	9	216
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20.5	6.4%	7	207
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	27.5	2.4%	6	151
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	26.3	3.7%	5	<100
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	21.8	5.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	27.7	7.5%	<5	154

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
5 / 9								

		批号：3387751			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	22	3.6%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.4	5%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.4	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20.1	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	26.5	1.5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12	12%	<5	<100
Maccura i3000 MAKER \ OTHER	pmol/L	24.2	2.5%	<5	<100
(FT4)游离T4					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.1	11%	9	190
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	70.5	4%	7	201
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	59	4.3%	6	150
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	56.4	6.4%	5	<100
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.5	6.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	55.8	7.9%	<5	190
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.7	1.2%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	43.9	2.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	49.3	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	628	2.4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	43.9	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	52.8	2.5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	41.2	12%	<5	<100
Maccura i3000 MAKER \ OTHER	pmol/L	49.1	2.3%	<5	<100
(hCG)人绒毛膜促性腺激素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	643	6.5%	5	112
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	449	2.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	417	2.3%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	516	10%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	621	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	369	2.1%	<5	<100
(HCY)同型半胱氨酸					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μmol/L	38.6	1.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	40.7	12%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	40	12%	<5	<100
(HGH)人生长激素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	10.3	7.7%	<5	<100
(IgE)免疫球蛋白 IgE					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	276	12%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	254	12%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	204	12%	<5	<100
(INS)胰岛素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	129	3.2%	6	150
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	mIU/L	114	4%	5	102
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	125	7.9%	<5	<100
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	159	16.2%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
6 / 9								

		批号：3387751			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	125	8.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	98.6	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	126	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	119	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.3	6.1%	<5	<100
(LH)促黄体生成素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	58	5.3%	9	198
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	71.6	6.8%	8	212
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	59.2	7.6%	<5	155
Tosoh AIA 900/1800 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	53.5	3%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	73.2	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	43.9	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	54.6	10%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	78.9	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	58.6	1.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	48	12%	<5	<100
(Myo / Mb)肌红蛋白					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	255	2.6%	6	150
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	μg/L	273	5.5%	<5	<100
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	387	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	214	17.8%	<5	125
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	233	1%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	292	12%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	281	2.9%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2366	2.5%	6	151
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	278	2.5%	6	151
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2385	8.3%	5	<100
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	280	8.3%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2382	3.6%	<5	122
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	280	3.6%	<5	122
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pg/mL	2654	1.7%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	312	1.7%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2232	12%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	262	12%	<5	<100
(PRL)催/ 泌乳素					
Siemens Immulite1000/2000 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	30.3	3.9%	8	212
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	35	3.6%	6	149
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	40.7	16.6%	<5	157
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	34.2	1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	35.8	12%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	30.4	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	35.7	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	38.6	12%	<5	<100

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
7 / 9								

批号：3387751					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	29.1	12%	<5	<100
(PROG)孕酮					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	17.1	6.4%	8	212
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22.6	3.9%	6	148
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.2	9.8%	<5	142
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.4	4.6%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.6	10%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.3	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.7	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.7	3.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.7	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.6	2.1%	<5	<100
(PSA)前列腺特异性抗原					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.8	5.8%	7	207
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.8	3%	6	149
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.4	9.3%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22.7	7.9%	<5	140
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28.8	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.6	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.6	1.4%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.9	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20.1	12%	<5	<100
(PTH)甲状旁腺素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	22.5	6.8%	8	179
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	24	9.8%	7	204
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	23.5	24.7%	<5	<100
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.5	2.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	20.5	8.6%	<5	144
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	3.8	11.5%	<5	<100
Autolumo A6000 Series Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	22.9	5.9%	<5	<100
(TES)睾酮					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.8	54.5%	9	203
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.27	7.6%	8	205
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.1	86.2%	<5	189
Tosoh AIA 900/1800 OTHER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.38	4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.49	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.02	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.16	12%	<5	<100
(Tg)甲状腺球蛋白					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.66	20.5%	7	206
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.6	66.7%	<5	138
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	16.4	10%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
8 / 9								

批号：3387751					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(THEO)茶碱					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	5.51	2.2%	<5	<100
(Tnl)肌钙蛋白I					
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.61	3.8%	<5	<100
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.22	26.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.73	1%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.34	0.3%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.63	2.6%	6	152
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.74	14.5%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.01	13.5%	<5	100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2	12%	<5	<100
(TSH)促甲状腺素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	26.2	7%	8	174
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.3	3.7%	7	208
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.7	2.8%	6	153
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	24.4	15.2%	<5	<100
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	25	2.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	35.5	6.4%	<5	208
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	37.6	2.1%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.5	4.5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29.2	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	26.2	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.8	1.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29.7	12%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.6	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	26.4	2.7%	<5	<100
Maccura i3000 MAKER \ OTHER	mIU/L	30.9	2%	<5	<100
(TT3)总三碘甲状原氨酸					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.27	4.8%	7	197
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.02	4%	6	147
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.82	6.2%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4	8.7%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.57	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.79	0.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.94	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.58	5.1%	<5	<100
(TT4)总T4					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	178	8.5%	7	208
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	155	3.3%	6	151
OTHER Technicon \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	161	17.8%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	157	7.4%	<5	128
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	131	2.5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	134	12%	<5	<100

REF	540-3	6×3mL	EXP	2026-09-23	LOT	3387751	更新日期：	2024-04
9 / 9								

批号：3387751					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	173	6.2%	<5	<100
(VB12)维生素B12					
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	800	6.8%	8	174
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	900	2.6%	6	150
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pg/mL	992	3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	927	6.3%	<5	124
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	430	3.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	986	12%	<5	<100
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	2798	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2838	12%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.502	4.7%	6	108
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.417	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.527	12%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.486	1.7%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。