

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
1 / 6								

		批号：33877B1			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Agilent aglient \ High performance liquid chromatography (HPLC)	ng/mL	5.38	10%	<5	<100
(25-OH-Vit D)25羟基维生素D					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	68.2	10%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	170	10%	<5	<100
(AFP)甲胎蛋白					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	229	3.8%	7	177
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	270	3.7%	6	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Electro chemi luminescence (ECL)	ng/mL	284	6%	<5	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	225	4.3%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	245	1.8%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	237	10.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	236	10%	<5	<100
chemclin LICA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	267	0.3%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	260	1.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	216	10%	<5	<100
Hotgen C2000 Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	244	2.6%	<5	<100
(ALD)醛固酮					
Agilent aglient \ Chemiluminescence Immunoassay (CLIA)	pg/mL	367	10%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.4	13.3%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	5.53	2.5%	<5	<100
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	IU/mL	5.42	7.6%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	35.9	1.5%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	8.69	15.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	18.2	12%	<5	<100
chemclin LICA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	27.6	2.2%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	18	30.4%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.2	3.8%	<5	<100
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	IU/mL	17.3	6.8%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	35.1	5.7%	<5	<100
chemclin LICA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	26.2	1.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.1	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.1	0.2%	<5	<100
(ASD)雄烯二酮					
Agilent aglient \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.96	10%	<5	<100
(CA 125)糖类抗原 125					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	123	3.8%	7	177
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	92.9	4.4%	6	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	93.4	2.9%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	144	2.2%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	104	9.9%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	80.3	1.6%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
2 / 6								

		批号：33877B1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	160	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	89.2	10%	<5	<100
Hotgen C2000 Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	69.5	1.4%	<5	<100
(CA 15-3)糖类抗原 15-3					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	84.1	5.6%	7	175
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.4	7.4%	6	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.9	5%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	33.8	3.7%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.1	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	70.6	1.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	59.7	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	84.7	10%	<5	<100
Hotgen C2000 Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	79	3.2%	<5	<100
(CA 19-9)糖类抗原 19-9					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	140	3.3%	6	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/mL	140	2.5%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	103	2.5%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	201	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	82.4	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	122	0.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	92.6	10%	<5	<100
(CEA)癌胚抗原					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28	3.8%	7	177
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	22.2	13.5%	6	<100
Siemens ADVIA Centaur XP/XPT Dedicated \ Electro chemi luminescence (ECL)	ng/mL	20.4	4.6%	<5	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.7	1.8%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	14.8	1.7%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19	10.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	42.3	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22.6	1.5%	<5	<100
Hotgen C2000 Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	14	2.3%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	22.2	0.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.3	10%	<5	<100
(CK)肌酸激酶					
TOSHIBA TBA TOSHIBA \ Rate -AMP buffer	U/L	131	13.5%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
TOSHIBA TBA TOSHIBA \ Immune suppression	U/L	23.6	1.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	45.3	10%	<5	<100
(Cortisol)皮质醇					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	424	5.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Enzyme immunoassay (EIA)	ng/mL	324	10%	<5	<100
(C-P)C-肽					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.76	5.4%	6	101
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.98	4.6%	<5	103
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.44	3.4%	<5	<100

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
3 / 6								

批号：33877B1					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.003	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.45	10%	<5	<100
(DGXN)地高辛					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.2	10%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Tosoh AIA 900/1800 TOSOH \ Enzyme immunoassay (EIA)	ng/mL	9359	10%	<5	<100
(E2)雌二醇					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1995	2.8%	6	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1711	7.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1714	10%	<5	<100
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	132	10%	<5	<100
(FA)叶酸					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.05	5%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20	10%	<5	<100
(Ferritin)铁蛋白					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	587	4.4%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	402	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	443	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1033	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	678	4.7%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.49	4.4%	7	173
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.36	9.1%	6	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.97	10%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.3	0.3%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.45	0.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.44	10%	<5	<100
(FSH)促卵泡刺激素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	51.4	3.6%	6	101
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	41.6	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	46.6	0.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	44	10%	<5	<100
(FT3)游离T3					
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20.1	5.9%	<5	158
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17.9	1.4%	<5	<100
Chemclin LiCA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	pmol/L	20.8	6.1%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.4	9.9%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	22.1	0.9%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	19.4	0.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	21.5	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	23.8	0.6%	<5	<100

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
4 / 6								

		批号：33877B1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	24.9	10%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	24	1.2%	<5	<100
(FT4)游离T4					
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.8	5.9%	<5	156
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	66.2	2.8%	<5	<100
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	pmol/L	40	4.2%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	43.9	10.5%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	41	1.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	51.4	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	47.9	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	51	1.3%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	69.2	1.2%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	56.7	5.2%	<5	<100
(hCG)人绒毛膜促性腺激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	538	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	342	10%	<5	<100
(HCY)同型半胱氨酸					
TOSHIBA TBA TOSHIBA \ Enzymatic Cycling Assays Method	μmol/L	41	2.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Enzyme immunoassay (EIA)	μmol/L	40.4	10%	<5	<100
(HGH)人生长激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	8.4	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	375	0.6%	<5	<100
(HS-CTnI) (高) 超敏肌钙蛋白I					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.74	10%	<5	<100
(IgE)免疫球蛋白 IgE					
Tosoh AIA 900/1800 TOSOH \ Enzyme linked immunosorbent assays (ELISA)	IU/mL	295	10%	<5	<100
(IGF1)胰岛素生长因子1					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	18.1	10%	<5	<100
(INS)胰岛素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	121	4.6%	6	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	129	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	119	10%	<5	<100
(LH)促黄体生成素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	76.7	10.3%	6	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	50.3	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	47.6	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	50.9	1.1%	<5	<100
(Myo / Mb)肌红蛋白					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	232	3.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	305	10%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2172	8.6%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	255	8.6%	<5	<100
(PRL)催/ 泌乳素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	41.3	2.6%	<5	<100

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
5 / 6								

		批号：33877B1			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	31.8	1.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	36.4	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	29	10%	<5	<100
(PROG)孕酮					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	27	8.2%	6	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	25.3	1.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	24.6	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	25.8	10%	<5	<100
(PSA)前列腺特异性抗原					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.4	4.7%	7	175
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	19.4	2.1%	<5	<100
Autolumo A6000 Series Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.6	1%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.6	10.1%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22.2	1.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.4	10%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	30.3	0.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	22.3	10%	<5	<100
(PTH)甲状旁腺素					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pmol/L	14.5	7.6%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17.7	2.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.3	1.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	21	10%	<5	<100
(TES)睾酮					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.63	11%	5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.4	1.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.45	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.09	10%	<5	<100
(Tg)甲状腺球蛋白					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.3	10.2%	<5	<100
(THEO)茶碱					
Roche cobas c 501/502/701/702 Roche \ Chemiluminescence Immunoassay (CLIA)	µg/mL	30.1	10%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.82	10%	<5	<100
(TSH)促甲状腺素					
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	24.3	4.1%	<5	158
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.3	2.3%	<5	<100
Chemclin LiCA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	mIU/L	46.5	2.8%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.5	10.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	34.2	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	33.8	1.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29.6	10%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	36.7	2.2%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	47.1	0.4%	<5	<100

REF	540-3	6×3mL	EXP	2027-05-18	LOT	33877B1	更新日期：	2024-04
6 / 6								

		批号：33877B1			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.9	1.5%	<5	<100
(TT3)总三碘甲状原氨酸					
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.44	4.7%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.95	1.7%	<5	<100
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.42	2%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.84	10.1%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.25	0.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.76	10%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.32	0.8%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.8	0.3%	<5	<100
(TT4)总T4					
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	159	4.3%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	165	3.5%	<5	<100
Chemclin LICA 800 CHIVD \ Chemiluminescence Immunoassay (CLIA)	ng/mL	139	3.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	185	10%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	171	10.9%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	126	0.3%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	184	1.6%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	160	1%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	162	2%	<5	<100
(VB12)维生素B12					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	882	3.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	919	10%	<5	<100
(β2M)β-2-微球蛋白					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2970	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。