

REF	540-3	6x3mL	EXP	2028-11-08	LOT	338781B	更新日期: 2025-03
							1 / 4

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(CK-MB)肌酸激酶-MB (活性)	(Myo / Mb)肌红蛋白	(TnT)肌钙蛋白T
(NT-proBNP)N-末端脑钠肽前体	(DGXN)地高辛	(THEO)茶碱
(TT3)总三碘甲状腺氨酸	(FT3)游离T3	(TT4)总T4
(FT4)游离T4	(TSH)促甲状腺素	(Cortisol)皮质醇
(FSH)促卵泡刺激素	(LH)促黄体生成素	(PROG)孕酮
(PRL)催 / 泌乳素	(TES)睾酮	(E2)雌二醇
(C-P)C-肽	(FA)叶酸	(INS)胰岛素
(VB12)维生素B12	(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白
(CEA)癌胚抗原	(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原
(CA 125)糖类抗原 125	(CA 15-3)糖类抗原 15-3	(CA 19-9)糖类抗原 19-9
(β2M)β-2-微球蛋白	(Ferritin)铁蛋白	(F-PSA)游离前列腺特异性抗原
(HCY)同型半胱氨酸	(NSE)神经元特异性烯醇化酶	(Cyfra 21-1)细胞角蛋白19片段
(CA 72-4)糖类抗原 72-4	(ASD)雄烯二酮	(ALD)醛固酮
(Anti-Tg)抗甲状腺球蛋白抗体	(Anti-TPO)抗甲状腺过氧化物酶抗体	(DHEA-S)硫酸脱氢表雄酮
(EPO)促红细胞生成素	(HGH)人生长激素	(IGF1)胰岛素生长因子1
(PTH)甲状旁腺素	(25-OH-Vit D)25羟基维生素D	(CT)降钙素
(IL6)白介素6	(HS-CTnl) (高) 超敏肌钙蛋白I	

REF	540-3	6x3mL	EXP	2028-11-08	LOT	338781B	更新日期:	2025-03
								2 / 4

批号: 338781B			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(25-OH-Vit D)25羟基维生素D			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	63.8	51-76.6
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	159	127-191
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	32	32-32
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	80.1	80.1-80.1
(AFP)甲胎蛋白			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	218	174-262
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	235	235-235
(ALD)醛固酮			
Agilent agilent \ Chemiluminescence Immunoassay (CLIA)	pg/mL	333	266-400
(Anti-Tg)抗甲状腺球蛋白抗体			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	4.71	3.77-5.65
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	1	1-1
(Anti-TPO)抗甲状腺过氧化物酶抗体			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	7.04	3.86-10.22
(ASD)雄烯二酮			
Agilent agilent \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.61	4.47-6.75
(CA 125)糖类抗原 125			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	148	118-178
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	110	110-110
(CA 15-3)糖类抗原 15-3			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	75	60-90
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	60.2	60.2-60.2
(CA 19-9)糖类抗原 19-9			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	75.3	60.2-90.4
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	95.8	95.8-95.8
(CA 72-4)糖类抗原 72-4			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	5.03	5.03-5.03
(CEA)癌胚抗原			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	37.8	30.2-45.4
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.4	17.4-17.4
(CK-MB)肌酸激酶-MB (活性)			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	42.3	33.9-50.7
(Cortisol)皮质醇			
Tosoh AIA Series TOSOH \ Enzyme immunoassay (EIA)	ng/mL	272	200-344
(C-P)C-肽			
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.82	2.22-3.42
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.75	2.75-2.75
(CT)降钙素			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	2.78	2.78-2.78
(Cyfra 21-1)细胞角蛋白19片段			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.68	1.68-1.68
(DGXN)地高辛			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.92	2.34-3.5

REF	540-3	6x3mL	EXP	2028-11-08	LOT	338781B	更新日期:	2025-03
								3 / 4

(DHEA-S)硫酸脱氢表雄酮			
Tosoh AIA Series TOSOH \ Enzyme immunoassay (EIA)	ng/mL	7848	6278-9418
(E2)雌二醇			
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1573	1258-1888
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1642	1563-1721
(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	105	84-126
(FA)叶酸			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.4	10.7-16.1
(Ferritin)铁蛋白			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	415	332-498
(F-PSA)游离前列腺特异性抗原			
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.8	5.56-6.04
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.56	6.56-6.56
(FSH)促卵泡刺激素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	39.2	31.4-47
(FT3)游离T3			
Tosoh AIA Series TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	19.3	15.4-23.2
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	24.5	22.8-26.2
(FT4)游离T4			
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	56.4	50.6-62.2
Tosoh AIA Series TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	44.5	35.6-53.4
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	56.6	56-57.2
(hCG)人绒毛膜促性腺激素			
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	IU/L	413	401-425
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	497	398-596
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	470	470-470
(HCY)同型半胱氨酸			
Tosoh AIA Series TOSOH \ Enzyme linked immunosorbent assays (ELISA)	μmol/L	39.6	31.7-47.5
(HGH)人生长激素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	8.1	6.48-9.72
(HS-CTnI) (高) 超敏肌钙蛋白I			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.99	3.99-5.99
(IgE)免疫球蛋白 IgE			
Tosoh AIA Series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	255	204-306
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	IU/mL	321	304-338
(IGF1)胰岛素生长因子1			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	17.1	13.7-20.5
(IL6)白介素6			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1.5	1.5-1.5
(INS)胰岛素			
Tosoh AIA Series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	98.8	79-118.6
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	131	131-131
(LH)促黄体生成素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	41.9	33.5-50.3
(Myo / Mb)肌红蛋白			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	283	226-340

REF	540-3	6x3mL	EXP	2028-11-08	LOT	338781B	更新日期:	2025-03
								4 / 4

(NSE)神经元特异性烯醇化酶			
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.97	0.97-0.97
(NT-proBNP)N-末端脑钠肽前体			
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1921	1537-2305
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	226	181-271
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	4230	4230-4230
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	498	498-498
(PRL)催乳素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	33.7	27-40.4
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	37.6	36.2-39
(PROG)孕酮			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.7	20.6-30.8
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.4	24.2-26.6
(PSA)前列腺特异性抗原			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20	16-24
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.2	20.6-21.8
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.9	16.9-16.9
(PTH)甲状旁腺素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.6	18.9-28.3
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.4	15.4-15.4
(TES)睾酮			
Tosoh AIA Series TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.82	2.26-3.38
(THEO)茶碱			
Roche cobas c 501/502/701/702 Roche \ Chemiluminescence Immunoassay (CLIA)	µg/mL	26.1	20.9-31.3
(TnT)肌钙蛋白T			
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.64	1.31-1.97
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.64	1.64-1.64
(TSH)促甲状腺素			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	27.6	22.1-33.1
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28	27.1-28.9
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	29.7	27.1-32.3
(TT3)总三碘甲状腺原氨酸			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.51	2.01-3.01
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.16	2.53-3.79
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.64	2.63-2.65
(TT4)总T4			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	165	132-198
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	161	129-193
OTHER OTHER \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	153	152-154
(VB12)维生素B12			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	859	689-1029
(β2M)β-2-微球蛋白			
Tosoh AIA Series TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2710	2168-3252

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。