

REF	530-1	6×3mL	EXP	2025-10-20	LOT	I28E6F2	更新日期：	2024-04
1 / 1								

		批号: I28E6F2			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CK-MB mass)肌酸激酶-MB质量					
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.54	2.3%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.89	10.3%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.75	5.4%	<5	<100
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.21	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.46	2.9%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.6	12%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Radiometer AQTSeries Radiometer \ Time-resolved fluorecence immunoassay (TRFIA)	mg/L	141	10%	<5	<100
Roche Cobas integra Roche \ Latex agglutination (LA)	mg/L	0.16	12.5%	<5	<100
(hs-CRP)超敏C反应蛋白					
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/L	1.25	8%	<5	<100
(Myo / Mb)肌红蛋白					
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	57.8	12%	<5	<100
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	45	1.9%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	57.2	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μg/L	56.1	1.2%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	59.4	10.1%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	102	12%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	11.9	12%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pg/mL	96.4	4.2%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	4.2%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	101	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	11.9	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	108	3.6%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.7	3.6%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	107	4.7%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.6	4.7%	<5	<100
(TnI)肌钙蛋白I					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.199	1.1%	8	<100
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.033	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.146	7%	<5	<100
Radiometer AQTSeries Radiometer \ Time-resolved fluorecence immunoassay (TRFIA)	ng/mL	11	10%	<5	<100
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.208	9.6%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。