

质控品名： 特定蛋白液体非定值质控品

REF	521-3	6×3mL	EXP	2024-04-03	LOT	31B7908	更新日期：	2024-04
1 / 3								

批号：31B7908					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Roche cobas 4000/6000/8000 Roche \ Bromocresol green	g/L	57.8	12%	<5	<100
(ASO)抗链球菌溶血素 O 抗体					
Beckman Immage 800 Beckman Coulter \ OTHER	IU/mL	188	5.9%	8	176
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	205	6%	31	678
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	IU/mL	206	4.3%	13	<100
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	IU/mL	199	4.3%	10	292
(AT-III)抗凝血酶III					
Stago STA series Dedicated \ Nephelometry	%	48	12%	<5	<100
(C3)补体 C3					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	205	5.1%	8	191
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	227	10.7%	5	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	231	3.2%	5	143
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	203	4.5%	31	662
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	216	3%	<5	<100
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	227	3%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	247	3%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	209	3.9%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	244	7.9%	16	444
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	213	3.6%	10	293
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	211	1.3%	<5	<100
(C4)补体 C4					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	55.4	5.3%	8	186
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	43.8	10.6%	7	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	48.4	4.9%	5	144
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	54.7	5.5%	31	675
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	48.9	2.6%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	56.5	3.6%	<5	<100
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	44.9	2.3%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	48.5	1.6%	<5	<100
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	55.4	4.5%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	54.6	7.2%	16	424
Beckman Immage 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	56.3	3.7%	10	291
(Cer)铜蓝蛋白					
Siemens BN II Siemens \ OTHER	g/L	0.778	9.1%	16	429
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	g/L	0.452	3.6%	10	<100
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	0.45	12%	<5	<100
(CRP)C反应蛋白					
Beckman Immage 800 Beckman Coulter \ OTHER	mg/dL	7.53	3.6%	8	188
Beckman Immage 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	7.49	5%	31	666
Beckman Immage 800 OTHER \ Immunonephelometry (INA)	mg/dL	7.71	3.5%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	6.37	7%	16	427
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	7.56	4.6%	10	294
(Ferritin)铁蛋白					
Siemens BN II Siemens \ OTHER	ng/mL	289	8.6%	16	432
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	ng/mL	235	1.5%	15	318
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	341	12%	<5	<100

质控品名： 特定蛋白液体非定值质控品

REF	521-3	6×3mL	EXP	2024-04-03	LOT	31B7908	更新日期：	2024-04
2 / 3								

批号：31B7908					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Hap (HPT))触珠蛋白 (结合珠蛋白)					
Siemens BN II Siemens \ OTHER	mg/L	2274	11.1%	16	442
(HbA1c)糖化血红蛋白A1c					
Roche cobas c 501/502/701/702 Dedicated \ immunoturbidimetric	%	10.2	1.5%	5	147
(hs-CRP)超敏C反应蛋白					
TOSHIBA TBA OTHER \ OTHER	mg/L	56.1	6%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	mg/L	60	8.5%	14	272
(IgA)免疫球蛋白 IgA					
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	377	15%	8	<100
Beckman Image 800 Beckman Coulter \ OTHER	mg/dL	338	4.3%	8	193
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	360	2.3%	5	142
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	341	4.5%	31	678
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	334	2.1%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	387	2.7%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	374	1.4%	<5	<100
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	330	3.6%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	447	7.4%	16	432
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	351	2.8%	10	291
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	363	2%	<5	<100
(IgE)免疫球蛋白 IgE					
Roche cobas c 501/502/701/702 MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	59	8.5%	5	148
Siemens BN II Siemens \ OTHER	IU/mL	76.6	5.8%	16	437
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	59.8	12%	<5	<100
(IgG)免疫球蛋白 IgG					
Beckman Image 800 Beckman Coulter \ OTHER	mg/dL	1855	4.6%	9	189
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	1731	13.2%	6	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	2002	2.9%	5	141
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1887	4.4%	31	683
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	2195	3.8%	<5	<100
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	2387	4%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	2172	3.4%	<5	<100
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	1884	3.6%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	2102	7.1%	16	433
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	1926	3%	10	292
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	1756	1.5%	<5	<100
(IgM)免疫球蛋白 IgM					
Beckman Image 800 Beckman Coulter \ OTHER	mg/dL	255	4.8%	8	188
Roche Hitachi/Hitachi Systems sysmex \ Immunoturbidimetry (ITA)	mg/dL	262	4.1%	5	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	253	2.7%	5	142
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	255	4.2%	31	677
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	175	1.9%	<5	<100
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	mg/dL	231	1.3%	<5	<100
Roche Hitachi/Hitachi Systems OTHER \ Immunoturbidimetry (ITA)	mg/dL	265	5.5%	<5	<100
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	255	2.6%	17	<100
Siemens BN II Siemens \ OTHER	mg/dL	293	6.2%	16	441
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	mg/dL	261	2.6%	10	293
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/dL	242	2.4%	<5	<100

质控品名： 特定蛋白液体非定值质控品

REF	521-3	6×3mL	EXP	2024-04-03	LOT	31B7908	更新日期：	2024-04
3 / 3								

批号：31B7908					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Kappa Light Chain)κ轻链抗体					
Beckman Image 800 OTHER \ Immunonephelometry (INA)	g/L	16.3	3.6%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	15.5	12%	<5	<100
(Lambda Light Chain)λ 轻链抗体					
Beckman Image 800 OTHER \ Immunonephelometry (INA)	g/L	8.31	5.4%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	8.13	12%	<5	<100
(MA)微量白蛋白					
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/L	133	5%	5	142
(PAB)前白蛋白					
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	g/L	0.476	2%	7	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	g/L	0.491	1.8%	30	593
Roche cobas 4000/6000/8000 DiaSys \ Immunonephelometry (INA)	g/L	0.48	12%	<5	<100
(RBP)视黄醇结合蛋白					
Roche cobas 4000/6000/8000 WanTai \ Immunonephelometry (INA)	mg/dL	6.56	12%	<5	<100
(RF)类风湿因子					
Beckman Image 800 Beckman Coulter \ OTHER	IU/mL	121	4%	8	142
TOSHIBA TBA OTHER \ OTHER	IU/mL	129	2.8%	6	<100
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	104	1.7%	6	<100
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	97.2	2.7%	5	141
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	111	6.5%	31	667
TOSHIBA TBA MedicalSystem \ Immunoturbidimetry (ITA)	IU/mL	136	3.5%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	89.1	4.8%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	90.9	2.8%	<5	<100
Siemens BN II Siemens \ OTHER	IU/mL	86.4	7.9%	16	430
Beckman Image 800 OTHER \ Immunonephelometry (INA)	IU/mL	112	4.1%	13	<100
Beckman Image 800 Dedicated \ Immunonephelometry (INA)	IU/mL	110	3.4%	10	293
(TRF)转铁蛋白					
Beckman Image 800 Beckman Coulter \ OTHER	mg/dL	353	6.1%	9	152
Siemens BN II Siemens \ OTHER	mg/dL	398	9.9%	16	429
Roche Hitachi/Hitachi Systems LEADMAN \ Immunoturbidimetry (ITA)	mg/dL	381	2.9%	15	318
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/dL	357	12%	<5	<100
(α-1-AT)α-1-抗胰蛋白酶					
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	mg/dL	201	10.4%	<5	<100
(α1M)α1微球蛋白					
Roche cobas c 501/502/701/702 Roche \ immunoturbidimetric	mg/L	90.7	4%	<5	<100
(β2M)β-2-微球蛋白					
TOSHIBA TBA OTHER \ OTHER	ng/mL	6332	1.8%	8	<100
Roche Hitachi/Hitachi Systems BSBE \ Immunoturbidimetry (ITA)	ng/mL	5017	2.6%	30	593
Roche Hitachi/Hitachi Systems OTHER \ OTHER	ng/mL	5422	4.8%	<5	<100
TOSHIBA TBA MedicalSystem \ OTHER	ng/mL	4659	4.2%	<5	<100
Mindray BS series Mindray \ OTHER	ng/mL	5962	5.5%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	ng/mL	5174	4%	14	282
Roche cobas 4000/6000/8000 JS (Shanghai) \ Immunonephelometry (INA)	ng/mL	5125	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。