

REF	670-2	6×1mL	EXP	2023-05-13	LOT	3FEACBB	更新日期：	2024-05
1 / 1								

		批号：3FEACBB			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	69.6	6.7%	8	487
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	73.3	2.6%	5	133
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	73.4	2.8%	5	123
Sysmex CS 5100 Siemens \ Light scattering	Seconds	69.7	3.8%	<5	<100
BE Thrombolyzer X/XRC/XRM Behnk Elektronik (BE) \ Light scattering	Seconds	67	1.7%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	87.5	10.2%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	67.2	10.1%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	%	38.8	6.4%	5	127
Sysmex CS-1300 SUNBIO \ OTHER	%	38.5	6.1%	5	112
Werfen ACL series Werfen \ Nephelometry	%	42.6	9.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	%	45.3	10.2%	<5	<100
Stago STA series Dedicated \ Nephelometry	%	46.8	9.8%	<5	<100
(FIB)纤维蛋白原					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	g/L	1.77	5.8%	8	483
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	g/L	1.88	3.2%	5	132
Sysmex CS-1300 SUNBIO \ OTHER	g/L	1.86	3.2%	5	116
Sysmex CS 5100 Siemens \ Light scattering	g/L	1.89	4.5%	<5	<100
BE Thrombolyzer X/XRC/XRM Behnk Elektronik (BE) \ Light scattering	g/L	1.88	3.4%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	1.74	9.2%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	1.66	9.6%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	27	6.6%	8	486
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	29.7	1.7%	5	123
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	29.7	2.7%	5	115
BE Thrombolyzer X/XRC/XRM Behnk Elektronik (BE) \ Light scattering	Seconds	32.5	1.8%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	26.2	19.2%	<5	102
Stago STA series Dedicated \ Magnetic Beads	Seconds	23.5	10.2%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	28.7	9.4%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	24.3	9.9%	<5	<100
(TT)凝血酶时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	17.5	4.9%	9	516
Sysmex CA 1500/6000/7000 SUNBIO \ OTHER	Seconds	18.5	2.7%	5	113
Sysmex CS-1300 SUNBIO \ OTHER	Seconds	18.5	2.1%	5	102
Sysmex CS 5100 Siemens \ Light scattering	Seconds	18.1	3.2%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	18.6	9.7%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。