

REF	521-1	6×3mL	EXP	2022-06-29	LOT	31B2A23	更新日期：	2019-09
								1 / 2

批号：31B2A23			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(a1-AG)a1 酸性糖蛋白			
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	41.8	33.4-50.2
(a2M)α2巨球蛋白			
Beckman Immage 800 OTHER \ immunoturbidimetric	mg/L	867	689-1,045
(Alb)白蛋白			
Beckman AU series Dedicated \ Bromocresol purple	mg/dL	1909	1,525-2,293
(ASO)抗链球菌溶血素 O 抗体			
Beckman AU series Reebio \ Immunoturbidimetry (ITA)	IU/mL	133	107-159
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	96.7	77.5-115.9
(C3)补体 C3			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	77.5	62.3-92.7
(C4)补体 C4			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	16.9	13.7-20.1
(Cer)铜蓝蛋白			
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	23.2	18.8-27.6
(CRP)C反应蛋白			
Beckman AU series KHB \ Immunoturbidimetry (ITA)	mg/dL	1.46	1.18-1.74
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	1.33	1.09-1.57
Perkinelmer SuperFlex Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/L	13.9	11.1-16.7
Siemens BN II Dedicated \ Immunonephelometry (INA)	mg/dL	1.36	1.12-1.60
(Ferritin)铁蛋白			
Beckman AU series Dedicated \ Immunonephelometry (INA)	ng/mL	49.6	39.8-59.4
(Hap)触珠蛋白			
Beckman AU series Dedicated \ Immunonephelometry (INA)	mg/L	655	525-785
(IgA)免疫球蛋白 IgA			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	116	92-140
(IgE)免疫球蛋白 IgE			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.7	18.3-27.1
(IgG)免疫球蛋白 IgG			
Beckman Immage 800 Dedicated \ Immunonephelometry (INA)	mg/dL	657	527-787
(IgM)免疫球蛋白 IgM			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	84.2	68.0-100.4
(Kappa Light Chain)κ 轻链抗体			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	6.39	5.13-7.65
(Lambda Light Chain)λ 轻链抗体			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	2.84	2.28-3.40
(PAB)前白蛋白			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	16.5	13.3-19.7
(RBP)视黄醇结合蛋白			
Siemens ADVIA 1800/2400 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	2.33	1.87-2.79

REF	521-1	6×3mL	EXP	2022-06-29	LOT	31B2A23	更新日期：	2019-09
								2 / 2

		批号：31B2A23	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(RF)类风湿因子			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	34.8	28.0-41.6
Roche Hitachi/Hitachi Systems Wako \ Immunoturbidimetry (ITA)	IU/mL	26.2	21.0-31.4
(TRF)转铁蛋白			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/dL	133	109-157
(α1 AT)a1 抗胰蛋白酶			
Beckman AU series OTHER \ Immunoturbidimetry (ITA)	mg/dL	67.6	54.2-81.0
(β2M)β-2-微球蛋白			
Beckman AU series Hyssen \ Immunonephelometry (INA)	mg/dL	0.21	0.17-0.25

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。