



REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
1 / 6								

		批号：3384F73	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.78	1.35-2.21
(25-OH-VD3)25羟基维生素D3			
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	19.5	14.8-24.2
(AFP)甲胎蛋白			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	88.9-145
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	106	80.6-131
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	92-150
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	133	101-165
(ALD)醛固酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	167	127-207
(Anti-Tg)抗甲状腺球蛋白抗体			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	9.86	7.49-12.2
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.6	49.1-80.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	56.5	42.9-70.1
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	20.9	15.9-25.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.4	17.8-29
(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.12	2.37-3.87
(CA 125)糖类抗原 125			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	85.7	65.1-106
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.9	49.3-80.5
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	52.7	40.1-65.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	61.4	46.7-76.1
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	62.3	47.3-77.3
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	71.2	54.1-88.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.3	34.4-56.2
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50.8	38.6-63
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	50.6	38.5-62.7
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.1	44.2-72
(CA 19-9)糖类抗原 19-9			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.4	29.9-48.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50.8	38.6-63
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	46.8	35.6-58
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	319	242-396

REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
2 / 6								

批号：3384F73			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(CEA)癌胚抗原			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.6	8.06-13.1
Abbott ARCHITECT i series Abbott \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	11.5	8.74-14.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.5	7.98-13
(CK)肌酸激酶			
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	57.6	43.8-71.4
(CK-MB)肌酸激酶-MB (活性)			
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.72	7.39-12.1
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.2	7.75-12.6
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	12.8-20.8
(Cortisol)皮质醇			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	228	173-283
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	225	171-279
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	234	178-290
(C-P)C-肽			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.74	1.32-2.16
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.22	1.69-2.75
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.46	1.11-1.81
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.35	1.79-2.91
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.28	1.73-2.83
(DGXN)地高辛			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.48	1.88-3.08
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.97	1.5-2.44
(DHEA-S)硫酸脱氢表雄酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μmol/L	17.2	13.1-21.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	16.3	12.4-20.2
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	15.5	11.8-19.2
(E2)雌二醇			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1040	790-1290
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	988	751-1220
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	968	736-1200
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	808	614-1000
(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	55.5	42.2-68.8
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	45.9	34.9-56.9
(FA)叶酸			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.73	5.87-9.59

REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
3 / 6								

		批号：3384F73	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.83	4.43-7.23
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.37	5.6-9.14
(Ferritin)铁蛋白			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µg/L	137	104-170
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	µg/L	153	116-190
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	µg/L	178	135-221
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	µg/L	210	160-260
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	129	98-160
(F-PSA)游离前列腺特异性抗原			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.81	4.42-7.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.92	5.26-8.58
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.43	4.89-7.97
(FSH)促卵泡刺激素			
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	38.9	29.6-48.2
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.8	20.4-33.2
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.8	20.4-33.2
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	28.4	21.6-35.2
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.8	26.4-43.2
(FT3)游离T3			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15	11.4-18.6
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	8.59-14
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	13.9	10.6-17.2
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.3	8.59-14
(FT4)游离T4			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.5	20.1-32.9
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.5	22.4-36.6
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	32.8	24.9-40.7
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	34.2	26-42.4
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	32.2	24.5-39.9
(hCG)人绒毛膜促性腺激素			
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	147	112-182
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	152	116-188
(HCY)同型半胱氨酸			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	µmol/L	24.3	18.5-30.1
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	µmol/L	24.8	18.8-30.8
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µmol/L	22.1	16.8-27.4
(HGH)人生长激素			
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	3.53	2.68-4.38

REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
4 / 6								

批号：3384F73			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxi 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	µg/L	3.53	2.68-4.38
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µg/L	3.97	3.02-4.92
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	4.59	3.49-5.69
(IgE)免疫球蛋白 IgE			
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	84.9	64.5-105
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.6	61.3-99.9
(IGF1)胰岛素生长因子1			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	22.9	17.4-28.4
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	211	160-262
(INS)胰岛素			
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.2	23.7-38.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.4	23.9-38.9
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	25.1	19.1-31.1
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.7	26.4-43
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.1	24.4-39.8
(LH)促黄体生成素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.3	14.7-23.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.3	16.9-27.7
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	16.6	12.6-20.6
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	17.9	13.6-22.2
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	28.5	21.7-35.3
(Myo / Mb)肌红蛋白			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µg/L	143	109-177
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	149	113-185
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	126	95.8-156
(NT-proBNP)N-末端脑钠肽前体			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1070	811-1320
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	127	96.5-157
(PRL)催/ 泌乳素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.2	11.6-18.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.3	12.4-20.2
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.1	10.7-17.5
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.9	12.1-19.7
(PROG)孕酮			
Siemens Immulite1000/2000 Siemens Dade Actin FS \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.23	6.25-10.2
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.6	8.82-14.4
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	8.05-13.2

REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
5 / 6								

批号：3384F73			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.7	8.14-13.3
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	8.51-13.9
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.7	9.65-15.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.5	9.5-15.5
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.4	10.2-16.6
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.5	8.74-14.3
(PTH)甲状旁腺素			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	90.2	68.6-112
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	83.3	63.3-103
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	80.4	61.1-99.7
(SHBG)性激素结合球蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	78.8	59.9-97.7
(TES)睾酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.59	0.448-0.732
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.82	0.623-1.02
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.47	0.358-0.582
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.42	0.319-0.521
(Tg)甲状腺球蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.3	5.55-9.05
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.2	9.27-15.1
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.37	1.8-2.94
(THEO)茶碱			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	24.1	18.3-29.9
(TnI)肌钙蛋白I			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.42	0.319-0.521
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.47	0.358-0.582
(TnT)肌钙蛋白T			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.19	0.906-1.47
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.03	0.782-1.28
(TSH)促甲状腺素			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.78	7.43-12.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	12.3	9.35-15.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.7	8.13-13.3
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.76	6.66-10.9
(TT3)总三碘甲状腺氨酸			
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	4.33	3.29-5.37

REF	540-2	6×3mL	EXP	2024-04-02	LOT	3384F73	更新日期：	2021-11
6 / 6								

批号：3384F73			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	3.78	2.87-4.69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.11	2.36-3.86
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.3	2.51-4.09
(TT4)总T4			
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	175	133-217
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	138	105-171
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	153	116-190
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	153	116-190
(VB12)维生素B12			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	659	501-817
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	567	431-703
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	680	517-843
(Vit D)维生素D			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	85.2	64.8-106
(β2M)β-2-微球蛋白			
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/L	1.39	1.06-1.72
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mg/L	1.44	1.09-1.79
(β-hCG)人绒毛膜促性腺激素β亚基			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	139	106-172
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	158	120-196

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。