

REF	541-2	6×3mL	EXP	2024-05-16	LOT	339D679	更新日期：	2022-06
1 / 1								

批号：339D679			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(ACTH)促肾上腺皮质激素			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	652	496-808
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pg/mL	595	476-714
(ALD)醛固酮			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	171	137-205
(Anti-Tg)抗甲状腺球蛋白抗体			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	1840	1470-2200
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	395	316-474
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	792	634-950
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	600	480-720
(CT)降钙素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	130	104-156
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/L	111	88.8-133
(PCT)降钙素原			
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.4	9.12-13.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.15	3.32-4.98
(Renin)肾素			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	569	455-683
(Tg)甲状腺球蛋白			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	110	83.6-136
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	100	80-120
(TRAb)促甲状腺素受体抗体			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.55	3.46-5.64

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。