



REF	662-3	6×4.5mL	EXP	2024-03-12	LOT	3F29FF8	更新日期：	2024-01
								1 / 3

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。



REF	662-3	6×4.5mL	EXP	2024-03-12	LOT	3F29FF8	更新日期： 2024-01
							2 / 3

批号：3F29FF8			
项目\仪器\试剂方法	单位	均值	+ / - 2 SD
(BASO %)嗜碱性粒细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series			
sysmex \ Calculate	%	52	32-72
Sysmex XS series			
sysmex \ Calculate	%	11	1-21
Sysmex XT series			
sysmex \ Calculate	%	16	1-31
Sysmex XN series			
sysmex \ Calculate	%	11	1-21
(EOS %)嗜酸性粒细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	16	6-26
Sysmex XE series			
sysmex \ Calculate	%	16	1-31
Sysmex XS series			
sysmex \ Calculate	%	16	6-26
Sysmex XT series			
sysmex \ Calculate	%	16	1-31
Sysmex XN series			
sysmex \ Calculate	%	21	1-41
(HCT)血细胞比容/压积			
Mindray BC-6000 series			
Mindray \ Light scattering	%	53.7	43.3-64.1
Sysmex XE series			
sysmex \ Light scattering	%	53.4	43-63.8
Sysmex XS series			
sysmex \ Light scattering	%	53.1	42.7-63.5
Sysmex XT series			
sysmex \ Light scattering	%	53.3	42.7-63.9
Sysmex XN series			
sysmex \ Light scattering	%	53.6	43.4-63.8
(HGB/ Hb)血红蛋白			
Mindray BC-6000 series			
Mindray \ Colorimetry	g/L	186	150-222
Sysmex XE series			
sysmex \ Light scattering	g/L	184	148-220
Sysmex XS series			
sysmex \ Light scattering	g/L	183	149-217
Sysmex XT series			
sysmex \ Light scattering	g/L	186	150-222
Sysmex XN series			
sysmex \ Light scattering	g/L	183	147-219
(LYM %)淋巴细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	20	10-30
Sysmex XE series			
sysmex \ Calculate	%	25	10-40
Sysmex XS series			
sysmex \ Calculate	%	20	10-30
Sysmex XT series			
sysmex \ Calculate	%	25	10-40
Sysmex XN series			
sysmex \ Calculate	%	25	10-40
(MCH)平均红细胞血红蛋白含量			
Mindray BC-6000 series			
Mindray \ Light scattering	pg	33.5	27.3-39.7
Sysmex XE series			
sysmex \ Light scattering	pg	32.4	26.2-38.6
Sysmex XS series			
sysmex \ Light scattering	pg	32.5	26.3-38.7
Sysmex XT series			
sysmex \ Light scattering	pg	32.6	26.2-39
Sysmex XN series			
sysmex \ Light scattering	pg	32.6	26.2-39
(MCHC)平均红细胞血红蛋白浓度			
Mindray BC-6000 series			
Mindray \ Light scattering	g/L	326	262-390
Sysmex XE series			
sysmex \ Light scattering	g/L	323	259-387
Sysmex XS series			
sysmex \ Light scattering	g/L	323	259-387
Sysmex XT series			
sysmex \ Light scattering	g/L	322	258-386
Sysmex XN series			
sysmex \ Light scattering	g/L	325	259-391
(MCV)平均红细胞体积			
Mindray BC-6000 series			
Mindray \ Light scattering	fL	95	77-113
Sysmex XE series			
sysmex \ Light scattering	fL	96	78-114

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。



REF	662-3	6×4.5mL	EXP	2024-03-12	LOT	3F29FF8	更新日期：	2024-01
3 / 3								

Sysmex XS series sysmex \ Light scattering	fL	96	78-114
Sysmex XT series sysmex \ Light scattering	fL	93	75-111
Sysmex XN series sysmex \ Light scattering	fL	95	77-113
(MONO %)单核细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	9	0-19
Sysmex XS series sysmex \ Calculate	%	9	0-19
Sysmex XT series sysmex \ Calculate	%	5	0-10
Sysmex XN series sysmex \ Calculate	%	5	0-10
(NEUT %)嗜中性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	61	41-81
Sysmex XE series sysmex \ Calculate	%	50	30-70
Sysmex XS series sysmex \ Calculate	%	52	32-72
Sysmex XT series sysmex \ Calculate	%	61	41-81
Sysmex XN series sysmex \ Calculate	%	52	32-72
(PLT)血小板			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	447	357-537
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E +9/L)	446	346-546
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E +9/L)	453	343-563
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	447	347-547
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	450	350-550
(RBC)红细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E +12/L)	5.41	4.35-6.47
Sysmex XE series sysmex \ Light scattering	T/L (1E +12/L)	5.44	4.38-6.5
Sysmex XS series sysmex \ Light scattering	T/L (1E +12/L)	5.48	4.42-6.54
Sysmex XT series sysmex \ Light scattering	T/L (1E +12/L)	5.45	4.33-6.57
Sysmex XN series sysmex \ Light scattering	T/L (1E +12/L)	5.45	4.33-6.57
(WBC)白细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	17.4	13.8-21
Sysmex XE series sysmex \ Light scattering	Giga/L (1E +9/L)	17.5	13.9-21.1
Sysmex XS series sysmex \ Light scattering	Giga/L (1E +9/L)	17.6	14-21.2
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	17.6	14-21.2
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	17.7	14.1-21.3

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的变化，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。