

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2023-05-06 | <b>LOT</b> | 376DF15 | 更新日期： | 2020-11 |
|            |       |       |            |            |            |         |       | 1 / 5   |

|                                                   |        |       |            |
|---------------------------------------------------|--------|-------|------------|
| 批号：376DF15                                        |        |       |            |
| 项目\仪器\试剂\方法                                       | 单位     | 均值    | + / - 2 SD |
| <b>(a-HBDH)α-羟丁酸脱氢酶</b>                           |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 280   | 222-338    |
| Siemens \ DGKC                                    |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 291   | 231-351    |
| Roche \ Enzymatic Cycling Assays Method           |        |       |            |
| Beckman UniCel Dx C 600/800                       | U/L    | 289   | 229-349    |
| Beckman Coulter \ Enzymatic Cycling Assays Method |        |       |            |
| <b>(Alb)白蛋白</b>                                   |        |       |            |
| Beckman AU series                                 | g/L    | 66.3  | 53.1-79.5  |
| Beckman Coulter \ Bromocresol green               |        |       |            |
| Roche cobas 4000/6000/8000                        | g/L    | 62.6  | 50.2-75    |
| Roche \ Bromocresol green                         |        |       |            |
| Siemens ADVIA 1800/2400                           | g/L    | 61.6  | 43.1-80.1  |
| Siemens \ Bromocresol green                       |        |       |            |
| <b>(ALP)碱性磷酸酶</b>                                 |        |       |            |
| Beckman AU series                                 | U/L    | 270   | 216-324    |
| Beckman Coulter \ Rate -AMP buffer                |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 237   | 189-285    |
| Roche \ Rate -AMP buffer                          |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 275   | 192-358    |
| Siemens \ Rate -AMP buffer                        |        |       |            |
| <b>(ALT/ SGPT)丙氨酸氨基转移酶</b>                        |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 139   | 97.3-181   |
| Siemens \ UV with P5P                             |        |       |            |
| Beckman AU series                                 | U/L    | 149   | 120-178    |
| Beckman Coulter \ UV with P5P                     |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 143   | 115-171    |
| Roche \ UV with P5P                               |        |       |            |
| <b>(AMY)淀粉酶</b>                                   |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 372   | 260-484    |
| Siemens \ PNP-G7                                  |        |       |            |
| Beckman AU series                                 | U/L    | 365   | 292-438    |
| Beckman Coulter \ Iodine colorimetry              |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 362   | 288-436    |
| Roche \ PNP-G7                                    |        |       |            |
| <b>(Apo-A I)载脂蛋白A I</b>                           |        |       |            |
| Siemens ADVIA 1800/2400                           | g/L    | 2.31  | 1.62-3     |
| Siemens \ immunoturbidimetric                     |        |       |            |
| Roche cobas 4000/6000/8000                        | g/L    | 1.89  | 1.53-2.25  |
| SekiSui \ immunoturbidimetric                     |        |       |            |
| <b>(Apo-B)载脂蛋白B</b>                               |        |       |            |
| Roche cobas 4000/6000/8000                        | g/L    | 1.07  | 0.83-1.31  |
| SekiSui \ immunoturbidimetric                     |        |       |            |
| Siemens ADVIA 1800/2400                           | g/L    | 1.09  | 0.763-1.42 |
| Siemens \ immunoturbidimetric                     |        |       |            |
| <b>(AST/ SGOT)天门冬氨酸氨基转移酶</b>                      |        |       |            |
| Beckman AU series                                 | U/L    | 239   | 191-287    |
| Beckman Coulter \ UV without P5P                  |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 209   | 167-251    |
| Roche \ UV with P5P                               |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 197   | 138-256    |
| Siemens \ UV without P5P                          |        |       |            |
| <b>(Ca)钙</b>                                      |        |       |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 3.02  | 2.11-3.93  |
| Siemens \ Arsenazo III                            |        |       |            |
| Beckman AU series                                 | mmol/L | 3.12  | 2.5-3.74   |
| Beckman Coulter \ BAPTA                           |        |       |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 3.17  | 2.53-3.81  |
| Roche \ BAPTA                                     |        |       |            |
| <b>(CHE)胆碱酯酶</b>                                  |        |       |            |
| Beckman AU series                                 | U/L    | 12000 | 8380-15600 |
| Beckman Coulter \ Butyrylthiocholine Method       |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 12700 | 8870-16500 |
| Roche \ Butyrylthiocholine Method                 |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 11800 | 8230-15300 |
| Siemens \ DTNB                                    |        |       |            |
| <b>(CHOL/ TC)总胆固醇</b>                             |        |       |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 6.33  | 4.43-8.23  |
| Siemens \ Cholesterol oxidase                     |        |       |            |
| Beckman AU series                                 | mmol/L | 5.81  | 4.65-6.97  |
| Beckman Coulter \ Cholesterol oxidase             |        |       |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 5.79  | 4.65-6.93  |
| Roche \ Cholesterol oxidase                       |        |       |            |
| <b>(CK)肌酸激酶</b>                                   |        |       |            |
| Beckman AU series                                 | U/L    | 324   | 259-389    |
| Beckman Coulter \ Urease-UV-rate                  |        |       |            |
| Roche cobas 4000/6000/8000                        | U/L    | 288   | 230-346    |
| Roche \ Rate -AMP buffer                          |        |       |            |
| Siemens ADVIA 1800/2400                           | U/L    | 296   | 207-385    |
| Siemens \ Rate -AMP buffer                        |        |       |            |

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2023-05-06 | <b>LOT</b> | 376DF15 | 更新日期： | 2020-11 |
|            |       |       |            |            |            |         |       | 2 / 5   |

|                                                         |        |       |             |
|---------------------------------------------------------|--------|-------|-------------|
| 批号：376DF15                                              |        |       |             |
| 项目\仪器\试剂\方法                                             | 单位     | 均值    | + / - 2 SD  |
| <b>(CK-MB)肌酸激酶-MB</b>                                   |        |       |             |
| Siemens ADVIA 1800/2400                                 | U/L    | 44.6  | 31.2-58     |
| Siemens \ DGKC                                          |        |       |             |
| Beckman AU series                                       | U/L    | 55.8  | 44.6-67     |
| Beckman Coulter \ Immune suppression                    |        |       |             |
| Roche cobas 4000/6000/8000                              | U/L    | 49.9  | 39.9-59.9   |
| Roche \ Immune suppression                              |        |       |             |
| <b>(Cl)氯</b>                                            |        |       |             |
| Beckman AU series                                       | mmol/L | 135   | 109-161     |
| Beckman Coulter \ ISE (direct)                          |        |       |             |
| Roche cobas 4000/6000/8000                              | mmol/L | 128   | 104-152     |
| Roche \ ISE (indirect)                                  |        |       |             |
| Siemens ADVIA 1800/2400                                 | mmol/L | 138   | 96.6-179    |
| Siemens \ ISE (indirect)                                |        |       |             |
| <b>(Cr)肌酐</b>                                           |        |       |             |
| Siemens ADVIA 1800/2400                                 | μmol/L | 240   | 168-312     |
| Siemens \ Enzymatic                                     |        |       |             |
| Beckman AU series                                       | μmol/L | 252   | 204-300     |
| Beckman Coulter \ Enzymatic                             |        |       |             |
| Roche cobas 4000/6000/8000                              | μmol/L | 258   | 210-306     |
| Roche \ Enzymatic                                       |        |       |             |
| <b>(CRP)C反应蛋白</b>                                       |        |       |             |
| Beckman Immage 800                                      | mg/dL  | 0.417 | 0.335-0.499 |
| Roche \ Immunonephelometry (INA)                        |        |       |             |
| Mindray BS series                                       | mg/dL  | 0.223 | 0.179-0.267 |
| SHSN \ Immunonephelometry (INA)                         |        |       |             |
| <b>(Cu)铜</b>                                            |        |       |             |
| Beckman AU series                                       | μmol/L | 33.9  | 27.2-40.6   |
| Beckman Coulter \ ICP-MS                                |        |       |             |
| <b>(DBIL)直接胆红素</b>                                      |        |       |             |
| Beckman AU series                                       | μmol/L | 16.9  | 13.7-20.1   |
| Beckman Coulter \ Diazotization -Blank Compensation     |        |       |             |
| Roche cobas 4000/6000/8000                              | μmol/L | 17.8  | 14.4-21.2   |
| Roche \ Diazotization -Blank Compensation               |        |       |             |
| Siemens ADVIA 1800/2400                                 | μmol/L | 20.8  | 14.6-27     |
| Siemens \ Vanadate                                      |        |       |             |
| <b>(DGXN)地高辛</b>                                        |        |       |             |
| Beckman UniCel Dxl 600/800                              | ng/mL  | 1.23  | 1.03-1.43   |
| Beckman Coulter \ Chemiluminescence Immunoassay (CLIA ) |        |       |             |
| <b>(E2)雌二醇</b>                                          |        |       |             |
| Roche cobas E601/602                                    | pmol/L | 218   | 176-260     |
| Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)   |        |       |             |
| <b>(FA)叶酸</b>                                           |        |       |             |
| Roche cobas E601/602                                    | ng/mL  | 1.35  | 1.08-1.62   |
| Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)   |        |       |             |
| <b>(Fe)铁</b>                                            |        |       |             |
| Beckman AU series                                       | μmol/L | 36.7  | 29.4-44     |
| Beckman Coulter \ Ferrozine colorimetry                 |        |       |             |
| Roche cobas 4000/6000/8000                              | μmol/L | 34.2  | 27.4-41     |
| Roche \ Ferrozine colorimetry                           |        |       |             |
| <b>(Ferritin)铁蛋白</b>                                    |        |       |             |
| Roche cobas 4000/6000/8000                              | μg/L   | 99.3  | 79.5-119    |
| Roche \ Immunoturbidimetry (ITA)                        |        |       |             |
| <b>(FT3)游离T3</b>                                        |        |       |             |
| Roche cobas E601/602                                    | pmol/L | 3.12  | 2.5-3.74    |
| Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)   |        |       |             |
| <b>(FT4 )游离T4</b>                                       |        |       |             |
| Roche cobas E601/602                                    | pmol/L | 17.5  | 13.9-21.1   |
| Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)   |        |       |             |
| <b>(GGT)L-γ-谷氨酰基转移酶</b>                                 |        |       |             |
| Beckman AU series                                       | U/L    | 154   | 123-185     |
| Beckman Coulter \ γ-glutamyl -4-nitroaniline            |        |       |             |
| Roche cobas 4000/6000/8000                              | U/L    | 146   | 118-174     |
| Roche \ γ-glutamyl -3-carboxy-4-nitroaniline            |        |       |             |
| Siemens ADVIA 1800/2400                                 | U/L    | 149   | 104-194     |
| Siemens \ γ-glutamyl -3-carboxy-4-nitroaniline          |        |       |             |
| <b>(GLDH)谷氨酸脱氢酶</b>                                     |        |       |             |
| Beckman AU series                                       | U/L    | 9.82  | 7.88-11.8   |
| Beckman Coulter \ DGKC                                  |        |       |             |
| Roche cobas 4000/6000/8000                              | U/L    | 8.73  | 6.99-10.5   |
| DiaSys \ Colorimetry                                    |        |       |             |

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2023-05-06 | <b>LOT</b> | 376DF15 | 更新日期： | 2020-11 |
|            |       |       |            |            |            |         |       | 3 / 5   |

|                                                   |        |      |            |
|---------------------------------------------------|--------|------|------------|
| 批号：376DF15                                        |        |      |            |
| 项目\仪器\试剂\方法                                       | 单位     | 均值   | + / - 2 SD |
| <b>(GLU)葡萄糖</b>                                   |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 10.3 | 7.21-13.4  |
| Siemens \ Hexokinase                              |        |      |            |
| Beckman AU series                                 | mmol/L | 9.89 | 7.97-11.8  |
| Beckman Coulter \ Dehydrogenase                   |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 9.76 | 7.84-11.7  |
| Roche \ Hexokinase                                |        |      |            |
| <b>(HCY)同型半胱氨酸</b>                                |        |      |            |
| Roche cobas 4000/6000/8000                        | μmol/L | 24.8 | 18.8-30.8  |
| Roche \ Enzymatic Cycling Assays Method           |        |      |            |
| Beckman UniCel Dx C 600/800                       | μmol/L | 18.4 | 12.9-23.9  |
| Beckman Coulter \ Enzymatic Cycling Assays Method |        |      |            |
| <b>(HDL-C)高密度脂蛋白胆固醇</b>                           |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 1.77 | 1.24-2.3   |
| Siemens \ Direct Measure                          |        |      |            |
| Beckman AU series                                 | mmol/L | 1.53 | 1.22-1.84  |
| Beckman Coulter \ Enzymatic                       |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 1.32 | 1.06-1.58  |
| Roche \ Enzymatic                                 |        |      |            |
| <b>(IgA )免疫球蛋白 IgA</b>                            |        |      |            |
| Beckman Immage 800                                | g/L    | 3.08 | 2.46-3.7   |
| Beckman Coulter \ Immunoturbidimetry (ITA)        |        |      |            |
| <b>(IgG)免疫球蛋白 IgG</b>                             |        |      |            |
| Beckman Immage 800                                | g/L    | 14.5 | 11.7-17.3  |
| Beckman Coulter \ Immunoturbidimetry (ITA)        |        |      |            |
| <b>(IgM)免疫球蛋白 IgM</b>                             |        |      |            |
| Beckman Immage 800                                | g/L    | 1.49 | 1.21-1.77  |
| Beckman Coulter \ Immunoturbidimetry (ITA)        |        |      |            |
| <b>(K)钾</b>                                       |        |      |            |
| Beckman AU series                                 | mmol/L | 6.03 | 4.85-7.21  |
| Beckman Coulter \ ISE (direct)                    |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 5.76 | 4.62-6.9   |
| Roche \ ISE (indirect)                            |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 6.11 | 4.28-7.94  |
| Siemens \ ISE (indirect)                          |        |      |            |
| <b>(LA)乳酸</b>                                     |        |      |            |
| Beckman AU series                                 | mmol/L | 5.88 | 4.72-7.04  |
| Beckman Coulter \ Lactate to pyruvate (L-P)       |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 5.68 | 4.56-6.8   |
| Roche \ Lactate to pyruvate (L-P)                 |        |      |            |
| <b>(LDH)乳酸脱氢酶</b>                                 |        |      |            |
| Siemens ADVIA 1800/2400                           | U/L    | 310  | 217-403    |
| Siemens \ Lactate to pyruvate (L-P)               |        |      |            |
| Beckman AU series                                 | U/L    | 325  | 260-390    |
| Beckman Coulter \ Lactate to pyruvate (L-P)       |        |      |            |
| Roche cobas 4000/6000/8000                        | U/L    | 328  | 262-394    |
| Roche \ Lactate to pyruvate (L-P)                 |        |      |            |
| <b>(LDL-C)低密度脂蛋白胆固醇</b>                           |        |      |            |
| Beckman AU series                                 | mmol/L | 3.58 | 2.86-4.3   |
| Beckman Coulter \ Enzymatic                       |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 3.75 | 3.03-4.47  |
| Roche \ Enzymatic                                 |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 3.64 | 2.55-4.73  |
| Siemens \ Direct Measure                          |        |      |            |
| <b>(Li)锂</b>                                      |        |      |            |
| Beckman AU series                                 | mmol/L | 1.79 | 1.43-2.15  |
| Beckman Coulter \ Colorimetry                     |        |      |            |
| <b>(LIP)脂肪酶</b>                                   |        |      |            |
| Beckman AU series                                 | U/L    | 119  | 97-141     |
| Beckman Coulter \ Glucose oxidase (GOD)           |        |      |            |
| Roche cobas 4000/6000/8000                        | U/L    | 93.7 | 75.3-112   |
| DiaSys \ Enzymatic Cycling Assays Method          |        |      |            |
| <b>(Mg)镁</b>                                      |        |      |            |
| Beckman AU series                                 | mmol/L | 1.43 | 1.15-1.71  |
| Beckman Coulter \ Xylidyl Blue                    |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 1.46 | 0-4.3      |
| Roche \ Xylidyl Blue                              |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 1.13 | 0.723-1.54 |
| Siemens \ Xylidyl Blue                            |        |      |            |
| <b>(Na)钠</b>                                      |        |      |            |
| Siemens ADVIA 1800/2400                           | mmol/L | 157  | 110-204    |
| Siemens \ ISE (indirect)                          |        |      |            |
| Beckman AU series                                 | mmol/L | 156  | 128-184    |
| Beckman Coulter \ ISE (indirect)                  |        |      |            |
| Roche cobas 4000/6000/8000                        | mmol/L | 153  | 125-181    |
| Roche \ ISE (indirect)                            |        |      |            |

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2023-05-06 | <b>LOT</b> | 376DF15 | 更新日期： | 2020-11 |
|            |       |       |            |            |            |         |       | 4 / 5   |

|                                                                                      |        |      |            |
|--------------------------------------------------------------------------------------|--------|------|------------|
| 批号：376DF15                                                                           |        |      |            |
| 项目\仪器\试剂\方法                                                                          | 单位     | 均值   | + / - 2 SD |
| <b>(NEFA)游离脂肪酸</b>                                                                   |        |      |            |
| Roche cobas 4000/6000/8000<br>DiaSys \ ACS - ACOD                                    | mmol/L | 0.72 | 0.56-0.88  |
| <b>(P)磷</b>                                                                          |        |      |            |
| Beckman AU series<br>Beckman Coulter \ Phosphomolybdate UV-rate                      | mmol/L | 3.12 | 2.54-3.7   |
| Roche cobas 4000/6000/8000<br>Roche \ Phosphomolybdate UV-rate                       | mmol/L | 3.09 | 2.47-3.71  |
| Siemens ADVIA 1800/2400<br>Siemens \ Phosphomolybdate UV-rate                        | mmol/L | 3.23 | 2.26-4.2   |
| <b>(PAMY)胰淀粉酶</b>                                                                    |        |      |            |
| Beckman AU series<br>Beckman Coulter \ Enzymatic                                     | U/L    | 184  | 147-221    |
| <b>(PRO/ TP)(总) 蛋白</b>                                                               |        |      |            |
| Beckman AU series<br>Beckman Coulter \ Biuret-endpoint                               | g/L    | 95.9 | 76.8-115   |
| Roche cobas 4000/6000/8000<br>Roche \ Biuret-endpoint                                | g/L    | 94.8 | 76-114     |
| Siemens ADVIA 1800/2400<br>Siemens \ Biuret-endpoint                                 | g/L    | 98.1 | 68.7-128   |
| <b>(PROG)孕酮</b>                                                                      |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | ng/mL  | 1.1  | 0.88-1.32  |
| <b>(TBA)总胆汁酸</b>                                                                     |        |      |            |
| Beckman AU series<br>Beckman Coulter \ Enzymatic                                     | μmol/L | 142  | 114-170    |
| Roche cobas 4000/6000/8000<br>DiaSys \ Enzymatic Cycling Assays Method               | μmol/L | 76.4 | 61.2-91.6  |
| Siemens ADVIA 1800/2400<br>Siemens \ Enzymatic Cycling Assays Method                 | μmol/L | 78.3 | 54.8-102   |
| <b>(TBIL)总胆红素</b>                                                                    |        |      |            |
| Siemens ADVIA 1800/2400<br>Siemens \ Vanadate                                        | μmol/L | 89.2 | 62.4-116   |
| Beckman AU series<br>Beckman Coulter \ Diazotization                                 | μmol/L | 74.1 | 59.3-88.9  |
| Roche cobas 4000/6000/8000<br>DiaSys \ Diazotization                                 | μmol/L | 73.5 | 58.9-88.1  |
| <b>(TES)睾酮</b>                                                                       |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | ng/mL  | 2.78 | 2.21-3.35  |
| <b>(TG)甘油三酯</b>                                                                      |        |      |            |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ GPO-POD                              | mmol/L | 5.12 | 4.1-6.14   |
| Roche cobas 4000/6000/8000<br>Roche \ Free glycerol removed                          | mmol/L | 4.88 | 3.9-5.86   |
| Siemens ADVIA 1800/2400<br>Siemens \ GPO-POD                                         | mmol/L | 4.97 | 3.48-6.46  |
| <b>(THEO)茶碱</b>                                                                      |        |      |            |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA) | μg/mL  | 21.4 | 17.4-25.4  |
| <b>(TRF)转铁蛋白</b>                                                                     |        |      |            |
| Roche cobas 4000/6000/8000<br>Roche \ Immunoturbidimetry (ITA)                       | mg/dL  | 389  | 313-465    |
| <b>(TSH)促甲状腺素</b>                                                                    |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | mIU/L  | 2.59 | 2.09-3.09  |
| <b>(TT3)总三碘甲状腺原氨酸</b>                                                                |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | nmol/L | 0.92 | 0.74-1.1   |
| <b>(TT4)总T4</b>                                                                      |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | nmol/L | 1.39 | 1.13-1.65  |
| <b>(UA)尿酸</b>                                                                        |        |      |            |
| Beckman AU series<br>Beckman Coulter \ Uricase colorimetry                           | μmol/L | 656  | 524-788    |
| Roche cobas 4000/6000/8000<br>Roche \ Uricase colorimetry                            | μmol/L | 669  | 537-801    |
| Siemens ADVIA 1800/2400<br>Siemens \ Uricase colorimetry                             | μmol/L | 688  | 482-894    |
| <b>(UIBC)不饱和铁结合力</b>                                                                 |        |      |            |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Ferrozine colorimetry                | μmol/L | 54.6 | 43.7-65.5  |



|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 581-2 | 6×5mL | EXP | 2023-05-06 | LOT | 376DF15 | 更新日期： | 2020-11 |
| 5 / 5 |       |       |     |            |     |         |       |         |

|                                                                               |        |      |            |
|-------------------------------------------------------------------------------|--------|------|------------|
| 批号：376DF15                                                                    |        |      |            |
| 项目\仪器\试剂\方法                                                                   | 单位     | 均值   | + / - 2 SD |
| Roche cobas 4000/6000/8000<br>Roche \ Ferrozine colorimetry                   | μmol/L | 53.5 | 42.9-64.1  |
| (Urea)尿素                                                                      |        |      |            |
| Siemens ADVIA 1800/2400<br>Siemens \ Urease-UV-rate                           | mmol/L | 13.7 | 9.59-17.8  |
| Beckman AU series<br>Beckman Coulter \ Urease-UV-rate                         | mmol/L | 12.4 | 9.92-14.9  |
| Roche cobas 4000/6000/8000<br>Roche \ Urease colorimetry                      | mmol/L | 12.9 | 10.5-15.3  |
| (VB12)维生素B12                                                                  |        |      |            |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA) | pg/mL  | 593  | 479-707    |
| (Zn)锌                                                                         |        |      |            |
| Beckman AU series<br>Beckman Coulter \ ICP-MS                                 | μmol/L | 33.2 | 26.6-39.8  |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。