

REF	539-2	6×3mL	EXP	2023-04-12	LOT	336C8D2	更新日期：	2021-07
1 / 4								

批号：336C8D2			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AFP)甲胎蛋白			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	88-136
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	126	102-150
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	127	101-153
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	93-141
(Anti-Tg)抗甲状腺球蛋白抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	30.2	24-36.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	66.2	52.4-80
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	6.15	4.67-7.63
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	33.8	25.7-41.9
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.3	23.5-35.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	21.7	16.9-26.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.1	23.5-36.7
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	59.6	47.8-71.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	61.2	48.6-73.8
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	54.3	42.7-65.9
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.2	61.8-98.6
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.6	31.8-47.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	24.9	19.7-30.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50.3	39.9-60.7
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	38.2	29.2-47.2
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	35.7	28.5-42.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	92.4	73.2-112
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	51.7	41.1-62.3
(CEA)癌胚抗原			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.92	7.96-11.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.3	7.9-12.7
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.2	8.6-13.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.88	5.98-9.78
(CK-MB)肌酸激酶-MB (活性)			
Roche cobas c 501/502/701/702 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.86	4.45-7.27

REF

539-2

6×3mL

EXP

2023-04-12

LOT

336C8D2

更新日期： 2021-07

2 / 4

批号：336C8D2

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(Cortisol)皮质醇			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	246	196-296
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	247	197-297
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	225	171-279
(C-P)C-肽			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.62	2.1-3.14
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.55	2.03-3.07
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.55	2.07-3.03
(DHEA-S)硫酸脱氢表雄酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μmol/L	17.2	13.6-20.8
(E2)雌二醇			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	947	751-1140
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1090	874-1300
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	929	733-1120
(FA)叶酸			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.26	6.62-9.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.17	4.11-6.23
(Ferritin)铁蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	201	159-243
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	143	115-171
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	226	174-278
(F-PSA)游离前列腺特异性抗原			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.82	3.86-5.78
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.96	3.98-5.94
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.04	4-6.08
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.21	3.97-6.45
(FSH)促卵泡刺激素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	32.6	25.6-39.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.2	27.2-41.2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	32.7	26.1-39.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	33.9	25.7-42.1
(FT3)游离T3			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.6	8.81-14.4
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.4	8-12.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.5	9.1-13.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	13.5	10.9-16.1

REF

539-2

6×3mL

EXP

2023-04-12

LOT

336C8D2

更新日期： 2021-07

3 / 4

		批号：336C8D2	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.5	8.74-14.3
(FT4)游离T4			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.3	23-37.6
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.6	22.8-36.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.4	24.4-36.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33.2	26-40.4
(Fβ-hCG)游离β-绒毛膜促性腺激素			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	108	84-132
(HCY)同型半胱氨酸			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	24.2	19.2-29.2
(HGH)人生长激素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.38	3.32-5.44
(INS)胰岛素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.5	27.3-41.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.8	28.6-43
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	35.2	27.8-42.6
(LH)促黄体生成素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.8	15.6-24
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.4	16.2-24.6
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	23.1	18.3-27.9
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	22.8	17.4-28.2
(PRL)催/ 泌乳素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.9	13.7-20.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.9	14.9-22.9
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.9	13.7-22.1
(PROG)孕酮			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.3	8.04-12.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.1	8.84-13.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.7	8.5-12.9
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.5	10.3-16.7
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.6	9.8-15.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.5	9.7-15.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.5	10.3-16.7
(PTH)甲状旁腺素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	78.4	61.8-95

REF	539-2	6×3mL	EXP	2023-04-12	LOT	336C8D2	更新日期：	2021-07
4 / 4								

批号：336C8D2			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	82.6	66-99.2
(TES)睾酮			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.56	0.45-0.67
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.531	0.421-0.641
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.43	0.33-0.53
(Tg)甲状腺球蛋白			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.86	7.9-11.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.26	5.7-8.82
(TSH)促甲状腺素			
Abbott ARCHITECT i series Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.68	7.72-11.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.1	7.7-12.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.2	8.6-13.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	12.3	9.35-15.3
(TT3)总三碘甲状腺氨酸			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.36	2.7-4.02
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.65	2.89-4.41
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	4.1	3.12-5.08
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.15	2.39-3.91
(TT4)总T4			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	142	116-168
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	133	107-159
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	143	113-173
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	nmol/L	161	122-200
(VB12)维生素B12			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	711	559-863
(β2M)β-2-微球蛋白			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/L	1.34	1.06-1.62
(β-hCG)人绒毛膜促性腺激素β亚基			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	128	104-152
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	135	107-163

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。