



REF	662-1	6×4.5mL	EXP	2024-01-21	LOT	3F251D6	更新日期：	2023-11
								1 / 3

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。



REF	662-1	6×4.5mL	EXP	2024-01-21	LOT	3F251D6	更新日期： 2023-11
							2 / 3

批号： 3F251D6			
项目\仪器\试剂方法	单位	均值	+ / - 2 SD
(BASO %)嗜碱性粒细胞 %			
Sysmex XS series sysmex \ Calculate	%	11	1-21
Sysmex XT series sysmex \ Calculate	%	53	33-73
Sysmex XN series sysmex \ Calculate	%	5	0-10
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	53	33-73
(EOS %)嗜酸性粒细胞 %			
Sysmex XE series sysmex \ Calculate	%	16	1-31
Mindray BC-6000 series Mindray \ Flow cytometry	%	21	1-41
Sysmex XS series sysmex \ Calculate	%	16	1-31
Sysmex XN series sysmex \ Calculate	%	21	1-41
Sysmex XT series sysmex \ Calculate	%	16	1-31
(HCT)血细胞比容/压积			
Sysmex XT series sysmex \ Light scattering	%	17.8	14.2-21.4
Sysmex XN series sysmex \ Light scattering	%	17.9	14.3-21.5
Sysmex XS series sysmex \ Light scattering	%	17.7	14.1-21.3
Mindray BC-6000 series Mindray \ Light scattering	%	17.7	14.1-21.3
Sysmex XE series sysmex \ Light scattering	%	17.9	14.5-21.3
(HGB/ Hb)血红蛋白			
Mindray BC-6000 series Mindray \ Colorimetry	g/L	64	51.8-76.2
Sysmex XE series sysmex \ Light scattering	g/L	64	51.8-76.2
Sysmex XS series sysmex \ Light scattering	g/L	64	51.6-76.4
Sysmex XT series sysmex \ Light scattering	g/L	64	51.8-76.2
Sysmex XN series sysmex \ Light scattering	g/L	64	51.6-76.4
(LYM %)淋巴细胞 %			
Sysmex XT series sysmex \ Calculate	%	15	0-30
Sysmex XN series sysmex \ Calculate	%	15	5-25
Sysmex XS series sysmex \ Calculate	%	15	5-25
Sysmex XE series sysmex \ Calculate	%	10	0-20
Mindray BC-6000 series Mindray \ Flow cytometry	%	15	0-30
(MCH)平均红细胞血红蛋白含量			
Sysmex XT series sysmex \ Light scattering	pg	25	20.4-29.6
Sysmex XN series sysmex \ Light scattering	pg	24.8	20-29.6
Sysmex XE series sysmex \ Light scattering	pg	24.5	19.7-29.3
Sysmex XS series sysmex \ Light scattering	pg	24.8	20.2-29.4
Mindray BC-6000 series Mindray \ Light scattering	pg	24.6	19.8-29.4
(MCHC)平均红细胞血红蛋白浓度			
Mindray BC-6000 series Mindray \ Light scattering	g/L	358	296-420
Sysmex XS series sysmex \ Light scattering	g/L	355	293-417
Sysmex XE series sysmex \ Light scattering	g/L	357	295-419
Sysmex XN series sysmex \ Light scattering	g/L	356	292-420
Sysmex XT series sysmex \ Light scattering	g/L	358	294-422
(MCV)平均红细胞体积			
Mindray BC-6000 series Mindray \ Light scattering	fL	85	69-101
Sysmex XS series sysmex \ Light scattering	fL	89	73-105

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。



REF	662-1	6×4.5mL	EXP	2024-01-21	LOT	3F251D6	更新日期：	2023-11
3 / 3								

Sysmex XE series sysmex \ Light scattering	fL	87	71-103
Sysmex XN series sysmex \ Light scattering	fL	89	73-105
Sysmex XT series sysmex \ Light scattering	fL	86	70-102
(MONO %)单核细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	5	0-10
Sysmex XN series sysmex \ Calculate	%	5	0-10
Sysmex XS series sysmex \ Calculate	%	10	0-20
Sysmex XT series sysmex \ Calculate	%	10	0-20
(NEUT %)嗜中性粒细胞 %			
Sysmex XT series sysmex \ Calculate	%	50	30-70
Sysmex XN series sysmex \ Calculate	%	50	30-70
Sysmex XE series sysmex \ Calculate	%	60	40-80
Sysmex XS series sysmex \ Calculate	%	60	40-80
Mindray BC-6000 series Mindray \ Flow cytometry	%	55	35-75
(PLT)血小板			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	50	30-70
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E +9/L)	55	35-75
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E +9/L)	55	35-75
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	55	35-75
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	50	30-70
(RBC)红细胞			
Sysmex XN series sysmex \ Light scattering	T/L (1E +12/L)	2.46	1.98-2.94
Sysmex XT series sysmex \ Light scattering	T/L (1E +12/L)	2.49	2.01-2.97
Sysmex XS series sysmex \ Light scattering	T/L (1E +12/L)	2.46	2-2.92
Sysmex XE series sysmex \ Light scattering	T/L (1E +12/L)	2.48	2.02-2.94
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E +12/L)	2.46	1.98-2.94
(WBC)白细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	3.38	2.72-4.04
Sysmex XE series sysmex \ Light scattering	Giga/L (1E +9/L)	3.36	2.7-4.02
Sysmex XS series sysmex \ Light scattering	Giga/L (1E +9/L)	3.4	2.72-4.08
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	3.38	2.72-4.04
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	3.37	2.71-4.03

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的变化，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。