

REF	540-3	6×3mL	EXP	2023-04-29	LOT	3387621	更新日期：	2021-03
1 / 6								

批号：3387621			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.46	3.54-5.38
(A II)血管紧张素II			
maglumi system Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/L	188	150-226
(AFP)甲胎蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	248	198-298
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	245	195-295
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	234	188-280
(ALD)醛固酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	351	280-422
maglumi system Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	264	212-316
(Anti-Tg)抗甲状腺球蛋白抗体			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	57.4	45.6-69.2
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	17.3	14.1-20.5
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	28.9	24.9-32.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	44.7	36.3-53.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	20.6	16.6-24.6
(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.53	4.38-6.68
(CA 125)糖类抗原 125			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	167	135-199
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	102	83-121
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	117	95-139
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	103	83-123
(CA 15-3)糖类抗原 15-3			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	66.9	53.7-80.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	64.4	54-74.8
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	67.3	54.1-80.5
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	69.5	54.1-84.9
(CA 19-9)糖类抗原 19-9			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	79.4	64.2-94.6
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	94.5	76.1-113
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	105	85.2-125
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	102	82-122

REF

540-3

6×3mL

EXP

2023-04-29

LOT

3387621

更新日期： 2021-03

2 / 6

批号：3387621

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(CEA)癌胚抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.6	20.4-30.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.9	20.7-31.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	39.7	31.3-48.1
(CK-MB)肌酸激酶-MB (活性)			
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	34.3	27.5-41.1
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	33.5	26.7-40.3
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	37.8	28.6-47
(Cortisol)皮质醇			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	355	283-427
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	335	263-407
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	283	219-347
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	267	215-319
(C-P)C-肽			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.76	3.12-4.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.12	2.5-3.74
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.89	3.09-4.69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.02	3.22-4.82
(DGXN)地高辛			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.38	2.69-4.07
(DHEA-S)硫酸脱氢表雄酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μmol/L	35.8	29.4-42.2
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μmol/L	35.6	29.2-42
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	38.7	31.1-46.3
(E2)雌二醇			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1720	1370-2070
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1630	1300-1950
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1470	1180-1760
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1510	1200-1810
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1690	1300-2070
(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	141	113-169
(FA)叶酸			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.2	13.2-19.2
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.9	15.3-22.5
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.6	11.7-17.5

REF	540-3	6×3mL	EXP	2023-04-29	LOT	3387621	更新日期：	2021-03
3 / 6								

		批号：3387621	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.8	9.6-14
(Ferritin)铁蛋白			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	452	362-542
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	415	321-509
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	654	530-778
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	μg/L	617	469-765
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	603	485-721
(F-PSA)游离前列腺特异性抗原			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.15	7.19-11.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.89	6.31-9.47
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.12	7.16-11.1
(FSH)促卵泡刺激素			
Siemens ADVIA Centaur XP Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	76.4	56.8-96
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	41.9	34.1-49.7
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	42.7	33.7-51.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	54.3	43.1-65.5
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	52.6	41.4-63.8
(FT3)游离T3			
Sysmex XP Series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	21.6	17.4-25.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	19.5	15.3-23.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	19.7	15.5-23.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	26.3	21.1-31.5
(FT4)游离T4			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	43.6	34.4-52.8
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	42.7	33.5-51.9
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	52.4	41.8-63
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	46.7	37.5-55.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	47.4	37.8-57
(HCY)同型半胱氨酸			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	49.4	39.4-59.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	48.9	39.1-58.7
(HGH)人生长激素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	9.23	7.39-11.1
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	8.3	6.57-10

REF

540-3

6×3mL

EXP

2023-04-29

LOT

3387621

更新日期：

2021-03

4 / 6

批号：3387621

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(IgE)免疫球蛋白 IgE			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	289	229-349
Beckman Image 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	224	174-274
(IGF1)胰岛素生长因子1			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	22.6	18.4-26.8
(INS)胰岛素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	127	103-151
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	119	95-143
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	132	106-158
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	119	93-145
(LH)促黄体生成素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	53.8	43-64.6
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	55.6	44.8-66.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	47.2	37.4-57
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	49.6	40.4-58.8
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	50.3	40.1-60.5
(Myo / Mb)肌红蛋白			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	266	210-322
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	287	233-341
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	228	184-272
(NT-proBNP)N-末端脑钠肽前体			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2060	1650-2460
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	243	195-291
(PRL)催/泌乳素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	37.5	29.9-45.1
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	34.6	27.4-41.8
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	36.3	28.9-43.7
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	36.3	29.7-42.9
(PROG)孕酮			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.6	19.1-28.1
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	27.1	21.5-32.7
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.6	20.6-30.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.5	19.5-29.5
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23.6	18.8-28.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.7	17.5-25.9

REF	540-3	6×3mL	EXP	2023-04-29	LOT	3387621	更新日期：	2021-03
5 / 6								

		批号：3387621	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.5	19.7-29.3
(PTH)甲状旁腺素			
Siemens Immulite1000/2000 Dedicated \ Enzyme linked immunosorbent assays (ELISA)	pg/mL	298	240-356
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	187	151-223
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	196	154-238
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	163	129-197
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	169	135-203
(Renin)肾素			
maglumi system Dedicated \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	17.8	14.4-21.2
(rT3)反三碘甲状腺原氨酸			
maglumi system Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.34	1.88-2.8
(SHBG)性激素结合球蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	68.1	51.8-84.4
(TES)睾酮			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.21	1.79-2.63
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.76	1.4-2.12
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.03	2.4-3.66
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.34	1.83-2.85
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.46	1.95-2.97
(Tg)甲状腺球蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.3	8.1-12.5
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.3	10.1-14.5
Siemens Immulite1000/2000 Dedicated \ Enzyme linked immunosorbent assays (ELISA)	ng/mL	4.23	3.31-5.15
(THEO)茶碱			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/mL	30.1	24-36.2
(TnI)肌钙蛋白I			
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.96	3.98-5.94
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.23	4.17-6.29
Sysmex XP Series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	19.3	15.5-23.1
(TnT)肌钙蛋白T			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.12	1.72-2.52
(TSH)促甲状腺素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.5	25.1-37.9
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	28.9	23.3-34.5
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.2	25.4-37
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	30.2	23.8-36.6

REF	540-3	6×3mL	EXP	2023-04-29	LOT	3387621	更新日期：	2021-03
6 / 6								

		批号：3387621	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxi 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.2	24.8-37.6
(TT3)总三碘甲状原氨酸			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	6.44	5.22-7.66
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	6.97	5.49-8.45
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	7.23	5.75-8.71
Sysmex XP Series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	7.21	5.77-8.65
(TT4)总T4			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	227	183-271
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	248	199-297
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	nmol/L	272	215-329
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	186	148-224
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	192	154-230
(VB12)维生素B12			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	846	675-1010
(β2M)β-2-微球蛋白			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mg/L	2.82	2.24-3.4
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/L	2.56	2.02-3.1
(β-hCG)人绒毛膜促性腺激素β亚基			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	402	330-474
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	404	324-484
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	416	332-500

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。