

REF	539-1	6×3mL	EXP	2022-11-15	LOT	336A293	更新日期：	2020-11
								1 / 3

批号：336A293			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.756	0.584-0.928
(AFP)甲胎蛋白			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.25	5.68-8.82
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.16	5.44-8.88
(ALD)醛固酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	79.3	63.4-95.2
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	175	133-217
(Anti-Tg)抗甲状腺球蛋白抗体			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	13.6	10.8-16.4
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	9.6	7.6-11.6
(ASD)雄烯二酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.916	0.71-1.12
(CA 125)糖类抗原 125			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	11.2	9-13.4
(CA 15-3)糖类抗原 15-3			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	13.6	11-16.2
(CA 19-9)糖类抗原 19-9			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	7.83	6.27-9.39
(CEA)癌胚抗原			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.55	3.61-5.49
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.42	2.74-4.1
(Cortisol)皮质醇			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	113	90.5-135
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	122	92.8-151
(C-P)C-肽			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.6	0.46-0.74
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.98	0.8-1.16
(DGXN)地高辛			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.618	0.478-0.758
(DHEA-S)硫酸脱氢表雄酮			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μmol/L	8.25	6.53-9.97
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	8.62	6.55-10.7
(E2)雌二醇			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	207	167-247
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	336	270-402

REF	539-1	6×3mL	EXP	2022-11-15	LOT	336A293	更新日期：	2020-11
								2 / 3

批号：336A293			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(FA)叶酸			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.99	3.18-4.8
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.31	3.45-5.17
(Ferritin)铁蛋白			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	78.8	62.8-94.8
Roche cobas c 501/502/701/702 Dedicated \ Immunoturbidimetry (ITA)	μg/L	180	158-202
(F-PSA)游离前列腺特异性抗原			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.11	0.0835-0.136
(FSH)促卵泡刺激素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	6.49	5.17-7.81
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.33	3.45-5.21
(FT3)游离T3			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.75	3.77-5.73
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.62	4.52-6.72
(FT4)游离T4			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.3	13.1-19.5
(hCG)人绒毛膜促性腺激素			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.57	3.63-5.51
(HCY)同型半胱氨酸			
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	16.6	12.6-20.6
(HGH)人生长激素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	1.62	1.29-1.95
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	2.25	1.71-2.79
(IgA)免疫球蛋白 IgA			
Siemens BN ProSpec Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	2.23	1.73-2.73
(IgE)免疫球蛋白 IgE			
Beckman Image 800 Roche \ Immunoturbidimetry (ITA)	IU/mL	93.2	70.8-116
(IgG)免疫球蛋白 IgG			
Siemens BN ProSpec Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	9.55	7.55-11.6
(IgM)免疫球蛋白 IgM			
Siemens BN ProSpec Dedicated \ Chemiluminescence Immunoassay (CLIA)	g/L	0.865	0.685-1.04
(INS)胰岛素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	7.75	6.15-9.35
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.18	7.34-11
(LH)促黄体生成素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.21	2.53-3.89
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.67	3.75-5.59

REF	539-1	6×3mL	EXP	2022-11-15	LOT	336A293	更新日期：	2020-11
								3 / 3

批号：336A293			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(PRL)催/泌乳素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.43	7.54-11.3
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.26	6.62-9.9
(PROG)孕酮			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.1	0.874-1.33
Roche cobas E601/602 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.969	0.775-1.16
(PSA)前列腺特异性抗原			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.73	0.59-0.87
(PTH)甲状旁腺素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	12.1	9.9-14.3
(SHBG)性激素结合球蛋白			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	283	227-339
(TES)睾酮			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.82	6.2-9.44
Genrui-PA120 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.81	5.48-8.14
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.12	5.68-8.56
(Tg)甲状腺球蛋白			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.3	8.1-12.5
(TSH)促甲状腺素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	1.55	1.21-1.89
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	1.91	1.53-2.29
(TT3)总三碘甲状腺原氨酸			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	1.79	1.43-2.15
(TT4)总T4			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	100	80-120
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	100	79.2-121
(VB12)维生素B12			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	432	348-516
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	267	215-319
(Vit D)维生素D			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	42.7	34.1-51.3
(β2M)β-2-微球蛋白			
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/L	0.753	0.572-0.934
(β-hCG)人绒毛膜促性腺激素β亚基			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.44	3.48-5.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	2.73	2.19-3.27

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。